Sentina® Ambidextrous

Sentina® Ivory

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en Instructions for Use

Dear: Valued Customer, This user information according to provisions of PPE Regulation (EU) 2016/425 is to be used in combination with the specific product related information on each product packaging.

Product: Latex Examination Gloves, Powder Free, Non Sterile, compliant with EN ISO 374 and EN ISO 21420

1. Cautions and Restriction

- Keep in clean, cool and dry area. Protect gloves from freezing and avoid excessive heat. Shield from direct sunlight, fluorescent lighting, x-rays, moisture and ozone.
- For single-use only.
- For latex gloves: This product contains natural latex which may cause allergic reactions.
- The product contains dithiocarbamate – not to be used in the case hypersensitivity to these substances.

2. Instructions before usage

- Before usage, inspect the gloves for any defect or imperfections and its suitability for the intended use. The permeation and penetration resistance are assessed under laboratory conditions and do not reflect the actual duration of protection in the workplace. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. The test results are only relate to the tested specimen and chemicals (mixtures may differ). When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. For more information contact [www.Lohmann-Rauscher.com!](#)

3. Chemical resistance

- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to over 400mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

4. Donning

- Remove all hand and wrist jewelry, and wash the hands before donning.
- Place the gloves on the prepared work surface.

- The user puts a glove on his/her dominant hand by grabbing it with the other hand, remembering to only touch the inside of the gloves, and slipping it over the dominant hand until it reaches finder level.
- The wearer uses the gloved dominant hand to slip the other glove onto the non-dominant hand.
- Once both gloves are on, the users can touch the outside of the gloves to ensure a proper fit.

5. Doffing

- Using the dominant hand, users start by grabbing the outside of the glove on the non-dominant hand on the palm side near the cuff.
- Pull the glove off the non-dominant hand and place it in the gloved hand, balling it up.
- Slip two fingers under the cuff of the other hand glove and carefully peel it off the hand without touching the wrist, turning the remaining glove inside out as it is removed and in turn encasing the first glove.
- The gloves can be disposed.

6. Disposal Method

- Throw it away to the special rubbish can after using.

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