

Confirmation of the antibacterial activity of a superabsorbent wound dressing¹ in vitro

Claudia Feldkamp¹, Martin Abel¹, Cornelia Wiegand²,

¹Lohmann & Rauscher GmbH & Co. KG (Clinical Regulatory Affairs, Neuwied, Germany)

²Department of Dermatology, University Hospital Jena, Jena, Germany



AIM:


The objective of this study was to confirm the antibacterial activity of a superabsorbent wound dressing¹ in vitro

METHOD:

The antibacterial activity against *E. coli* was investigated in vitro according to JIS L 1902:2015 with polyester material as growth control. Furthermore, an experimental *E. coli* biofilm was treated with the dressing for 24h. Biomass was then quantified by determination of the bacterial respiratory activity after removal of the dressings and compared to an untreated control without dressing. The release of endotoxins during treatment of planktonic and biofilm *E. coli* as well as the binding of endotoxins were investigated via absorbance measurements on a plate photometer and compared to a control.

RESULTS:

The superabsorbent¹ wound dressing **reduced the growth of *E. coli* bacteria significantly** (2.65 log cfu vs. control: 0.00 log cfu reduction).



GRAPH 1: Antimicrobial activity

Reduced growth of *Escherichia coli* under the influence of the superabsorbent¹ wound dressing over 24 hours (a) and strong reduction of microbial growth (2.65 log cfu vs. control: 0.00 log cfu) according to JIS L 1902 (b).

THE JOINT EWMA-GNEAUPP 2025 CONFERENCE


BARCELONA, SPAIN • 26 – 28 MARCH 2025

¹Vliwazell® Pro – Lohmann & Rauscher GmbH & Co. KG

Confirmation of the antibacterial activity of a superabsorbent wound dressing¹ in vitro

RESULTS:

Treatment of and experimental E. coli biofilm with the dressing **significantly reduced biomass** (87.8% reduction compared to untreated control).



GRAPH 2: Biofilm
The superabsorbent¹ wound dressing significantly reduced E. coli biofilm biomass compared to the control without dressing (87.8% reduction).

THE JOINT EWMA-GNEAUPP 2025 CONFERENCE

BARCELONA, SPAIN • 26 – 28 MARCH 2025


¹Vliwazell® Pro – Lohmann & Rauscher GmbH & Co. KG

Confirmation of the antibacterial activity of a superabsorbent wound dressing¹ in vitro

RESULTS:

Release of endotoxins was **low** after treatment of planktonic E. coli (1.40 EU/mL; polyester control: 0.05 EU/mL) and an E. coli biofilm (1.38 EU/mL; untreated control: 1.46 EU/mL).

GRAPH 3: Endotoxin release after treatment of planktonic E. coli
Release of endotoxins was low after treatment of planktonic E. coli with the dressing (1.40 EU/mL; polyester control: 0.05 EU/mL).



GRAPH 4: Endotoxin release after treatment of Biofilm E. coli
Release of endotoxins was low after treatment of biofilm E. coli with the dressing (1.38 EU/mL; untreated control: 1.46 EU/mL). No significant differences between endotoxin release after treatment with the dressing and untreated biofilm were found.

THE JOINT EWMA-GNEAUPP 2025 CONFERENCE

BARCELONA, SPAIN • 26 – 28 MARCH 2025

¹Vliwazell® Pro – Lohmann & Rauscher GmbH & Co. KG


Confirmation of the antibacterial activity of a superabsorbent wound dressing¹ in vitro

RESULTS:

The superabsorbent¹ wound dressing **bound significant amounts of E. coli endotoxin** and **reduced its concentration in the simulated wound fluid**. After 1h the endotoxin reduction was already >90% compared to the Control.

CONCLUSION:

Test results indicate that the superabsorbent¹ wound dressing **significantly reduces microbial growth and biomass** in vitro. The dressing **binds endotoxins rapidly and efficiently**, with a **low release of endotoxins** after treatment of planktonic E. coli and an E. coli biofilm.



GRAPH 5: Binding of E. coli endotoxin

The amount of unbound endotoxin in simulated wound fluid was significantly reduced after treatment with the dressing compared to the control.

THE JOINT EWMA-GNEAUPP 2025 CONFERENCE

BARCELONA, SPAIN • 26 – 28 MARCH 2025

¹Vliwazell® Pro – Lohmann & Rauscher GmbH & Co. KG

Confirmation of the antibacterial activity of a superabsorbent wound dressing¹ in vitro

Claudia Feldkamp¹, Martin Abel¹, Cornelia Wiegand²,

¹Lohmann & Rauscher GmbH & Co. KG (Clinical Regulatory Affairs, Neuwied, Germany)

²Department of Dermatology University Hospital Jena, Jena, Germany



AIM:


The objective of this study was to confirm the antibacterial activity of a superabsorbent wound dressing¹ in vitro.

METHOD:

The antibacterial activity against *E. coli* was investigated in vitro according to JIS L 1902:2015 with polyester material as growth control. Furthermore, an experimental *E. coli* biofilm was treated with the dressing for 24h. Biomass was then quantified by determination of the bacterial respiratory activity after removal of the dressings and compared to an untreated control without dressing. The release of endotoxins during treatment of planktonic and biofilm *E. coli* as well as the binding of endotoxins were investigated via absorbance measurements on a plate photometer and compared to a control.


RESULTS:

The superabsorbent¹ wound dressing **reduced the growth of *E. coli* bacteria significantly** (2.65 log cfu vs. control: 0.00 log cfu reduction). Treatment of and experimental *E. coli* biofilm with the dressing **significantly reduced biomass** (87.8% reduction compared to untreated control). **Release of endotoxins was low** after treatment of planktonic *E. coli* (1.40 EU/mL; polyester control: 0.05 EU/mL) and an *E. coli* biofilm (1.38 EU/mL; untreated control: 1.46 EU/mL) and the dressing **bound significant amounts of *E. coli* endotoxin** and **reduced its concentration in the simulated wound fluid**. After 1h the endotoxin reduction was already >90% compared to the Control.




GRAPH 1: Antimicrobial activity

Reduced growth of *Escherichia coli* under the influence of the superabsorbent¹ wound dressing over 24 hours (a) and strong reduction of microbial growth (2.65 log cfu vs. control: 0.00 log cfu) according to JIS L 1902 (b).



GRAPH 2: Biofilm

The superabsorbent¹ wound dressing significantly reduced *E. coli* biofilm biomass compared to the control without dressing (87.8% reduction).




GRAPH 3: Endotoxin release after treatment of planktonic *E. coli*

Release of endotoxins was low after treatment of planktonic *E. coli* with the dressing (1.40 EU/mL; polyester control: 0.05 EU/mL).



GRAPH 4: Endotoxin release after treatment of Biofilm *E. coli*

Release of endotoxins was low after treatment of biofilm *E. coli* with the dressing (1.38 EU/mL; untreated control: 1.46 EU/mL). No significant differences between endotoxin release after treatment with the dressing and untreated biofilm were found



GRAPH 5: Binding of *E. coli* endotoxin

The amount of unbound endotoxin in simulated wound fluid was significantly reduced after treatment with the dressing compared to the control

CONCLUSION:

Test results indicate that the superabsorbent¹ wound dressing **significantly reduces microbial growth and biomass** in vitro. The dressing **binds endotoxins rapidly and efficiently**, with a **low release of endotoxins** after treatment of planktonic *E. coli* and an *E. coli* biofilm.

THE JOINT EWMA-GNEAUPP 2025 CONFERENCE

BARCELONA, SPAIN • 26 – 28 MARCH 2025