

# Protective effect of polihexanide on HaCaT keratinocytes in co-culture with *Staphylococcus aureus*

Friedrich-Schiller-Universität Jena



seit 1558

C. Wiegand<sup>1</sup>, M. Abe<sup>2</sup>, P. Ruth<sup>2</sup>, U.-C. Hipler<sup>1</sup>

<sup>1</sup>Department of Dermatology, University Medical Center, Jena, Germany

<sup>2</sup>Lohmann & Rauscher GmbH & Co. KG, Rengsdorf, Germany

## Introduction

*Staphylococcus aureus* is one of the most important pathogens of nosocomial infections and is a common complication during the treatment of chronic wounds. It can exhibit a range of antibiotic resistance (MRSA, methicillin-resistant *Staphylococcus aureus*). Polihexanide (PHMB) is regarded first choice for the antimicrobial treatment of critical colonized or infected chronic wounds because of its good skin tolerability beside its antimicrobial effects. It possesses a specific mechanism of action against acidic lipids of the bacterial membrane and has only little effect on the neutral lipids of the human cell membrane [1]. In fact, Kramer et al. showed that polihexanide promotes wound healing in an animal model [2]. Furthermore, we have investigated the effect of polihexanide on human keratinocytes and found that polihexanide in low concentrations stimulates cell proliferation [3]. Hence, we have used a co-culture system of HaCaT keratinocytes and *Staphylococcus aureus* to test the capacity of polihexanide to protect the cells from the bacterial damage.




Fig. 1: Negative effect of increasing *Staphylococcus aureus* concentrations on HaCaT keratinocyte viability and proliferation (measurement of keratinocyte ATP content, mean  $\pm$  SE).

## Material & Methods

HaCaT keratinocytes were cultured with increasing concentrations of *Staphylococcus aureus* and with or without the addition of polihexanide as well as the extract of a polihexanide-containing HydroBalanced biocellulose based wound dressing (PHWD)\* in different concentrations. Viability and proliferation of HaCaT keratinocytes was investigated by means of the ATPLite™-M kit (Perkin Elmer). Viable *Staphylococcus aureus* cells were quantified via staining with SYTO-9 (Molecular Probes).




Fig. 2: Addition of 1  $\mu$ g/ml polihexanide to a co-culture of *Staphylococcus aureus* and HaCaT keratinocytes protects the cells against bacterial damage (quantification of keratinocytes by measurement of cellular ATP, mean  $\pm$  SE).

\* PHWD = Suprasorb® X+PHMB, Lohmann & Rauscher




Fig. 3: Polihexanide reduces the number of living *Staphylococcus aureus* cells and inhibits bacterial growth (staining with SYTO-9, mean  $\pm$  SE)

## Results

*Staphylococcus aureus* had a concentration-dependent negative effect on HaCaT cell viability and proliferation (Fig. 1). The addition of polihexanide (1  $\mu$ g/ml) prevented damage to the HaCaT cells and restored normal cell proliferation (Fig. 2). In accordance, the addition of 0.4  $\mu$ g/ml and 1  $\mu$ g/ml of polihexanide, respectively, reduced the number of viable bacterial cells as determined via SYTO-9 staining (Fig. 3). Because polihexanide is often used in wound dressings we have investigated the effect of an extract of PHWD\* in our co-culture system and observed a significant reduction of *Staphylococcus aureus* growth (data not shown). Hence, the extract of PHWD\* was able to protect the keratinocytes from bacterial damage and, furthermore, had a positive influence on the cell proliferation (Fig. 4).




Fig. 4: PHWD\* extract protects HaCaT keratinocytes in co-culture with *Staphylococcus aureus* against bacterial damage (determination of HaCaT cell viability via measurement of ATP content, mean  $\pm$  SE).

## Conclusions

Polihexanide is a highly potent antimicrobial agent, which possesses low cytotoxicity and very good skin tolerability. In addition, it is able to induce cell proliferation *in vitro* [3] as well as in an animal model [2]. Therefore, polihexanide seems to be an ideal antimicrobial substance in wound dressings for treating chronic wounds. Furthermore, we have been able to proof the antimicrobial activity of polihexanide and PHWD\* extract in a co-culture of HaCaT keratinocytes and *Staphylococcus aureus in vitro*. It protects the cells from the bacterial damage and allows normal cell growth.

## References

- Ikeda T. et al. *Biochem Biophys Acta* 1983; 735: 380-6
- Kramer A. et al. *Skin Pharmacol Physiol* 2004; 17: 141-6
- Wiegand C. et al. *GMS Krankenhaushyg Interdisziplinär* 2007; 2(2): Doc43