



People.Health.Care.

Raucolast® Conforming bandage

Field of Application

- Compression bandaging for lymphedema
- Fingers, Toes and Genital lymphedema
- Fixing foam padding
- Dressing retention

Features/Benefits

- Soft, comfortable, closed-weave gauze bandage
- Ideal for bandaging fingers and toes
- Conforms well to body contours
- Air permeable
- Wider size is excellent to secure foam padding from moving under bandages
- Extensibility up to approximately 80% lengthwise

Product composition

100% polyester

Not made with natural rubber latex



Raucolast® conforming bandage

rolled, bulk packed in a cardboard box

Size	Length stretched	Item No.	Shipping Units per box/case
4 cm / 1.6 in.	4 m / 4.4 yd	11493	20/320
6 cm / 2.4 in.		11496	20/200

