



People.Health.Care.

Suprasorb[®] P Suprasorb[®] P Border Foam dressings you can rely on.

Soft and flexible

Comfortable to wear, easy to use.¹

- Easy to mold² — conforms to body contours
- Low potential for dressing induced sore spots²

Perforated wound contact layer

High rate of patient acceptance.²

- Perforated wound contact layer minimizes skin irritation (e.g., maceration, hematoma)²

Reliable exudate management



Improved wound healing.^{1,2}

- Fewer dressing changes lead to longer rest periods for wound healing
- Longer dressing change intervals result in lower costs



People.Health.Care.

Suprasorb® P and Suprasorb® P Border



- 1 Quick absorption of exudate
- 2 Good contact with wound bed
- 3 Bacteria and water repellent
- 4 Excess water vapor evaporates

<p>Exudate levels Wounds with moderate exudate levels</p>	<p>Wound depth For superficial wounds and for deep wounds when combined with an appropriate wound filler</p>	<p>Indications</p> <ul style="list-style-type: none"> ■ Venous or arterial leg ulcers ■ Pressure ulcers ■ Diabetic ulcers ■ First and second degree burns 	<p>Products</p> <p>Suprasorb® P</p> <p>Suprasorb® P Border</p> <p>Suprasorb® P Sacrum</p>
--	---	--	--

Suprasorb® P Foam dressings

Size	Item No.	HCPCS Code	Shipping Units (per box/case)
non-border			
5 x 5 cm (2 x 2 in.)	20405	A6209	10/100
7.5 x 7.5 cm (3 x 3 in.)	20406	A6209	10/100
10 x 10 cm (4 x 4 in.)	20407	A6209	10/100
15 x 15 cm (6 x 6 in.)	20408	A6210	5/50
15 x 20 cm (6 x 8 in.)	20409	A6210	5/50
border			
7.5 x 7.5 cm (3 x 3 in.)	20416	A6212	10/100
10 x 10 cm (4 x 4 in.)	20417	A6212	10/100
15 x 15 cm (6 x 6 in.)	20418	A6212	5/50
15 x 20 cm (6 x 8 in.)	20419	A6213	5/50
sacrum			
18 x 20.5 cm (7.2 x 8.1 in.)	20415	A6213	5/50

Sources

1. Bernhard Lange-Asschenfeldt, Annette Schätzle, Christiane Zell, Gabriele Stern, Ulrike Wagner, Reinhard Schaupp, Martin Abel: "Evaluation of ergonomic aspects in the daily routine of a new foam with a wound contact layer" (2013) EWMA 2013 Copenhagen; WundDach 2013 Friedrichshafen
2. Michael Schmitz: "Haptiktest Suprasorb P WCL" ASS GmbH, Bocholt, test report (2012)

