



People.Health.Care.

Suprasorb® P Heel **NEW**

Easy to apply, soft, flexible, and conformable.



Atraumatic dressing changes



High fluid absorption capacity



Specialty shape offers protection from friction, shear, and pressure



People.Health.Care.

Suprasorb® P Heel
Polyurethane foam dressing.

Three-layer material:

Your benefit:

Outer protective film

Reduces friction and shear. High Moisture Vapor Transmission Rate (MVTR) extends wear time.

Absorbent PU foam

Redistributes pressure and manages exudate.

Perforated non-adherent wound contact layer

Reduced pain and trauma when changing dressing.



Indications

- Superficial acute and chronic wounds, such as: pressure ulcers, diabetic foot ulcers, abrasions and surgical wounds
- Protects heels and elbows from friction, shear, and pressure
- For moderate to highly exuding wounds or to protect intact eschar



Suprasorb® P Heel Non-adhesive PU foam dressing
sterile, individually-sealed

Size	HCPCS Code	REF	Shipping Units (per box/case)
13 x 10 cm (5 x 4 in.)	A6010	180040	15/120

