



People.Health.Care.

Zorflex®

Removes odor, naturally.





People.Health.Care.

Zorflex® Activated Carbon Cloth Wound Dressing


Zorflex® activated carbon wound dressing is indicated as a contact layer for use in the management of acute and chronic, infected, and non-healing wounds.

Benefits

- Removes and sequesters odor-causing biomolecules
- Non-cytotoxic and safe to use
- Can be left on the wound for up to 7 days between changes
- Pain-free in application and use
- Easy to use with no special procedures or techniques required

Indications

- Acute and Chronic wounds including:
 - ◆ Venous leg ulcers
 - ◆ Diabetic foot ulcers
 - ◆ Pressure sores
 - ◆ Surgical wounds
 - ◆ Trauma wounds
 - ◆ Graft donor and recipient sites
 - ◆ Fungating carcinomas



- A 100% pure activated carbon cloth, derived from a natural source
- B Maintains direct wound contact to protect the wound bed
- C Removes and sequesters odor-causing biomolecules

Zorflex® activated carbon cloth wound contact dressing sterile, individually-sealed

Size	HCPCS Code	Item No.	Shipping Units (per box/case)
10 x 10 cm (4 x 4 in.)	A6206	140535	10/100
10 x 20 cm (4 x 8 in.)	A6207	140536	10/100
15 x 25 cm (6 x 10 in.)	A6208	140537	10/50



L&R USA INC.
3880 W Wheelhouse Road
Milwaukee, WI 53208 USA

Phone: 855-892-4140, 414-892-4140
Fax: 414-892-4150
Email: inquiries@us.LRmed.com



Chemviron Cloth Division, Rainton Bridge Industrial Estate, Houghton le Spring, Tyne and Wear, DH4 5PP, UK