



Wound Care Products

- **Wound Dressing Strips**
- **Precuts**
- **Wound Closure Strips**
- **Sterile Film Dressings and Adhesive Wound Dressings**
- **Absorbent Gauze, Gauze Dressings**
- **Gauze Balls**
- **Nonwoven Dressings**
- **Special Wound Dressings**



Curaplast® sensitive Wound Dressing Strips

Soft, supple nonwoven plaster for sensitive skin; first-aid dressing for hospital and doctor's surgery as well as for minor injuries at home and on the move.

The nonwoven material, which is slightly elastic widthwise, readily conforms to body contours and holds the dressing pad securely on the skin. A

micro-perforated separating film prevents adhesion in the wound area. The healing process is supported by dressing changes that are painless and gentle.



Curaplast® sensitive Wound Dressing Strips

In various sizes skin coloured nonwoven, 100% polyester, with continuous dressing pad, skin-friendly polyacrylate adhesive, free from colophonium and colophonium derivatives.

Curaplast® Sensitive Wound Dressing Strips			
for the practice and outpatient department	Size	REF	Shipping Units
5 m roll in cardboard dispenser	4 cm x 5 m	30 616	1 pc. / 36 pcs.
	6 cm x 5 m	30 617	1 pc. / 40 pcs.
	8 cm x 5 m	30 618	1 pc. / 24 pcs.
5 dressings, 10 cm each in cardboard box	6 cm x 0.5 m	17 076	1 pc. / 56 pcs.
10 dressings, 10 cm each	4 cm x 1 m	17 077	1 pc. / 47 pcs.
	6 cm x 1 m	17 078	1 pc. / 47 pcs.
	8 cm x 1 m	17 079	1 pc. / 36 pcs.
for household and travel			
Curaplast® Sensitive Finger Dressing			
for household and travel			
10 pcs. in cardboard box	2 x 18 cm	17 086	1 pack / 30 packs

Curaplast® **sensitive** Adhesive Wound Dressings



Curaplast® **sensitive** Injection Plasters

Skin coloured precut adhesive plasters off the roll (cardboard dispenser); quick and easy to handle for blood tests, injections, etc.

Curaplast® **sensitive** Precuts

Skin coloured, in different sizes, with central wound pad, each plaster sealed individually and hygienically.

Curaplast® Sensitive Injection Plaster			
in cardboard dispenser	Size	REF	Shipping Units
dispenser with 250 dressings	2 cm x 4 cm	30 625	1 pack / 16 packs
Curaplast® Sensitive Precuts			
10 pcs., assorted sizes in transparent pouch 6 dressings 4 dressings	1.9 cm x 7.2 cm 2.5 cm x 7.2 cm	17 082	1 pack / 340 packs
20 pcs., assorted sizes in cardboard box 10 dressings 4 dressings 4 dressings 2 dressings	1.9 cm x 7.2 cm 2.5 cm x 7.2 cm Ø 2.2 cm 3.8 cm x 3.8 cm	17 083	1 pack / 45 packs
500 pcs., individually sealed in cardboard box	2.5 x 7.2 cm	30 631	1 pack / 38 packs
Curaplast® Dressings, waterproof		Skin coloured, precut adhesive plasters, 2-layer composite of polyurethane film and polyester nonwoven, skin-friendly polyacrylate adhesive.	
10 pcs., assorted sizes in transparent pouch 8 dressings 2 dressings	2.5 cm x 7.2 cm 1.9 cm x 7.2 cm	17 084	1 pack / 30 packs
Curaplast® Kids		Adhesive plaster with children's motifs; water and dirt repellent, sealed individually and hygienically.	
15 dressings in cardboard box	1.7 cm x 6 cm	30 630	1 pack / 48 packs
Curaplast® Sensitive round		Elastic, skin coloured nonwoven, with skin-friendly polyacrylate adhesive, free from colophonium and colophonium derivatives.	
in cardboard dispenser roll of 100 pcs.	Ø 2.3 cm	30 626	1 pack / 24 packs



Curapor® Island Dressing, sterile

Ready-for-use island dressing for initial sterile wound care and for changing dressings on wounds of all types – in the operating theatre, outpatient department and doctor's surgery. This dressing with its extremely high absorbing capacity has a polyester/polyethylene layer on the wound side; this allows discharge to diffuse into the absorptive part, accelerates wound

cleansing, yet does not stick to the wound. This has a positive influence on the healing process by ensuring that the wound is not disturbed when the dressing is changed. Reliable all-round closure provided by particularly adherent skin-friendly polyacrylate adhesive; the Curapor wound dressing can be removed from the skin painlessly and without leaving any residue.



Curapor® Island Dressing, sterile

Elastic nonwoven carrier material, white, 100% polyester, with skin-friendly polyacrylate adhesive (free from colophonium and colophonium derivatives); dressing pad with an absorbent layer made from viscose/polyester and a non-adhesive polyester/poly-ethylene wound contact layer, individually sealed.

Curapor® Island Dressing, sterile			
Size	Wound Pad	REF	Shipping Units
7 x 5 cm	4 x 2.5 cm	30 504	5 pcs. / 100 pcs.
10 x 8 cm	6 x 4 cm	30 505	5 pcs. / 100 pcs.
7 x 5 cm	4 x 2.5 cm	22 120	100 pcs. / 600 pcs.
10 x 8 cm	6 x 4 cm	22 121	50 pcs. / 300 pcs.
10 x 15 cm	5 x 11 cm	22 122	50 pcs. / 300 pcs.
10 x 20 cm	5 x 16 cm	22 123	50 pcs. / 300 pcs.
10 x 25 cm	5 x 20 cm	22 124	50 pcs. / 250 pcs.
10 x 30 cm	5 x 25 cm	22 125	50 pcs. / 250 pcs.
10 x 34 cm	5 x 29 cm	22 126	50 pcs. / 250 pcs.

Curapor® Transparent Wound Dressing, sterile

Curapor transparent is a ready-for-use wound dressing for sterile care of postoperative wounds, accident injuries and small cuts and abrasions. For the protection against secondary infections. Allows patients to take a shower.


The transparent film holds safely and fits accurately like a second skin – even after longer wearing periods and while moving. Both the carrier film and the adhesive are skin-friendly – and tolerated well even by sensitive skin.

Curapor® Transparent Wound Dressing, sterile


Transparent film carrier, 100% polyurethane, with skin-friendly polyacrylate adhesive (free of colophonium or colophonium derivates). Dressing pad with an absorbent layer made from viscose/polyester and a non-adhesive polyester/polyethylene wound contact layer. Individually sealed, sterilised by EO.




Curapor transparent wound dressing is easy to apply:




Simply open the paper cover, ...



... pull apart, ...



... gently press onto the cleaned skin ...



... and pull off the support film – done!

Curapor® Transparent Wound Dressing, sterile

Size	Wound Pad	REF	Shipping Units
7 x 5 cm	4 x 2,5 cm	13 099	5 pcs. / 120 pcs.
10 x 8 cm	6 x 4 cm	13 100	5 pcs. / 170 pcs.
7 x 5 cm	4 x 2,5 cm	13 101	50 pcs. / 500 pcs.
10 x 8 cm	6 x 4 cm	13 102	25 pcs. / 200 pcs.
10 x 15 cm	5 x 11 cm	13 103	25 pcs. / 225 pcs.
10 x 20 cm	5 x 16 cm	13 104	25 pcs. / 225 pcs.
10 x 25 cm	5 x 20 cm	13 105	25 pcs. / 225 pcs.



Curapont® • Porofix® Wound Closure Strips

Porofix® Wound Closure Strips

Windowed first-aid dressing, purpose-designed for adapting the wound edges in superficial stab wounds and incised wounds, or in lacerations; also suitable for relieving pressure on sutures.

In the windowed area the plaster has an aluminized film backing on the wound side, which prevents the plaster from sticking to the wound. Discharge is able to drain away, and the wound can be checked visually.

Unstretchable fabric with a high tensile strength, 100% viscose, skin coloured, coated with a skin-compatible synthetic rubber adhesive.



Curapont® Wound Closure Strips

Plaster strips with high adhesive power, made from polyurethane reinforced with non-woven material, supple, tear-resistant and non-elastic; coated with hypoallergenic polyacrylate adhesive; for dependable, stitch-free wound closure; also for relieving pressure on sutures.

By adapting the wound edges cleanly and over a wide area, the resulting scar is perfect in cosmetic terms. The local anaesthetic

that is often required for a suture can be avoided; puncture channel infections, blood supply disorders in the suture area as well as the characteristic stitch marks are avoided.

Rapid and safe access can be gained at either end of the strip if necessary. This is particularly useful where half-length wound closure strips are sufficient; they are simply cut in half, complete with carrier paper, before being peeled off.



Porofix® Wound Closure Strip

2 pcs. in a sealed bag

REF

Shipping Units

2.5 x 12.5 cm

30 100

10 pcs. / 1560 pcs.

Curapont® Wound Closure Strips, sterile

sealed bags with

5 pcs.	3 x 75 mm	30 058	100 pcs. / 1000 pcs.
3 pcs.	6 x 75 mm	30 059	120 pcs. / 1200 pcs.
10 pcs.	6 x 100 mm	30 056	100 pcs. / 2000 pcs.
6 pcs.	12 x 100 mm	30 060	120 pcs. / 1200 pcs.
5 pcs.	12 x 140 mm	30 057	100 pcs. / 1000 pcs.

Suprasorb® F

Film Wound Dressing

Suprasorb® F Protects and holds (in)visibly!

Suprasorb F offers transparent help for protection and dressing retention. With the breathability of a semi-permeable polyurethane film.

Suprasorb F enables a moist wound climate conducive to wound healing by means of gas

and heat exchange, repels fluids and protects the wound against colonisation by germs from outside.


Indications: lightly exuding, superficial, non-infected wounds in the epithelialisation phase, e.g.: abrasions, dermabrasion, ulcers, decubitus, diabetic foot syndrome (DFS), second degree burns.



Similar to skin and extremely fine.

The Suprasorb F microporous film wound dressing offers a variety of "clear" benefits: visual wound inspection, very easy application and use, effective protection from external influences, allows showering and washing, flexible, elastic and tear-resistant.

Suprasorb F is available as a sterile film wound dressing and also as a non-sterile film dressing off the roll.



- 1 Transparent, semi-permeable polyurethane film.
- 2 Repels bacteria and fluids.
- 3 Surplus water vapour evaporates.

Interval between dressing changes

Depending on the degree of exudation the dressing should be changed when there is an accumulation of exudate over the wound size.

Suprasorb® F Film Wound Dressing, sterile, individually sealed

Adhesive Surface	REF	Shipping Units
5 x 7 cm	20 460	10 pcs. / 380 pcs.
5 x 7 cm	20 461	100 pcs. / 1000 pcs.
10 x 12 cm	20 462	10 pcs. / 100 pcs.
10 x 12 cm	20 463	50 pcs. / 200 pcs.
10 x 25 cm	20 464	10 pcs. / 140 pcs.
15 x 20 cm	20 465	10 pcs. / 400 pcs.
20 x 30 cm	20 466	10 pcs. / 180 pcs.

Suprasorb® F Film Wound Dressing, off the roll, non-sterile

10 cm x 1 m	20 467	1 pc. / 16 pcs.
10 cm x 10 m	20 468	1 pc. / 16 pcs.
15 cm x 10 m	20 469	1 pc. / 10 pcs.



Gazin® Gauze Dressings

Gazin® Absorbent Gauze

Gazin® Dressings

Sterile dressings made from absorbent gauze type 17, all cut edges inside.



Gazin® Gauze Dressings, sterile			
	Size	REF	Shipping Units
2 dressings per pouch, 8-ply	5 x 5 cm	13 600	10 pcs. / 600 pcs.
	7.5 x 7.5 cm	13 601	10 pcs. / 540 pcs.
	10 x 10 cm	13 602	10 pcs. / 560 pcs.
	10 x 20 cm	13 603	10 pcs. / 320 pcs.
2 dressings per pouch, 8-ply	5 x 5 cm	13 610	50 pcs. / 600 pcs.
	7.5 x 7.5 cm	13 611	100 pcs. / 1800 pcs.
	10 x 10 cm	13 612	50 pcs. / 900 pcs.
	10 x 20 cm	13 613	50 pcs. / 500 pcs.
2 dressings per pouch, 8-ply	5 x 5 cm	13 620	100 pcs. / 1800 pcs.
	7.5 x 7.5 cm	13 621	100 pcs. / 1800 pcs.
	10 x 10 cm	13 622	100 pcs. / 2000 pcs.
	10 x 20 cm	13 623	100 pcs. / 800 pcs.
2 dressings per pouch, 12-ply	5 x 5 cm	13 690	80 pcs. / 1440 pcs.
	7.5 x 7.5 cm	13 691	80 pcs. / 1440 pcs.
	10 x 10 cm	13 692	80 pcs. / 1600 pcs.
	10 x 20 cm	13 694	80 pcs. / 640 pcs.
5 dressings per pouch, 12-ply	5 x 5 cm	13 695	100 pcs. / 2100 pcs.
	7.5 x 7.5 cm	13 696	100 pcs. / 800 pcs.
	10 x 10 cm	13 697	100 pcs. / 600 pcs.
	10 x 20 cm	13 699	100 pcs. / 400 pcs.
5 dressings per pouch, 8-ply	7.5 x 7.5 cm	11 624	100 pcs. / 1200 pcs.
	10 x 10 cm	11 586	100 pcs. / 1000 pcs.
10 dressings per pouch, 8-ply	7.5 x 7.5 cm	13 646	100 pcs. / 1200 pcs.
	10 x 10 cm	11 583	300 pcs. / 1500 pcs.
10 dressings per pouch, 12-ply	7.5 x 7.5 cm	23 617	100 pcs. / 1200 pcs.
	10 x 10 cm	11 574	100 pcs. / 1000 pcs.
10 dressings per pouch, 16-ply	10 x 10 cm	13 683	100 pcs. / 1000 pcs.
Dressings with Y-cut	7.5 x 7.5 cm	13 521	100 pcs. / 1800 pcs.
	10 x 10 cm	13 522	100 pcs. / 2000 pcs.
2 dressings per pouch			

Gazin® Gauze Dressings		Absorbent gauze, all cut edges inside, no fraying threads; in paper bag, sterilisable (A 121°C).			
		Size	REF	Shipping Units	
gauze type 17	8-ply	5 x 5 cm	18 500	100 pcs. / 1000 pcs.	
		7.5 x 7.5 cm	18 503	100 pcs. / 1000 pcs.	
		10 x 10 cm	18 506	100 pcs. / 1000 pcs.	
	12ply	5 x 5 cm	18 501	100 pcs. / 1000 pcs.	
		7.5 x 7.5 cm	18 504	100 pcs. / 1000 pcs.	
		10 x 10 cm	18 507	100 pcs. / 1000 pcs.	
		10 x 12.5 cm	18 509	100 pcs. / 1000 pcs.	
		10 x 20 cm	18 510	100 pcs. / 1000 pcs.	
	16-ply	5 x 5 cm	18 502	100 pcs. / 1000 pcs.	
		7.5 x 7.5 cm	18 505	100 pcs. / 1000 pcs.	
		10 x 10 cm	18 508	100 pcs. / 1000 pcs.	
		10 x 12.5 cm	10 956	100 pcs. / 800 pcs.	
gauze type 20	32-ply	10 x 10 cm	10 942	50 pcs. / 500 pcs.	
X-ray detectable					
gauze type 17	8-ply	10 x 10 cm	10 947	100 pcs. / 1000 pcs.	
	12-ply	7.5 x 7.5 cm	18 511	100 pcs. / 1000 pcs.	
		10 x 10 cm	18 514	100 pcs. / 1000 pcs.	
		10 x 12.5 cm	18 516	100 pcs. / 1000 pcs.	
		10 x 20 cm	18 517	100 pcs. / 1000 pcs.	
	16-ply	5 x 5 cm	10 972	100 pcs. / 2000 pcs.	
		7.5 x 7.5 cm	18 512	100 pcs. / 1000 pcs.	
		10 x 10 cm	18 515	100 pcs. / 1000 pcs.	
		10 x 12.5 cm	10 957	100 pcs. / 800 pcs.	
	gauze type 20	32-ply	10 x 10 cm	10 943	50 pcs. / 500 pcs.
	gauze type 17	32-ply	10 x 10 cm	18 603	50 pcs. / 500 pcs.
Gazin® Absorbent Gauze		Absorbent gauze type 20, width: 80 cm, 8-ply folded to 10 cm width, rolled			
in plastic dispenser	8-ply	10 cm x 1 m	12 002	1 pack / 80 packs	
		10 cm x 2 m	12 003	1 pack / 60 packs	
		10 cm x 5 m	12 004	1 pack / 36 packs	
		10 cm x 10 m	12 005	1 pack / 24 packs	
in a pouch refill (paper)		2 x 20 m = 40 m	12 007	1 pack / 6 packs	
		10 m	12 006	1 pack / 22 packs	
Gazin® Absorbent Gauze		5 m on a roll, in special pouch and cardboard box, ready for sterilisation			
	4-ply	6 cm x 5 m	12 010	1 pc. / 30 pcs.	
		8 cm x 5 m	12 011	1 pc. / 30 pcs.	
		10 cm x 5 m	12 012	1 pc. / 43 pcs.	
Dispenser		Polystyrol, transparent, suitable for rolled dressing materials, in surgeries and outpatient departments, e.g. Gazin, Zemuko, Zelletten, etc. max. width of contents 10 cm			
empty max. width of contents		10 cm	23 900	1 pc.	

Gazin® Gauze Balls

Gazin® Gauze Balls

Made from DIN EN 14079 absorbent gauze 20 threads; optionally with or without X-ray contrast thread; for absorbing blood and discharge, for wound cleansing, also for disinfection of the skin.

Gazin® Dissecting Swabs

Made from DIN EN 14079 absorbent gauze 24 threads; with X-ray contrast thread; firm construction for preparing layers of tissue and anatomical structures for surgical interventions.



Gazin® Gauze Balls	non-sterile without X-ray detectable thread in a pouch	
	REF	Shipping Units
walnut size	15 175	1000 pcs. / 5000 pcs.
plum size	15 176	1000 pcs. / 5000 pcs.
egg size	15 177	500 pcs. / 2500 pcs.
extra large	15 179	500 pcs. / 2000 pcs.
fist size	15 180	250 pcs. / 1000 pcs.
	non-sterile with X-ray detectable thread in a pouch	
walnut size	15 183	1000 pcs. / 5000 pcs.
plum size	15 184	1000 pcs. / 5000 pcs.
egg size	15 185	500 pcs. / 2500 pcs.
extra large	15 186	500 pcs. / 2000 pcs.
fist size	15 223	250 pcs. / 1000 pcs.
Gazin® Gauze Balls	sterile in blister pack, packs with 2+3 or 2+2 pcs., without X-ray detectable thread	
plum size, 2+2 pcs.	14 961	100 pcs. (25x 4) / 1000 pcs. (250x 4)
plum size, 2+3 pcs.	14 962	125 pcs. (25x 5) / 1250 pcs. (250x 5)
egg size, 2+3 pcs.	14 968	100 pcs. (20x 5) / 1000 pcs. (200x 5)
extra large, 2+3 pcs.	14 969	100 pcs. (20x 5) / 1000 pcs. (200x 5)
dialysis, 2+2 pcs.	14 965	100 pcs. (25x 4) / 1000 pcs. (250x 4)
dialysis, 2+3 pcs.	14 966	125 pcs. (25x 5) / 1250 pcs. (250x 5)
	with X-ray detectable thread	
plum size, 2+3 pcs.	14 967	125 pcs. (25x 5) / 1250 pcs. (250x 5)
	sterile in blister pack, without X-ray detectable thread	
plum size, 3 pcs.	15 501	90 pcs. (30x 3) / 720 pcs. (240x 3)
plum size, 10 pcs.	15 500	300 pcs. (30x10) / 1500 pcs. (150x10)
plum size, 20 pcs.	15 499	400 pcs. (20x20) / 2000 pcs. (100x20)
egg size, 10 pcs.	15 498	300 pcs. (30x10) / 1500 pcs. (150x10)
Gazin® Dissecting Swabs	non-sterile with X-ray detectable thread in a pouch	
small	15 155	1000 pcs. / 3000 / 45000 pcs.
medium	15 156	1000 pcs. / 2000 / 30000 pcs.
large	15 157	500 pcs. / 1000 / 15000 pcs.

Soft, economical universal dressing, made from a combination of a fine nonwoven layer and bleached cellulose; for wounds with heavy discharge, as an additional ab-

sorbent layer, for ointment dressings and moist dressings; also for padding.

The wound-friendly nonwoven material and the absorbent lay-

er are joined together without any creases. Zemuko dressings can thus be cut off the roll with well-defined edges and no fluff.

Zemuko®

Cellulose/nonwoven combination. Wound side: nonwoven, 67% viscose, 33% polyester; absorbent layer: bleached cellulose (in accordance with German Pharmacopoeia DAB); the two layers are joined together by embossing; sterilised, or sterilisable (A 121°C).

Zemuko® seamless

Alternative universal dressing; made from bleached cellulose DAB with gauze **covering on all sides**.



Zemuko®			
	Size	REF	Shipping Units
rolled, in plastic dispenser	10 cm x 1 m	23 512	1 pc. / 60 pcs.
	10 cm x 2 m	23 515	1 pc. / 30 pcs.
rolled, in a pouch for the practice and outpatient department	10 cm x 10 m	23 502	1 pc. / 20 pcs.
	15 cm x 10 m	23 503	1 pc. / 16 pcs.
	28 cm x 10 m	23 602	1 pc. / 10 pcs.
Zemuko® seamless			
rolled, in a pouch	10 cm x 10 m	23 702	1 pc. / 20 pcs.
	15 cm x 10 m	23 703	1 pc. / 16 pcs.
Zemuko® Wound Dressings sterile, individually sealed			
	10 x 10 cm	23 802	10 pcs. / 200 pcs.
Dispenser		Polystyrol, transparent, suitable for rolled dressing materials, in surgeries and outpatient departments, e.g. Gazin, Zemuko, Zelletten, etc. max. width of contents 10 cm	
empty, max. width of contents: 10 cm		23 900	1 pc.

Solvaline® N

Solvaline N dressings consist of highly absorbent cotton wool, covered on both sides by a finely perforated polyester film which is interconnected with the cotton by heat-treatment.

Wound exudates are quickly absorbed through the perforations,

and the dressing does not stick to the wound. This accelerates the healing process and enables an atraumatic dressing change.

Because there is a film cover on both sides of the dressing, either side of Solvaline N may be applied to the wound. This makes

mistakes of application virtually impossible and Solvaline N really safe – to the satisfaction of patients and carers.

Solvaline® N

– ideally suited for

- superficial wounds in post-operative treatment
- first aid, e.g. abrasions and cuts
- 1st and 2nd degree burns



Solvaline® N			
	Size	REF	Shipping Units
sterile, individually sealed	5 x 5 cm	23 230	100 pcs. / 7.200 pcs.
	10 x 10 cm	23 231	100 pcs. / 2.400 pcs.
	10 x 20 cm	23 232	50 pcs. / 1.000 pcs.
	20 x 30 cm	23 233	25 pcs. / 300 pcs.
	5 x 5 cm	23 234	25 pcs. / 900 pcs.
	10 x 10 cm	23 235	10 pcs. / 300 pcs.
Solvaline® N			
non-sterile, in cardboard box	5 x 5 cm	23 220	200 pcs. / 20.000 pcs.
	10 x 10 cm	23 221	150 pcs. / 3.600 pcs.
	10 x 20 cm	23 222	50 pcs. / 1.200 pcs.
	20 x 30 cm	23 223	25 pcs. / 300 pcs.



Vliwasoft®

Soft and supple universal dressing friendly with low adhesion tendency. Vliwasoft is economical for use in the field of wound cleansing and wound care.



The Vliwasoft drainage sponge is used for covering drains and cannulae in dialysis, traction following fractures, etc.

Vliwasoft®			
	Size	REF	Shipping Units
sterile, 2 dressings sealed, 75 x 2 dressings per pack	5 x 5 cm	12 076	150 pcs. / 3000 pcs.
	7.5 x 7.5 cm	12 082	150 pcs. / 3000 pcs.
	10 x 10 cm	12 065	150 pcs. / 3000 pcs.
	10 x 20 cm	12 071	150 pcs. / 1500 pcs.
4-ply	5 x 5 cm	12 079	100 pcs. / 2000 pcs.
	7.5 x 7.5 cm	12 085	100 pcs. / 2000 pcs.
	10 x 10 cm	12 068	100 pcs. / 2000 pcs.
	10 x 20 cm	12 073	100 pcs. / 1000 pcs.
non-sterile, paper bag with 100 dressings	5 x 5 cm	12 075	100 pcs. / 5000 pcs.
	7.5 x 7.5 cm	12 081	100 pcs. / 2500 pcs.
	10 x 10 cm	12 064	100 pcs. / 2000 pcs.
	10 x 20 cm	12 070	100 pcs. / 1000 pcs.
4-ply	5 x 5 cm	12 078	100 pcs. / 5000 pcs.
	7.5 x 7.5 cm	12 084	100 pcs. / 2000 pcs.
	10 x 10 cm	12 067	100 pcs. / 2000 pcs.
	10 x 20 cm	12 072	100 pcs. / 1000 pcs.
6-ply			
Vliwasoft® Drainage Sponges with Y-cut			
sterile, 2 dressings sealed, 50 x 2 dressing per pack	7.5 x 7.5 cm	12 099	100 pcs. / 2000 pcs.
	10 x 10 cm	12 098	100 pcs. / 2000 pcs.

Vliwasorb®

Superabsorbent Wound Dressing

Wound dressing with an absorbent core consisting of superabsorbent polymer particles, for the management of heavily exuding wounds, e.g. pressure ulcers, leg ulcers, etc.

The Vliwasorb wound dressing has a particularly high absorption capacity, as superabsorbent polymer particles are capable of absorbing several times their own weight in liquids (e.g. 20 times as much blood). The exudate is safely retained even when mechanical pressure is exerted. Vliwasorb is therefore perfectly safe to use under compression dressings.

Vliwasorb allows an economical management of heavily exuding wounds, as fewer dressing changes are required due to Vliwasorb's high absorption capacity. This allows longer peri-

ods of wound rest and reduces the risk of pathogen transmission. The absorption and retention of exudate helps to prevent wound edge maceration.

Another major benefit of Vliwasorb is its supreme wearing comfort. The elastic non-adherent wound contact layer is gentle on the skin, protecting it from

damage during dressing change. A blue nonwoven clothing protection layer keeps wound exudates from seeping through and makes the correct dressing application easier and more convenient. With its slim design and high flexibility, Vliwasorb provides superior wearing comfort even underneath compression dressing and in areas that are difficult to dress.



Vliwasorb®		
Size	REF	Shipping Units
10 x 10 cm	24 501	10 pcs. / 360 pcs.
10 x 20 cm	24 502	10 pcs. / 340 pcs.
20 x 20 cm	24 503	10 pcs. / 300 pcs.



Vliwazell®

Absorbent Wound Dressing

Fully-enclosed, supple dressing with very high absorption capacity; for the care of wounds with heavy discharge. Vliwazell dressings are wound-friendly

through their particularly soft covering of nonwoven material. Wound exudate is quickly channelled into the absorptive cellulose layer, whose capacity is

many times greater than that of a gauze dressing. A moisture inhibiting layer of nonwoven material on the back protects the clothing.



Vliwazell®	Size	REF	Shipping Units
non-sterile, 25 pcs. in plastic bag	10 x 10 cm	30 430	25 pcs. / 900 pcs.
	10 x 20 cm	30 431	25 pcs. / 450 pcs.
	15 x 20 cm	30 432	25 pcs. / 200 pcs.
	15 x 25 cm	30 433	25 pcs. / 200 pcs.
	20 x 20 cm	30 434	25 pcs. / 200 pcs.
	20 x 25 cm	30 435	25 pcs. / 150 pcs.
	20 x 40 cm	30 436	25 pcs. / 200 pcs.
sterile, individually sealed	10 x 10 cm	30 450	60 pcs. / 360 pcs.
	10 x 20 cm	30 451	30 pcs. / 300 pcs.
	15 x 20 cm	30 452	30 pcs. / 180 pcs.
	15 x 25 cm	30 453	30 pcs. / 240 pcs.
	20 x 20 cm	30 454	30 pcs. / 120 pcs.
	20 x 25 cm	30 455	30 pcs. / 150 pcs.
	20 x 40 cm	30 456	30 pcs. / 180 pcs.

Highly absorbent, soft, supple universal dressing for the care of wounds with heavy discharge. The polypropylene tissue on the wound side allows discharge to pass rapidly to the absorbent core without becoming adherent, so that dressings can be changed without damaging tissue.

Vliwin®

Absorbent core:

A layer of nonwoven material transfers discharge quickly into the absorbent core (cellulose fluff) and allows large quantities of exudate to be absorbed and distributed evenly.

Covering:

2 moisture repellent, breathable layers of nonwoven material inhibit the outward passage of discharge, and clothing is protected against soiling. At the same time, it is easy to recognise when the dressing needs to be changed.



Vliwin®			
	Size	REF	Shipping Units
non-sterile, 100 pcs.	10 x 10 cm	12 815	100 pcs. / 900 pcs.
	10 x 20 cm	12 816	100 pcs. / 600 pcs.
	15 x 25 cm	12 818	100 pcs. / 200 pcs.
	20 x 20 cm	12 817	100 pcs. / 300 pcs.
50 pcs.	20 x 40 cm	12 819	50 pcs. / 150 pcs.
sterile, individually sealed	10 x 10 cm	12 810	25 pcs. / 400 pcs.
	10 x 20 cm	12 811	25 pcs. / 250 pcs.
	15 x 25 cm	12 813	15 pcs. / 135 pcs.
	20 x 20 cm	12 812	25 pcs. / 200 pcs.
	20 x 40 cm	12 814	12 pcs. / 96 pcs.

Vliwaktiv®

Vliwaktiv® Ag

Vliwaktiv® Activated Charcoal Dressing

Ideal combination of a highly absorbent pad with a wound dressing that absorbs both odour and bacteria, for the treatment of malodorous wounds such as leg ulcers, decubital ulcers, anal fistulae, etc.

A skin-friendly, breathable non-woven cover speeds up the passage of exudate into the absorbent core. Bacteria and odours are effectively bound in Vliwaktiv. The clothing protection prevents the seeping of wound exudate.

Vliwaktiv® Ag Activated Charcoal Dressing with Silver

Activated charcoal dressing impregnated with silver, in a high-quality nonwoven cover that does not stick to the wound; for the treatment of malodorous infected and critically colonised wounds, such as leg ulcers, decubital ulcers, exulcerative tumors, etc. Presentations: rope and absorbent dressings.

Vliwaktiv Ag combines the antimicrobial action of silver with the benefits of activated charcoal which adsorbs odours and binds bacteria.

When Vliwaktiv Ag makes contact with exudate, it releases silver ions and destroys the retained microorganisms. Thus, microbial counts are effectively reduced without undue stress to the wound.

Unpleasant wound odours which are often associated with infected wounds are effectively adsorbed by the silver/activated

charcoal fabric, which enhances the well-being of patients and carers alike. Vliwaktiv Ag also effectively binds endotoxins which may impair the wound healing process.

Deep wounds and wound cavities can be treated easily with the flexible, soft Vliwaktiv activated charcoal rope. For heavily exuding superficial wounds, the Vliwaktiv Ag activated charcoal dressing is made with a highly absorbent core and integrated clothing protection which prevents the seeping of exudate.



Vliwaktiv® sterile			
	Size	REF	Shipping Units
individually sealed	10 x 10 cm	20 254	20 pcs. / 400 pcs.
	10 x 20 cm	20 255	20 pcs. / 200 pcs.
Vliwaktiv® Ag Activated Charcoal Rope with Silver, sterile			
individually sealed	6.5 x 10 cm	20 700	10 pcs. / 300 pcs.
	10 x 10 cm	20 701	10 pcs. / 300 pcs.
	10 x 15 cm	20 702	10 pcs. / 340 pcs.
Vliwaktiv® Ag Absorbent Activated Charcoal Dressing with Silver, sterile			
individually sealed	10 x 10 cm	20 705	10 pcs. / 150 pcs.
	10 x 20 cm	20 706	10 pcs. / 110 pcs.
	20 x 20 cm	20 707	10 pcs. / 130 pcs.

Metalline® Wound Dressings

The aluminized wound pad of the Metalline dressings is a supple and soft, extremely fine non-woven material with good drainage. The smooth surface does not stick to the wound. This ensures that the wound is not disturbed, and the dressing can be changed without trauma, thus promoting the healing process. **Metalline Drainage Dressing:** the central slit and those radiating from the punch hole ensure a good fit around small calibre drains.



Metalline®

Multi-layered special dressing material.

Wound side: fine nonwoven veil, with vacuum-aluminized fibres.

Absorbent layer: high-bulk non-woven material.

Rear side: needle-punched non-woven material. Layers joined through needle punching and impervious dispersion adhesive, resulting in good absorptive characteristics.



Metalline® Dressings

For all types of wounds – in traumatology, surgery, dermatology and phlebology, especially for extensive wounds, abrasions, burns; also for first aid use.

Metalline® Drainage Dressings

For care of wounds when small calibre drains or infusion catheters are involved, as well as for wire or nail extensions.

Hole diameter:
from 4.5 mm - 11 mm
for Charrière 9 - 30

Metalline® Tracheo Dressings

For care with tracheostomies or large-diameter drains.

Hole diameter:
from 12 mm - 19 mm
for Charrière 31 - 50

Metalline® Wound Dressings			
	Size	REF	Shipping Units
sterile, individually sealed	8 x 10 cm	23 083	10 pcs. / 560 pcs.
	10 x 12 cm	23 084	10 pcs. / 300 pcs.
sterile, individually sealed; for OR, outpatient department and practice, in dispenser box	8 x 10 cm	23 081	50 pcs. / 600 pcs.
	10 x 12 cm	23 082	50 pcs. / 300 pcs.
	10 x 20 cm	23 087	50 pcs. / 300 pcs.
non-sterile, rolled, in foil pouch	5 m x 10 cm	23 080	1 pc. / 20 pcs.
Metalline® Drainage Dressing			
sterile, individually sealed, in dispenser box	6 x 7 cm	23 093	50 pcs. / 600 pcs.
Metalline® Tracheo Dressing			
sterile, individually sealed, in dispenser box	8 x 9 cm	23 094	50 pcs. / 600 pcs.



Metalline® Non-adherent Dressing Kits and Dressing Cloths

Aluminized nonwoven material, sometimes combined with absorbent material, for the treatment of large area wounds and open wounds. The aluminized surface of Metalline dressing does not stick to the wound.

Metalline® Non-adherent Dressing Kit

For first aid and emergency care.

Combination of Metalline dressing 35 x 45 cm and gauze bandage 6 cm x 3 m.

Metalline® Non-adherent Bed Sheet

Aluminized nonwoven combined with absorbent material, for use as an underlay in the treatment of large-scale burns, chemical burns and abrasions.

Aluminized nonwoven material

combined with absorbent nonwoven material.

Metalline® Non-adherent Cloth

Aluminized nonwoven material off the roll, for use as an underlay or for covering large-scale wounds, for the cautious bedding of patients.

Aluminized nonwoven material off the roll. Metalline non-adherent cloth and any other absorbent layers must be sterilised before use.

Metalline® Non-adherent Dressing Cloth

For the care of large-scale injuries, burns, chemical burns and abrasions, in first aid and for emergency use.

Aluminized nonwoven combined with absorbent nonwoven material. The cloths are provided with ties to make them easier to handle during positioning and securing.



Metalline® Non-adherent Dressing Kit		
	REF	Shipping Units
sterile , in foil pouch		
consisting of: 1 dressing 35 cm x 45 cm 1 gauze bandage 6 cm x 3 m	23 000	1 pc. / 150 pcs.
Metalline® Non-adherent Bed Sheet		
sterile , in foil pouch		
2.5 m x 73 cm	23 050	1 pc. / 20 pcs.
Metalline® Non-adherent Cloth		
non-sterile , rolled, in cardboard box		
10 m x 73 cm	23 070	1 pc. / 10 pcs.
Metalline® Non-adherent Dressing Cloth		
sterile , in foil pouch		
40 x 60 cm	23 012	1 pc. / 100 pcs.
60 x 80 cm	23 010	1 pc. / 100 pcs.
80 x 120 cm	23 011	1 pc. / 100 pcs.