



Support and Relief

- Dressings for Compression Therapy

Porotape® Haftan®

Porotape®

Non-elastic adhesive tape with high tensile strength, which can be torn both lengthwise and widthwise, with an effective zinc oxide rubber adhesive.

Porotape can be applied to the skin directly, but is generally used in combination with adhesive bandages, for partial immobilization with a support and

relieving function on joints, muscles and ligaments; also as a prophylactic dressing aimed at avoiding overloading.

The adhesive tape is dispensed easily from the reel; the serrated edge makes tearing easy and prevents frayed edges.

Fine fabric, 100% viscose, zinc oxide rubber adhesive, raw white.



Haftan®

Extremely thin bandage (0.8 mm) made from polyurethane foam, open-pored, supple. Used as skin protection under adhesive bandages and tape dressings, as well as for securing pads which provide extra compression. The material's surface makes bandage turns adhere well to each other.

Haftan is also used for the smooth application of padding material underneath P.o.P. and synthetic casts.

Polyurethane foam, white, resistant to ointments, sweat and cleaning agents.

Porotape®		REF	Shipping Units
length: 10 m	2.5 cm	30 081	6 pcs. / 120 pcs.
	3.8 cm	30 084	1 pc. / 96 pcs.
	3.8 cm	30 082	12 pcs. / 96 pcs.
	5.0 cm	30 083	6 pcs. / 72 pcs.
Haftan®			
thickness: 0.8 mm	7 cm x 27.5 m	22 315	6 pcs. / 60 pcs.
	10 cm x 27.5 m	22 316	6 pcs. / 60 pcs.

Porodress®

Adhesive Bandage

Widthwise elastic adhesive bandage, lengthwise inelastic, with strong adhesive power.

Porodress adhesive bandages are used chiefly as a stirrup for lesions on capsules or ligaments and for muscle and tendon injuries or disorders; also for redressement and traction,

especially in combination with Porelast applied in circular turns.

The fabric with porously applied adhesive is permeable to air and water vapour, thus offering good ventilation even with overlapping bandage turns.



Porodress®

Elastic textile fabric, with cut edges, skin coloured, 58% cotton, 42% viscose, widthwise extensibility approx. 60%, strong synthetic rubber adhesive, generally well tolerated by the skin, X-ray-permeable.

Porodress®			
individually in cardboard box	Width	REF	Shipping Units
2.5 m (lengthwise inelastic)	6 cm	30 240	1 pc. / 50 pcs.
	8 cm	30 241	1 pc. / 50 pcs.
	10 cm	30 242	1 pc. / 50 pcs.



Panelast®

Panelast® PrO₂ Adhesive Bandages

Lengthwise and widthwise elastic adhesive bandages, offering a choice of two types of adhesive, for compression of the extremities (permanent bandages) in the field of phlebology; as a foundation layer or anchorage for tape applications in traumatology and sports medicine.

For support and relief with distortions, contusions and dislocations; also used for aftercare with fractures. The fabric with porously applied adhesive is permeable to air and water vapour, thus offering good ventilation even with overlapping bandage turns.

Panelast®

Lengthwise and widthwise elastic adhesive bandage with selvages, skin coloured, 100% cotton, lengthwise extensibility approx. 80%, widthwise approx. 45%, strong synthetic rubber adhesive, generally well tolerated by the skin, X-ray-permeable.



Panelast® PrO₂

Lengthwise and widthwise elastic adhesive bandage with cut edges, skin coloured, 100% cotton, lengthwise extensibility approx. 70%, widthwise approx. 45%, hypoallergenic adhesive, with high air permeability, X-ray-permeable.



Panelast®			
individually in cardboard box	Width	REF	Shipping Units
2.5 m (unstretched length)	6 cm	30 270	1 pc. / 50 pcs.
	8 cm	30 271	1 pc. / 50 pcs.
	10 cm	30 272	1 pc. / 50 pcs.
Panelast® PrO₂			
individually in cardboard box	Width	REF	Shipping Units
2.5 m (unstretched length)	6 cm	33 690	1 pc. / 50 pcs.
	8 cm	33 691	1 pc. / 50 pcs.
	10 cm	33 692	1 pc. / 50 pcs.

Porelast® • Porelast Superforte

Porelast® PrO₂ Adhesive Bandages

Lengthwise elastic adhesive bandages, with a choice of two types of adhesive, for permanent compression bandages in the field of phlebology and for functional bandages in sports medicine; for support and relief with distortions, contusions and

dislocations; also for redressment bandages, e.g. for malposition of the foot.

The fabric with porously applied adhesive is permeable to air and water vapour, thus offering good ventilation even with overlapping bandage turns.



Porelast® Porelast® Superforte

Elastic textile fabric, with selvedges, 100% cotton, skin coloured, lengthwise extensibility approx. 80%, strong synthetic rubber adhesive, generally well tolerated by the skin, X-ray-permeable.

The difference between the two bandages is the higher adhesive power of the Superforte version.

Porelast® PrO₂

Elastic textile fabric with cut edges, 100% cotton, skin coloured, lengthwise extensibility approx. 70%, hypoallergenic adhesive, with high air-permeability, X-ray-permeable.

Porelast®	Width	REF	Shipping Units
2.5 m (unstretched length) individually in cardboard box	6 cm	30 210	1 pc. / 50 pcs.
	8 cm	30 211	1 pc. / 50 pcs.
	10 cm	30 212	1 pc. / 50 pcs.
2.5 m (unstretched length) unwrapped	6 cm	30 213	10 pcs. / 50 pcs.
	8 cm	30 214	10 pcs. / 50 pcs.
	10 cm	30 215	10 pcs. / 50 pcs.
Porelast® Superforte			
2.5 m (unstretched length) individually in cardboard box	6 cm	95 410	1 pc. / 50 pcs.
	8 cm	95 411	1 pc. / 50 pcs.
	10 cm	95 412	1 pc. / 50 pcs.
Porelast® PrO₂			
2.5 m (unstretched length) individually in cardboard box	6 cm	30 220	1 pc. / 50 pcs.
	8 cm	30 221	1 pc. / 50 pcs.
	10 cm	30 222	1 pc. / 50 pcs.
2.5 m (unstretched length) unwrapped	6 cm	30 223	10 pcs. / 50 pcs.
	8 cm	30 224	10 pcs. / 50 pcs.
	10 cm	30 225	10 pcs. / 50 pcs.



Fortelast® • Fortelast® Plus

Fixona® Elastic Adhesive Bandages

Lengthwise elastic bandages, made of skin-friendly fabric, with excellent extensibility as a result of the porously applied

synthetic adhesive. Both bandages are permeable to air and moisture. They are indicated for compression and support bandaging in traumatology, functional bandaging in sports medicine, and in phlebology for compression.

bandaging in traumatology, functional bandaging in sports medicine, and in phlebology for compression.



Fortelast®

White, 85% cotton, 15% viscose, lengthwise elastic, extensibility approx. 70%, hypoallergenic porously applied synthetic adhesive.

Fixona®

White, 100% cotton, lengthwise elastic, extensibility approx. 80%, synthetic adhesive with increased adhesive power.

Fortelast®	Width	REF	Shipping Units
2.5 m (unstretched length) individually in cardboard box	3 cm	94 740	1 pc. / 75 pcs.
	6 cm	94 741	1 pc. / 90 pcs.
	8 cm	94 742	1 pc. / 90 pcs.
	10 cm	94 743	1 pc. / 72 pcs.
	15 cm	94 745	1 pc. / 24 pcs.
Fortelast® Plus			
rolled inwards	2.5 cm	16 752	1 pc. / 84 pcs.
	5 cm	16 753	1 pc. / 60 pcs.
	7.5 cm	16 754	1 pc. / 60 pcs.
	10 cm	16 755	1 pc. / 72 pcs.
	15 cm	16 756	1 pc. / 24 pcs.
rolled outwards	2.5 cm	16 757	1 pc. / 84 pcs.
	5 cm	16 758	1 pc. / 60 pcs.
	7.5 cm	16 759	1 pc. / 60 pcs.
	10 cm	16 760	1 pc. / 72 pcs.
Fixona®			
4.5 m (stretched length) individually in cardboard box adhesive coating on the outside	5 cm	94 751	1 pc. / 50 pcs.
	7.5 cm	94 752	1 pc. / 50 pcs.
	10 cm	94 753	1 pc. / 50 pcs.
	15 cm	94 755	1 pc. / 50 pcs.

Rosidal® haft short stretch

Raucodur Cohesive

Cohesive short stretch bandages with strong compression; for compression of the extremities with disorders of the veins and lymphatic system, as permanent bandages, for the pre- and postoperative thromboprophylaxis, and for compression of the thigh to complement the lower leg bandage. The compression remains therapeutically effective over a prolonged period.

These bandages are also suitable for use as functional bandages, for immobilising, supporting and relieving capsules, ligaments, muscles and joints.

The stretch characteristics are similar to those of Rosidal K (see page 89), and – despite their adhesive properties – the bandages are permeable to air and water vapour.



Skin coloured fabric, 100% cotton, elastic textile with twisted warp threads, low stretch, cohesive through microfine latex dots; can be washed at 40°C.

Rosidal haft: Extensibility approx. 60%.

Raucodur Cohesive: Extensibility approx. 75%.



Rosidal® haft			
individually in cardboard box	Width	REF	Shipping Units
5 m (stretched length)	6 cm	22 270	1 pc. / 60 pcs.
	8 cm	22 271	1 pc. / 40 pcs.
	10 cm	22 272	1 pc. / 40 pcs.
	12 cm	22 273	1 pc. / 30 pcs.
Raucodur Cohesive			
individually in cardboard box			
5 m (stretched length)	8 cm	20 173	1 pc. / 50 pcs.
	10 cm	20 171	1 pc. / 50 pcs.
	12 cm	20 172	1 pc. / 48 pcs.

Perfekta® Cohesive

long stretch


Long stretch bandage with durable elasticity, cohesive on both sides; a support and relief bandage for the prophylactic relief of capsules and ligaments, prevention of haematomas, and for contusions, distortions, dislocations and tendon injuries; also used as a compressive bandage.

Its cohesivity makes it especially suitable for conically shaped or frequently moved parts of the body.



Perfekta® Cohesive

Long stretch bandage with durable elasticity, skin coloured, 96% cotton, 4% elastane, lengthwise extensibility approx. 140%, cohesive on both sides due to microfine latex dots, washable at 40°C.



Perfekta® Cohesive

individually in cardboard box	Width	REF	Shipping Units
5 m (stretched length)	6 cm	20 106	1 pc. / 50 pcs.
	8 cm	20 107	1 pc. / 50 pcs.
	10 cm	20 104	1 pc. / 50 pcs.
	12 cm	20 105	1 pc. / 50 pcs.

Mollelast® haft

Mollelast® haft color

Cohesive conforming bandage with special stretch characteristics.

Mollelast haft is lengthwise elastic and widthwise extensible. This makes the bandage perfectly suitable for application on cone-shaped parts of the body, and offers a comfortable

and secure fit even in joint areas.

The new Mollelast haft was further optimised and provided with safe woven edges. The side threads are firmly anchored in the knit, so that they cannot become untied or fray.

Mollelast haft is also suitable as a simple support for sporting activities, as well as for cannula retention. If the dressing is fitted correctly, it will remain in position without cutting in, and is both air – permeable and well tolerated by the skin.



Mollelast® haft

Knitted bandage, white, 66% viscose, 34% polyamide, lengthwise elasticity approx. 65%, widthwise extensibility approx. 50%, made cohesive by a microfine latex coating.

Mollelast® haft color

Knitted bandage, red and blue, 56% viscose, 44% polyamide, lengthwise elasticity approx. 55%, widthwise extensibility approx. 50%, made cohesive by a microfine latex coating.

Mollelast® haft			
individually in cardboard box	Width	REF	Shipping Units
white 4 m (stretched length)	4 cm	30 063	1 pc. / 120 pcs.
	6 cm	30 064	1 pc. / 100 pcs.
	8 cm	30 065	1 pc. / 120 pcs.
	10 cm	30 066	1 pc. / 120 pcs.
	12 cm	30 067	1 pc. / 100 pcs.
white 20 m (stretched length)	4 cm	30 068	1 pc. / 18 pcs.
	6 cm	30 069	1 pc. / 18 pcs.
	8 cm	30 070	1 pc. / 18 pcs.
	10 cm	30 071	1 pc. / 18 pcs.
	12 cm	30 072	1 pc. / 18 pcs.
Mollelast® haft color			
individually in cardboard box			
red 20 m (stretched length)	6 cm	30 073	1 pc. / 18 pcs.
	8 cm	30 074	1 pc. / 18 pcs.
	10 cm	30 075	1 pc. / 18 pcs.
blue 20 m (stretched length)	6 cm	30 076	1 pc. / 18 pcs.
	8 cm	30 077	1 pc. / 18 pcs.
	10 cm	30 078	1 pc. / 18 pcs.



Velpeau® crêpe F • Velpeau® crêpe H

Crepe Bandages

These bandages with their special crepe character and high elasticity are slip-free and remain securely in position. They are suitable for retention, support, and for light compression, as well as for the treatment of distortions and contusions.

The bandages consist of 100% unbleached cotton, with an

open fabric structure, and are thus highly skin-friendly.

Velpeau® crêpe F

With blue edge threads, 100% unbleached cotton, lengthwise elastic, extensibility approx. 150%.

Velpeau® crêpe H

Same as Velpeau crêpe F, but with a higher surface weight.



Velpeau® crêpe F		Width	REF	Shipping Units
3 m individually packed		5 cm	17 797	30 pcs. / 840 pcs.
		7 cm	17 798	40 pcs. / 800 pcs.
		10 cm	17 799	40 pcs. / 480 pcs.
		15 cm	17 805	25 pcs. / 300 pcs.
		20 cm	17 806	20 pcs. / 240 pcs.
		25 cm	17 807	15 pcs. / 180 pcs.
	30 cm	17 808	10 pcs. / 140 pcs.	
3 m unwrapped		5 cm	17 809	30 pcs. / 840 pcs.
		7 cm	17 810	40 pcs. / 800 pcs.
		10 cm	17 811	40 pcs. / 480 pcs.
		15 cm	17 812	25 pcs. / 300 pcs.
		20 cm	17 813	20 pcs. / 240 pcs.
4 m individually packed		5 cm	17 814	25 pcs. / 700 pcs.
		7 cm	17 815	25 pcs. / 500 pcs.
		10 cm	17 816	30 pcs. / 360 pcs.
		15 cm	17 817	25 pcs. / 300 pcs.
		20 cm	17 818	15 pcs. / 180 pcs.
		25 cm	17 819	15 pcs. / 180 pcs.
	30 cm	17 820	10 pcs. / 140 pcs.	
5 m individually packed		20 cm	17 821	15 pcs. / 180 pcs.
4 m unwrapped		5 cm	17 822	25 pcs. / 700 pcs.
		7 cm	17 823	25 pcs. / 500 pcs.
		10 cm	17 824	30 pcs. / 360 pcs.
		15 cm	17 825	25 pcs. / 300 pcs.
		20 cm	17 826	15 pcs. / 180 pcs.
		30 cm	17 827	10 pcs. / 140 pcs.
Velpeau® crêpe H				
4 m (stretched length) individually packed		5 cm	15 297	25 pcs. / 700 pcs.
		7 cm	15 298	25 pcs. / 500 pcs.
		10 cm	15 299	30 pcs. / 360 pcs.
		15 cm	15 300	25 pcs. / 300 pcs.
		20 cm	15 301	15 pcs. / 180 pcs.
unwrapped	25 per box	5 cm	15 404	25 pcs. / 700 pcs.
	25 per box	7 cm	15 405	25 pcs. / 500 pcs.
	30 per box	10 cm	15 406	30 pcs. / 360 pcs.
	25 per box	15 cm	15 408	25 pcs. / 300 pcs.
	15 per box	20 cm	15 409	15 pcs. / 180 pcs.

Lenkelast®

Elastra Ideal medium stretch

Firm-grip, all purpose bandages with durable elasticity and with a wide range of applications: everyday retention bandages with compressive wound care, or splint dressings, for support

and pressure relief; also for light compression.

The extensibility of these bandages makes them conform well to all body contours, and the dressing remains in posi-

tion. The fine fabric is air-permeable and highly skin-friendly.

Lenkelast®

White, 64% cotton, 35% polyamide, 1% elastane, medium stretch, extensibility approx. 130%.

Elastra Ideal

White, 58% cotton, 41% polyamide, 1% elastane, medium stretch, extensibility approx. 135%, washable, sterilisable (A 121°C).



Lenkelast®	Width	REF	Shipping Units
5 m (stretched length) individually sealed	6 cm	19 460	10 pcs. / 250 pcs.
	8 cm	19 461	10 pcs. / 200 pcs.
	10 cm	19 462	10 pcs. / 240 pcs.
	12 cm	19 463	10 pcs. / 200 pcs.
	15 cm	19 464	10 pcs. / 160 pcs.
	20 cm	19 465	10 pcs. / 100 pcs.
5 m (stretched length) unwrapped 10 per box	6 cm	19 490	10 pcs. / 250 pcs.
	8 cm	19 491	10 pcs. / 200 pcs.
	10 cm	19 492	10 pcs. / 240 pcs.
	12 cm	19 493	10 pcs. / 200 pcs.
	15 cm	19 494	10 pcs. / 160 pcs.
	20 cm	19 495	10 pcs. / 100 pcs.
Elastra Ideal			
5 m (stretched length) unwrapped 10 per box	6 cm	10 122	10 pcs. / 160 pcs.
	8 cm	10 124	10 pcs. / 120 pcs.
	10 cm	10 108	10 pcs. / 120 pcs.
	12 cm	10 110	10 pcs. / 100 pcs.
	15 cm	10 116	10 pcs. / 80 pcs.
	20 cm	10 118	10 pcs. / 60 pcs.



Lenkideal®

short stretch

Short stretch bandage with durable elasticity, very similar to the DIN "Ideal" bandage, but with longer-lasting elasticity and more durable, easy-care and economical – a bandage for all situations: for compressive

retention of wound dressings; for the application of ointment dressings; for compression of the extremities in phlebology, for pre-, intra- and postoperative thrombo-prophylaxis; for support and relief in cases of

distortions and contusions as well as for sport supports; for treating tendosynovitis, and for securing splints.

Lenkideal with its flat fabric does not become bulky, is air-permeable and skin-friendly.



Lenkideal®

White, dense fabric, 64% cotton, 35% polyamide, 1% elastane, with durable elasticity, short stretch, extensibility approx. 90%.



Lenkideal®			
	Width	REF	Shipping Units
5 m (stretched length) with bandage clips individually sealed	6 cm	14 581	1 pc. / 100 pcs.
	8 cm	14 582	1 pc. / 140 pcs.
	10 cm	14 583	1 pc. / 120 pcs.
	12 cm	14 584	1 pc. / 120 pcs.
	15 cm	14 585	1 pc. / 120 pcs.
5 m (stretched length) without bandage clips unwrapped 10 per box	4 cm	19 580	10 pcs. / 140 pcs.
	6 cm	19 581	10 pcs. / 100 pcs.
	8 cm	19 582	10 pcs. / 140 pcs.
	10 cm	19 583	10 pcs. / 120 pcs.
	12 cm	19 584	10 pcs. / 120 pcs.
	15 cm	19 585	10 pcs. / 120 pcs.
	20 cm	19 586	10 pcs. / 60 pcs.

Idealbinde

short stretch

Bandage type “Ideal” conforming to DIN 61632, textile-elastic, the classic bandage for everyday use in hospital and doctor’s surgery; for powerful compression of the extremities in phlebology and lymphology; for pre-, intra-, and postoperative thrombo-prophylaxis; for sup-

port and relief in cases of distortions and contusions and for sport supports; for treating tendosynovitis, and for securing splints. “Ideal” bandages are air-permeable and skin-friendly; their firm-grip fabric structure ensures that the bandage turns do not slip.

Idealbinde

White fabric, conforming to DIN 61632, at least 85% cotton, viscose, textile-elastic through twisted cotton warp threads; short stretch, extensibility approx. 90%, can be washed and boiled.



Idealbinde			
	Width	REF	Shipping Units
5 m (stretched length) with bandage clips individually sealed	4 cm	14 100	1 pc. / 80 pcs.
	6 cm	14 101	1 pc. / 100 pcs.
	8 cm	14 102	1 pc. / 100 pcs.
	10 cm	14 103	1 pc. / 60 pcs.
	12 cm	14 104	1 pc. / 80 pcs.
	15 cm	14 105	1 pc. / 60 pcs.
	20 cm	14 106	1 pc. / 35 pcs.
	30 cm	14 108	1 pc. / 15 pcs.
5 m (stretched length) without bandage clips unwrapped 10 per box	4 cm	14 110	10 pcs. / 80 pcs.
	6 cm	14 111	10 pcs. / 100 pcs.
	8 cm	14 112	10 pcs. / 100 pcs.
	10 cm	14 113	10 pcs. / 60 pcs.
	12 cm	14 114	10 pcs. / 80 pcs.
	15 cm	14 115	10 pcs. / 60 pcs.
	20 cm	14 116	10 pcs. / 60 pcs.
	30 cm	14 118	10 pcs. / 40 pcs.



Durelast®

Durelast® Combi very short stretch

Bandage with very short stretch, exceptionally high compressive power and durable elasticity; for strong compression of the extremities in phlebology as well

as for the treatment of acute and chronic lymphoedemas; for support and relief with distortions and contusions in traumatology and sports medicine.

Durelast® Bandage

Firm-grip fabric, skin coloured, 66% cotton, 34% polyamide, with durable elasticity due to textured polyamide warp threads, very short stretch, extensibility approx. 45%, can be washed and boiled, sterilisable (A 121 °C).



Durelast® Bandage			
	Width	REF	Shipping Units
5 m (stretched length) with bandage clips individually in cardboard box	6 cm	22 100	1 pc. / 108 pcs.
	8 cm	22 101	1 pc. / 108 pcs.
	10 cm	22 102	1 pc. / 90 pcs.
	12 cm	22 103	1 pc. / 54 pcs.
5 m (stretched length) with bandage clips, unwrapped, 10 per box	6 cm	22 150	10 pcs. / 160 pcs.
	8 cm	22 151	10 pcs. / 80 pcs.
	10 cm	22 152	10 pcs. / 120 pcs.
	12 cm	22 153	10 pcs. / 80 pcs.
Durelast® Combi		Durelast bandages in varying widths in pairs, in cardboard box	
5 m (stretched length) with bandage clips	8 cm / 10 cm	22 110	1 box / 30 boxes

Rosidal® K • Elko®-Rosidal® K

Raucodur Kräftig (strong) short stretch

Short stretch bandages with high compressive power and durable elasticity; for strong compression of the extremities in phlebology as well as for the treatment of acute and chronic

lymphoedemas; for support and relief in traumatology and sports medicine.

The above bandages are permeable to air and well tolerated by the skin. The firm-grip fabric

texture ensures that the bandage turns hold securely.

The cohesive short stretch bandage Rosidal haft (see page 81) offers even greater protection against slipping.



Rosidal® K Raucodur Kräftig (strong)

Firm-grip fabric, skin coloured, 100% cotton, textile-elastic due to twisted warp threads, short stretch, extensibility approx. 90%, washable up to 95°C.

Elko®-Rosidal® K

For compressive dressings using the double bandaging technique, two Rosidal bandages of 10 cm width are available under the name of Elko-Rosidal K.

Rosidal® K			
	Width	REF	Shipping Units
5 m (stretched length) with bandage clips (except the bandage with 4 cm) individually in cardboard box	4 cm	22 199	1 pc. / 125 pcs.
	6 cm	22 200	1 pc. / 108 pcs.
	8 cm	22 201	1 pc. / 108 pcs.
	10 cm	22 202	1 pc. / 90 pcs.
	12 cm	22 203	1 pc. / 54 pcs.
10 m (stretched length) without bandage clips individually in cardboard box	10 cm	22 206	1 pc. / 60 pcs.
	12 cm	22 205	1 pc. / 40 pcs.
5 m (stretched length) with bandage clips unwrapped, 10 per box	8 cm	22 250	10 pcs. / 140 pcs.
	10 cm	22 251	10 pcs. / 140 pcs.
	12 cm	22 252	10 pcs. / 120 pcs.
Elko®-Rosidal® K 2 Rosidal K bandages packed in pairs in cardboard box			
5 m (stretched length) with bandage clips	10 cm	22 225	1 box / 30 boxes
Raucodur Kräftig (strong)			
5 m (stretched length) with bandage clips individually in cardboard box	6 cm	20 168	1 pc. / 50 pcs.
	8 cm	20 169	1 pc. / 50 pcs.
	10 cm	20 164	1 pc. / 50 pcs.
	12 cm	20 166	1 pc. / 50 pcs.
unwrapped 10 per box	8 cm	20 170	10 pcs. / 100 pcs.
	10 cm	20 165	10 pcs. / 100 pcs.
	12 cm	20 167	10 pcs. / 100 pcs.
10 m	10 cm	20 268	30 pcs.

Application	Lymphology Set 1 for leg, small REF 59 795	Lymphology Set 2 for leg, large REF 59 796	Lymphology Set 3 for arm, small REF 59 797	Lymphology Set 4 for arm, large REF 59 798	Lymphology Set 5 for leg, small REF 16 838	Lymphology Set 6 for leg, large REF 16 839
tg® Tubular Bandage For skin protection (hygiene) under decongestive compression dressings.	Size 7 (1x)	Size 9 (1x)	Size 5 (1x)	Size 7 (1x)	Size 9 (1x)	Size 9 (1x)
Mollelast® Bandage Comfortably soft, highly elastic conforming bandage for the compression of fingers and toes in lymphological compression dressings.	4 cm x 4 m (20x)	4 cm x 4 m (20x)	4 cm x 4 m (20x)	4 cm x 4 m (20x)	4 cm x 4 m (20x)	6 cm x 4 m (20x)
Komprex® Foam Rubber Bandage To cut out U- and L-shaped pieces (edges should be chamfered) for padding and increasing pressure in the ankle joint. Additionally, as basic material for padding packs (strips or balls) and for the local softening of tissue.	10 cm x 1 m Thickness 1 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)	8 cm x 1 m Thickness 0.5 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)
Cellona® Synthetic Padding Padding under compression bandages for the even distribution of pressure and for the prevention of constrictions and pressure spots.	6 cm x 3 m (8x) 10 cm x 3 m (8x)	10 cm x 3 m (8x) 15 cm x 3 m (8x)	6 cm x 3 m (4x) 10 cm x 3 m (4x)	6 cm x 3 m (4x) 10 cm x 3 m (4x)		
Rosidal® soft Skin-friendly, air-permeable foam bandage for padding under compression dressings. Prevents the dressing from slipping.					10 cm x 2.5 m Thickness 0.3 cm (2x) 15 cm x 2.5 m Thickness 0.4 cm (2x)	15 cm x 2.5 m Thickness 0.3 cm (2x) 15 cm x 2.5 m Thickness 0.4 cm (3x)
Rosidal® K Bandage (strong) Short stretch bandage, 100% cotton, for strong compression of extremities in lymphological treatment.	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (4x)	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (4x) 12 cm x 5 m (6x)	6 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (2x)	6 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (3x)	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (4x)	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (4x) 12 cm x 5 m (6x)
Idealbinde For strong compression and dressing stabilisation on the thigh and abdomen.					20 cm x 5 m (1x)	20 cm x 5 m (2x)

Dauerbinde® • Dauerbinde® biforte Perfekta® long stretch

Dauerbinde® K and F

Long stretch bandages with strong compression, with a characteristic firm-grip fabric texture in Dauerbinde K and a finer weave in Dauerbinde F.

Due to their long stretch properties, the compressive power of these bandages can be easily adjusted. They are used as support and relief bandages on ligaments and joints.

Because of the high resting pressure the bandage must be removed at night or when the patient is in a resting position.

The Dauerbinde can be fitted easily even without much experience in bandaging techniques; the bandages conform well to body contours, have little tendency to bulk up, are durable and easy-care.

Skin coloured, durable elasticity, long stretch, extensibility approx. 180%, washable at 95°C.



Dauerbinde® biforte

Two-way stretch elastic bandage. The bandage is indicated in sports injuries, in supporting and relieving sprains and contusions, for after-care in traumatology, and for strong and permanent compression in phlebology.

Skin coloured, bi-elastic long stretch, 144% in length, 84% in width, washable at 95°C.



Perfekta® Super (extra fine)

Long stretch bandage with durable elasticity, with low power of return; particularly well tolerated by the skin due to its high cotton content.

Support and relief bandage for prophylactic relief of capsules and ligaments, prevention of haematomas and for contusions, distortions, dislocations and tendon injuries.

Perfekta Super can also be used as a compressive bandage.

Skin coloured, resistant to ointments, long stretch, extensibility approx. 160%, washable at 95°C.

Perfekta® Kräftig (strong)

Long stretch bandage with durable elasticity and strong power of return; particularly well tolerated by the skin due to its high cotton content.

Compression bandages for varicosis, oedemas, angiodysplasia, sclerotherapy, following surgical interventions on veins and for prevention of thrombosis. Support and relief bandages for phlebitis or existing injuries.

Perfekta® Fein (fine)

Same as Perfekta Kräftig, but with a ribbed surface and medium power of return.

Skin coloured, resistant to ointments, long stretch, extensibility approx. 170%, washable at 95°C.



Dauerbinde® K		86% cotton, 7% polyamide, 7% elastane		
	Width	REF	Shipping Units	
7 m (stretched length) with bandage clips individually in cardboard box	6 cm	22 000	1 pc. / 100 pcs.	
	8 cm	22 001	1 pc. / 100 pcs.	
	10 cm	22 002	1 pc. / 100 pcs.	
	12 cm	22 003	1 pc. / 50 pcs.	
	20 cm	22 004	1 pc. / 20 pcs.	
14 m (stretched length) with bandage clips individually sealed	10 cm	22 007	1 pc. / 25 pcs.	
	12 cm	22 008	1 pc. / 25 pcs.	
7 m (stretched length) with bandage clips unwrapped 10 per box	6 cm	22 020	10 pcs. / 200 pcs.	
	8 cm	22 021	10 pcs. / 180 pcs.	
	10 cm	22 022	10 pcs. / 120 pcs.	
	12 cm	22 023	10 pcs. / 100 pcs.	
	20 cm	22 024	10 pcs. / 60 pcs.	
Dauerbinde® F		82% cotton, 13% polyamide, 5% elastane		
7 m (stretched length) with bandage clips individually in cardboard box	6 cm	22 010	1 pc. / 100 pcs.	
	8 cm	22 011	1 pc. / 100 pcs.	
	10 cm	22 012	1 pc. / 100 pcs.	
	12 cm	22 013	1 pc. / 50 pcs.	
7 m (stretched length) with bandage clips unwrapped 10 per box	6 cm	22 025	10 pcs. / 160 pcs.	
	8 cm	22 026	10 pcs. / 80 pcs.	
	10 cm	22 027	10 pcs. / 120 pcs.	
	12 cm	22 028	10 pcs. / 80 pcs.	
Dauerbinde® biforte		69% viscose, 24% cotton, 7% elastodiene		
7 m (stretched length) individually sealed	8 cm	16 980	1 pc. / 40 pcs.	
	10 cm	16 981	1 pc. / 30 pcs.	
	12 cm	16 982	1 pc. / 28 pcs.	
Perfekta® Super (extra fine)		96% cotton, 4% elastane		
7 m (stretched length) unwrapped 8 per box	8 cm	20 136	8 pcs. / 80 pcs.	
	10 cm	20 131	8 pcs. / 80 pcs.	
	12 cm	20 133	8 pcs. / 80 pcs.	
Perfekta® Kräftig (strong)		94% cotton, 6% elastane		
7 m (stretched length) unwrapped 8 per box	8 cm	20 129	8 pcs. / 80 pcs.	
	10 cm	20 109	8 pcs. / 80 pcs.	
	12 cm	20 117	8 pcs. / 80 pcs.	
Perfekta® Fein (fine)		94% cotton, 6% elastane		
7 m (stretched length) individually in cardboard box	8 cm	20 103	1 pc. / 50 pcs.	
	10 cm	20 099	1 pc. / 50 pcs.	
	12 cm	20 100	1 pc. / 50 pcs.	
10 m (stretched length)	10 cm	20 098	1 pc. / 50 pcs.	
	12 cm	93 185	1 pc. / 50 pcs.	



tg[®] grip Elasticated Tubular Support Bandage

Tubular bandage with durable elasticity, supple, exerting even pressure; highly practical and handy for dressing retention with light compression, for support and relief on ligaments,

tendons and joints and for wide-area (thorax) dressings.

The bandage is air-permeable, slightly padding and very skin-friendly through its high cotton

content. tg grip can be fitted, removed and re-applied easily and without creasing. The compression effect can be increased with a second layer or by selecting a smaller size.



tg[®] grip Elasticated Tubular Support Bandage

Seamlessly knitted tube, raw white, 85% cotton, 10% elastodiene (rubber), covered with 5% polyamide yarn, widthwise elastic, can be washed (60°C) and sterilised (A 134°, EO, irradiation).



tg[®] grip Elasticated Tubular Support Bandage

10 m roll in cardboard dispenser	Size	REF	Shipping Units
for children's hands, children's forearms	A (4.6 cm)	24 120	1 roll / 24 rolls
for children's arms, children's lower legs	B (6.0 cm)	24 121	1 roll / 16 rolls
for slim arms, children's legs	C (6.5 cm)	24 122	1 roll / 16 rolls
for medium-sized arms, lower legs	D (7.5 cm)	24 123	1 roll / 16 rolls
for large arms, legs, slim thighs	E (8.75 cm)	24 124	1 roll / 12 rolls
for legs, medium-sized thighs	F (10 cm)	24 125	1 roll / 12 rolls
for legs, for large thighs	G (12 cm)	24 126	1 roll / 8 rolls
for small trunks	J (17.5 cm)	24 127	1 roll / 6 rolls
for medium-sized trunks	K (21.5 cm)	24 128	1 roll / 4 rolls
for large trunks	L (32.5 cm)	24 129	1 roll / 4 rolls
for extra-large trunks	M (37.5 cm)	24 130	1 roll / 2 rolls

Varicex® Zinc Paste Bandages (Unna's Boot) Starched Gauze Bandages

Bandages with different fabric texture, impregnated with zinc paste; ready-for-use, time-saving application of semi-rigid permanent bandages for the treatment of thrombophlebitis, oedemas, and chronic venous insufficiency. Also used for leg ulcers in the healing phase and for the

aftercare of fractures. Whether they are applied in circular turns or with the semi-figure-of-eight technique – Varicex zinc paste bandages do not bulk up.

Varicex® S

Bi-elastic knitted bandage for application in circular turns. It is par-

ticularly easy to apply, because of the material's stretch effect.

Varicex® E

Elastic creped bandage for application in circular turns.

Varicex® F

Produced from gauze bandages, with a low moisture content, for rapid drying.

Varicex® T

Produced from gauze bandages, with a high moisture content, for slow drying.

Starched Gauze Bandage

20-thread, starched, for light rigid bandaging or for the reinforcement of bandages with an immobilising function.



Varicex® S		
individually in aluminium pouch	REF	Shipping Units
10 cm x 5 m (stretched length)	20 630	1 pc. / 80 pcs.
10 cm x 7 m (stretched length)	20 631	1 pc. / 60 pcs.
10 cm x 10 m (stretched length)	20 632	1 pc. / 30 pcs.
individually in aluminium pouch, 20 pcs. in a cardboard box		
10 cm x 5 m (stretched length)	20 635	20 pcs.
10 cm x 7 m (stretched length)	20 636	20 pcs.
10 cm x 10 m (stretched length)	20 637	20 pcs.
Varicex® E individually in aluminium pouch		
8 cm x 5 m (stretched length)	20 620	1 pc. / 60 pcs.
Varicex® F individually in aluminium pouch		
10 cm x 5 m	20 600	1 pc. / 80 pcs.
10 cm x 7 m	20 601	1 pc. / 60 pcs.
10 cm x 10 m	20 602	1 pc. / 50 pcs.
Varicex® T individually in aluminium pouch		
10 cm x 10 m	20 612	1 pc. / 50 pcs.
Starched Gauze Bandage unwrapped, 10 per box		
6 cm x 4 m	14 701	10 pcs. / 100 pcs.
8 cm x 4 m	14 702	10 pcs. / 100 pcs.
10 cm x 4 m	14 703	10 pcs. / 100 pcs.
12 cm x 4 m	14 704	10 pcs. / 50 pcs.
15 cm x 4 m	14 705	10 pcs. / 50 pcs.

In ready-for-use formats of different sizes, with lasting vertical elasticity; for compression in the treatment of oedemas in phlebology and lymphology; also used for functional bandages in sports medicine.

The fine-pored material with its velvety surface can also be easily cut to the required shape.

Free of natural latex.



Komprex® Foam Rubber Pads

Synthetic latex (styrene-butadiene), permeable to air and water vapour.

Komprex® Foam Rubber Bandages

Komprex bandages are indicated for severe oedemas, thrombophlebitis and after extensive sclerotherapy of varices. The bandages are applied in spiral turns, without stretching, and edge-to-edge – i.e. without overlapping; the increased compression effect is achieved by applying elastic bandages on top. If a

dressing is cut to a specific shape from a sheet of Komprex or from a bandage (e.g. U- and L-shaped pieces), the edges should be chamfered to prevent pressure spots.

Pads, bandages or material cut to size can also be used as partial padding with dressings and splints.

Synthetic latex (styrene-butadiene), permeable to air and water vapour.

Komprex® II Foam Plate, wavy

Narrow strips of foam material, embedded in a skin-friendly nonwoven, offer an extremely breathable, non-slip padding material. The special, wavy surface structure of Komprex II creates a firm yet gentle micro-massaging effect.

Komprex® Foam Rubber Pads				
Size	Thickness	Size	REF	Shipping Units
individually in cardboard box				
00 oval	10 mm	9 x 6.5 cm	22 300	1 pc. / 60 pcs.
0 kidney-shaped	10 mm	9 x 5 cm	22 301	1 pc. / 75 pcs.
1 kidney-shaped	12 mm	12 x 6.5 cm	22 302	1 pc. / 30 pcs.
2 rectangular	12 mm	17 x 13 cm	22 303	1 pc. / 25 pcs.
4 trapezoid	5 mm	39/22 x 33 cm	22 305	1 pc. / 20 pcs.
Komprex® Foam Rubber Bandages				
	5 mm	8 cm x 1 m	22 310	1 pc. / 30 pcs.
rolled	10 mm	8 cm x 1 m	22 311	1 pc. / 20 pcs.
individually	10 mm	10 cm x 1 m	22 312	1 pc. / 20 pcs.
in cardboard box	5 mm	8 cm x 2 m	22 313	1 pc. / 30 pcs.
	10 mm	8 cm x 2 m	22 314	1 pc. / 15 pcs.
Komprex® Foam Rubber Sheet				
rolled, in plastic pouch	10 mm	100 x 50 cm	22 320	1 pc. / 10 pcs.
Komprex® II Foam Plate, wavy				
individually packed in a pouch		65 x 65 cm	16 985	1 pc. / 6 pcs.

Rosidal® sys

Short Stretch Compression System

For an effective compression therapy of venous ulcers.

The use of this set in long-term treatment is particularly effective and economical. All components are carefully adjusted to

each other, are of high and lasting quality and are combined in a handy set.

The reusability of the padding and compression bandages results in a considerable cost-ef-

fectiveness compared with disposable products. The padding bandage Rosidal soft and the short stretch bandage Rosidal K are washable up to 50 times at up to 95°C.

A system solution also saves handling time and offers highest therapy safety. The administration of individual components is eliminated.

For wound treatment, the dressing can be adequately added, so that the wound is optimally cared for.

It is possible to schedule bandage changes individually for each patient. Thus, a compression bandage can be worn for up to 7 days.



Rosidal® sys			
Components	Size		REF
Set	Shipping Units 1 piece/4 pieces		16 984
All components may be ordered individually (except those marked with *).			
tg® tubular bandage For skin protection under the decongesting compression dressing	Size 7	1x	24 005
Rosidal® soft bandage Skin-friendly, latex-free foam bandage made from 100% polyurethane foam; washable up to 95°C	10 cm x 2 m x 0.2 cm	4x	13 016
Rosidal® K bandage Tightly compressing short stretch bandage made from 100% cotton; washable up to 95°C	8 cm x 5 m	2x	22 201
	10 cm x 5 m	2x	22 202
Porofix® Skin coloured adhesive tape on the roll; coated with synthetic rubber adhesive	2.5 cm x 5 m	1x	34 311
Mollelast® haft bandage Soft cohesive conforming bandage, lengthwise and widthwise elastic	8 cm x 20 m	1x	30 070
Laundry bag*		1x	

Rosidal® soft

Rosidal soft is ideally suited for padding underneath lymphological and phlebological compression bandages. The padding bandage leads to an even pressure distribution of short stretch bandages (e.g. Rosidal K) over the extremity, this makes the compression treatment safe by preventing constrictions. The bandage will not slip as the open pores of Rosidal soft are interlocking.

The skin-friendly, latex-free bandage made from polyurethane foam is highly permeable to air and moisture.

Rosidal soft is particularly tear-resistant and does not lose its elasticity even after repeated washing.

(We recommend a washing temperature of 40 - 60°C. If the bandage is extremely dirty, it may be boiled.)



100% polyurethane foam, highly extensible, washable up to 95°C, sterilisable (A 121°C, A 134°C).



Rosidal® soft			
	Size	REF	Shipping Units
in bulk (unwrapped, with bandage clips)	10 cm x 0.3 cm x 2.5 m	23 100	1 box / 36 pcs.
	10 cm x 0.4 cm x 2.5 m	23 101	24 pcs.
	12 cm x 0.4 cm x 2.5 m	23 102	18 pcs.
	15 cm x 0.4 cm x 2.5 m	23 103	14 pcs.
individually packed (in cardboard box, with bandage clips)	10 cm x 0.3 cm x 2.5 m	23 110	1 pc. / 30 pcs.
	10 cm x 0.4 cm x 2.5 m	23 111	1 pc. / 24 pcs.
	12 cm x 0.4 cm x 2.5 m	23 112	1 pc. / 16 pcs.
	15 cm x 0.4 cm x 2.5 m	23 113	1 pc. / 14 pcs.
two pieces in a cardboard box (with bandage clips)	10 cm x 0.2 cm x 2 m	13 016	1 box / 23 boxes

Dauerlastic® TPS

Thrombo-prophylactic Stockings

Seamless, circular knitted stocking for thrombo-prophylaxis.

The stocking is widthwise elastic and shaped anatomically in such a way that the compression applied to the leg and the veins in the leg falls gradually from 18 mmHg at the ankle to a

maximum of 8 mmHg at the thigh. The foot and ankle are firmly enclosed.

Dauerlastic TPS are provided with an inspection opening below the toes. Dauerlastic TPS stockings assist the blood to flow back through the veins in

immobilised patients before, during and after operations, and in obstetrics.

Seamless, circular knitted, white, 80% polyamide, 20% elastane; flat garter strip with silicone rubber coating, can be boiled and sterilised (A 121°C).



The compression values are unaffected even by repeated washing.

The stocking sizes fit approx. 90% of patients; if the given values fall below or above the specified dimensions, compression bandages (e. g. Dauerbinde K or F) should be used instead.

Contraindications: arterial disorders which act haemodynamically; severe cardiac insufficiency; serious oedemas in the leg; dimensions of the thigh which exceed the maximum dimensions of the stockings; care should be exercised with venous ligatures.

Dauerlastic® TPS

Thigh length, one pair per pack

Size	Colour code foot	Thigh length to gluteal fold	Calf circumference	Thigh circumference (max.)	REF	Shipping Units
A1	1 blue stripe	up to 75 cm	up to 34.5 cm	up to 65 cm	23 570	1 pack / 26 packs
A2	2 blue stripes	75-85 cm	up to 34.5 cm	up to 65 cm	23 571	1 pack / 26 packs
A3	3 blue stripes	over 85 cm	up to 34.5 cm	up to 65 cm	23 572	1 pack / 26 packs
B1	1 green stripe	up to 75 cm	34.5-40.5 cm	up to 70 cm	23 573	1 pack / 26 packs
B2	2 green stripes	75-85 cm	34.5-40.5 cm	up to 70 cm	23 574	1 pack / 26 packs
B3	3 green stripes	over 85 cm	34.5-40.5 cm	up to 70 cm	23 575	1 pack / 26 packs
C1	1 yellow stripe	up to 75 cm	over 40.5 cm	up to 75 cm	23 580	1 pack / 26 packs
C2	2 yellow stripes	75-85 cm	over 40.5 cm	up to 75 cm	23 581	1 pack / 26 packs
C3	3 yellow stripes	over 85 cm	over 40.5 cm	up to 75 cm	23 582	1 pack / 26 packs

Knee length, one pair per pack

Size	Colour code foot	Knee length to hollow of the knee	Calf circumference	REF	Shipping Units
A	1 blue stripe	up to 40 cm	up to 34.5 cm	22 530	1 pack / 30 packs
B	1 green stripe	40-43 cm	34.5-40.5 cm	22 531	1 pack / 30 packs
C	1 yellow stripe	up to 47 cm	over 40.5 cm	22 532	1 pack / 30 packs