



Full Service in the Plaster Room



Cellona® Plaster of Paris Bandages

Ready-to-dip Plaster of Paris bandages with exceptionally fine creaminess giving enhanced mouldability; for casts applied in circular turns, for im-

mobilisation after fractures, operations and orthopaedic corrections; for treating joint and bone disorders. Cellona Plaster of Paris bandages are fully



soaked after brief immersion. The finished cast has an early weight bearing after approx. 30 minutes; it is fully set after approx. 24 hours.

Cellona® P.o.P. Bandages

Gauze bandage, porously coated with plaster compound made from extremely fine α and β semi-hydrate.

The bandages are sealed in packaging impermeable to water vapour.

Cellona® P.o.P. Bandages		REF	Shipping Units
2 m 5 bandages in a bag	6 cm	25 400	60 pcs.
	8 cm	25 401	60 pcs.
	10 cm	25 402	40 pcs.
	12 cm	25 403	40 pcs.
	15 cm	25 404	40 pcs.
	20 cm	25 405	40 pcs.
3 m	8 cm	25 411	40 pcs.
	10 cm	25 412	40 pcs.
	12 cm	25 413	40 pcs.
	15 cm	25 414	40 pcs.
	20 cm	25 415	30 pcs.
4 m	8 cm	25 421	40 pcs.
	10 cm	25 422	40 pcs.
	12 cm	25 423	30 pcs.
	15 cm	25 424	20 pcs.
	20 cm	25 425	20 pcs.
1 m, individually sealed	5 cm	24 999	10 pcs. / 200 pcs.
2.7 m individually in a bag	5 cm	91 380	12 pcs. / 96 pcs.
	7.5 cm	91 381	12 pcs. / 72 pcs.
	10 cm	91 382	12 pcs. / 48 pcs.
	12 cm	91 383	12 pcs. / 48 pcs.
	15 cm	91 384	12 pcs. / 48 pcs.
	20 cm	91 385	12 pcs. / 36 pcs.

Cellona® Plaster of Paris Splints Wide Material

Cellona P.o.P., ready-to-dip, folded in 4 layers, or as wide material off the roll; for splints and shell-shaped plaster casts and cast sections, e.g. casts for thorax and combined thorax and

extremities, pelvis/leg casts, plaster corsets and plaster beds; also plaster casts for orthotic devices and supports. The splinting technique shortens the application time, re-

duces material consumption and is also suitable for supplementing or reinforcing casts applied in circular turns.

Cellona® P.o.P. Splints Cellona® Wide Material

Gauze bandage, porously coated with plaster compound made from extremely fine α and β semi-hydrate; 4-fold splints, zigzag folded, in PE bag impermeable to water vapour, in cardboard dispenser.



Cellona® Splints		REF	Shipping Units
4-ply			
20 m in a foil pouch and cardboard dispenser	10 cm	25 500	1 pc.
	12 cm	25 501	1 pc.
	15 cm	25 502	1 pc.
	20 cm	25 503	1 pc.
Cellona® Wide Material			
5 m roll individually in a foil pouch	40 cm width	20 301	1 pc. / 10 pcs.
	60 cm width	20 302	1 pc. / 10 pcs.
	80 cm width	20 303	1 pc. / 10 pcs.



Cellona® Xtra P.o.P. Bandages and Splints

The main difference between Cellona and Cellona Xtra Plaster of Paris is the shorter setting time of Cellona Xtra.

Simple dip technique: with easy-soak through properties and low plaster loss.

Enhanced mouldability: extra plaster with extra creaminess. Guaranteed 2.5 minutes setting time.

Windowing and bivalving after 10 minutes demonstrate excellent early strength.

Partial weight bearing after half an hour, with assistance of walking aid.

Full weight bearing after 24 hours.

Cellona Xtra Plaster of Paris bandages and splints are indicated for the immobilisation of fractures, postoperatively after orthopaedic surgical procedures, and in the treatment of bone and joint diseases.



Cellona® Xtra Plaster of Paris Bandages

	Size	REF	Shipping Units
individually sealed	5 cm x 2.7 m	91 180	1 box / 60 pcs.
	7.5 cm x 2.7 m	91 181	1 box / 60 pcs.
	10 cm x 2.7 m	91 182	1 box / 60 pcs.
	15 cm x 2.7 m	91 184	1 box / 36 pcs.
	20 cm x 2.7 m	91 185	1 box / 36 pcs.
	6 cm x 3 m	25 050	1 box / 60 pcs.
	8 cm x 3 m	25 051	1 box / 60 pcs.
	10 cm x 3 m	25 052	1 box / 60 pcs.
	12 cm x 3 m	25 053	1 box / 48 pcs.
	15 cm x 3 m	25 054	1 box / 36 pcs.
	20 cm x 3 m	25 055	1 box / 30 pcs.

Cellona® Xtra Plaster of Paris Splints, 4-ply

each splint sealed in a protective PE bag and packed to cardboard dispenser	10 cm x 20 m	25 056	1 box / 1 splint
	12 cm x 20 m	25 057	1 box / 1 splint
	15 cm x 20 m	25 058	1 box / 1 splint
	20 cm x 20 m	25 059	1 box / 1 splint

Cellacast® Xtra Casting Tape and Splints

Synthetic material for exceptionally lightweight supports with high strength; for immobilisation after fractures, operations and orthopaedic corrections; for treatment of joint and bone disorders; ideal for pro-



moting mobilisation and for bandaging child patients.

Cellacast Xtra bandages and splints are economical in material consumption. They have widthwise and diagonal stretch for easier bandaging and moulding.

All the techniques familiar from Plaster of Paris apply equally to Cellacast Xtra, e.g. for opening, splitting, windowing, combining with walking heels, hinging etc.

Cellacast Xtra bandages are lightweight, permeable to air, X-ray transparent and waterproof once hardened.



Cellacast® Xtra

Knitted fibreglass fabric impregnated with polyurethane resin; bandages wound on a reel core (PP), splints packed in 4 layers and individually sealed in aluminium composite foil; unpadded pre-cuts for applying splints and for the reinforcement of circular casts.

Cellacast® Xtra Casting Tape			
individually sealed	Size	REF	Shipping Units
Creme	2.5 cm x 1.8 m	25 205	10 pcs. / 40 pcs.
	5.0 cm x 3.6 m	25 201	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 202	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 203	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 204	10 pcs. / 40 pcs.
Blue	5.0 cm x 3.6 m	25 211	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 212	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 213	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 214	10 pcs. / 40 pcs.
Orange	5.0 cm x 3.6 m	25 235	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 236	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 237	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 238	10 pcs. / 40 pcs.
Yellow	5.0 cm x 3.6 m	25 239	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 240	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 241	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 242	10 pcs. / 40 pcs.
Green	5.0 cm x 3.6 m	25 221	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 222	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 223	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 224	10 pcs. / 40 pcs.
Cellacast® Xtra Splints			
4-ply, individually sealed, Creme	5 cm x 20 cm	25 660	1. pc. / 15 pcs. / 60 pcs.
	7.5 cm x 38 cm	25 261	1. pc. / 10 pcs. / 40 pcs.
	10.0 cm x 38 cm	25 262	1. pc. / 10 pcs. / 40 pcs.
	10.0 cm x 96 cm	25 263	1. pc. / 5 pcs. / 20 pcs.
	12.5 cm x 76 cm	25 264	1. pc. / 5 pcs. / 20 pcs.

Cellacast® Splints

Synthetic casting tape, padded on one side. For a clean, time saving and economical splinting technique; for functional and partially immobilising casts, e.g.

after distortions, fibular ligament rupture; in cases of tendovaginitis, epicondylitis humeri or distal fracture of the forearm.

Cellacast® Splints

7 layers of knitted fibreglass, impregnated with polyurethane resin to form the mouldable core of the splint, are integrated between fine padding felt (100% polypropylene) on the skin contact side and a film/padding cover. Cellacast splints are sealed in a moisture-proof aluminium composite foil.



Cellacast® Splints			
in a cardboard dispenser		REF	Shipping Units
folded length 3.6 m	5.0 cm width	25 039	1 box / 2 boxes
	7.5 cm width	25 040	1 box / 2 boxes
	10.0 cm width	25 041	1 box / 2 boxes
	12.5 cm width	25 042	1 box / 2 boxes
	15.0 cm width	25 043	1 box / 2 boxes
Cellacast® Pre-cut Splints			
individually sealed	7.5 x 30 cm	25 250	10 pcs. / 40 pcs.
	10.0 x 38 cm	25 251	10 pcs. / 40 pcs.
	10.0 x 76 cm	25 252	5 pcs. / 20 pcs.
	12.5 x 76 cm	25 253	5 pcs. / 20 pcs.

Cellacast® Active Bandages

Cellacast® Hook and Loop Tapes

Cast without fibreglass for focussed immobilisation. Cellacast Active allows an ideal combination of the necessary immobilisation and controlled movement. For immobilisation in primary and secondary treatment after fractures, operations, tendon injuries and overload

conditions. Cellacast Active is particularly suitable to promote mobilisation, also for children's casts.

Cellacast Active bandages can be handled fast and easily. The thin carrier material makes the bandage highly conformable to body contours. All the usual

manipulations used in casting techniques – from cast application to bivalving, windowing, and a combination with Shoe-cast – apply equally to Cellacast Active.

Cellacast Active casts are thin, light, 100% X-ray transparent, and can be easily removed and worked on with cast scissors.



Cellacast® Active

Polyester carrier material, impregnated with polyurethane resin; Bandages are rolled on a PP core and individually sealed in an laminated aluminium foil.

Cellacast® Hook and Loop Tapes

For splint fixation in focussed immobilisation, available in two versions.

Both tapes made of polyamide, dark blue, off the roll, in a dispenser.

Cellacast® Active			
individually sealed	Size	REF	Shipping Units
Creme	2.5 cm x 1.8 m	23 034	10 pcs. / 40 pcs.
	5.0 cm x 3.6 m	23 035	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 036	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 037	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 038	10 pcs. / 40 pcs.
Blue	5.0 cm x 3.6 m	23 040	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 041	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 042	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 043	10 pcs. / 40 pcs.
Orange	5.0 cm x 3.6 m	23 045	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 046	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 047	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 048	10 pcs. / 40 pcs.
Cellacast® Hook Tape self -adhesive			
dark blue	2.5 cm x 25 m	18 040	1 pc. / 24 pcs.
Cellacast® Loop Tape non -adhesive			
dark blue	2.5 cm x 25 m	18 041	1 pc. / 24 pcs.

Cellacast® Soft Casting Tape

Semi-rigid casting tape for focussed immobilisation. This casting tape is suitable for making semi-rigid casts with high wearer comfort. In combination with a Cellacast Xtra splint, it allows a precise, selective increase in stability. For immobilisation in primary and secondary care of bone fractures, in postoperative care, or in the treatment of ligament injuries and conditions due to excessive physical stress. Cellacast Soft is particularly appropriate for use in the functional

mobilisation of patients and for children's casts.

Cellacast Soft tapes are quick and easy to apply. The widthwise and diagonal stretch of the material makes the application easy and convenient; its fine structure allows smooth and accurate moulding to body contours.

All common P.o.P. techniques—from application to splitting, windowing, or the combination with Cellona Shoecast—can be applied to Cellacast Soft as well.

Cellacast Soft casts are thin, lightweight, and X-ray transparent; they can be removed and trimmed to size by use of cast scissors.

Cellacast® Soft

Knitted fibreglass fabric impregnated with polyurethane resin; rolled on a reel core (PP), individually sealed in aluminium composite foil.



Cellacast® Soft Casting Tape			
	Size	REF	Shipping Units
Individually sealed, creme	2.5 cm x 1.8 m	23 721	10 pcs. / 40 pcs.
	5.0 cm x 3.6 m	23 722	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 723	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 724	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 725	10 pcs. / 40 pcs.
blue	5.0 cm x 3.6 m	23 726	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 727	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 728	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 729	10 pcs. / 40 pcs.

tg[®] soft Tubular Padding Bandage

tg soft is the fast and safe padding alternative for P.o.P. and synthetic casts – particularly for focussed immobilisation. The soft terry cloth bandage adheres to the cast. This makes it possible to cut the cast easily to individual needs. The highly extensible bandage can be applied crease-free, does not slip and conforms to any body contour.

The padding properties remain unimpaired over the whole period of wear.

High water vapour and air permeability as well as a low absorption volume of the tightly-meshed bandage guarantee a convincing wearer comfort. The high cotton content makes tg soft very skin-friendly.

Only three sizes to cover all areas of application for upper and lower extremities.



tg[®] soft

85% cotton, 15% elastane, unbleached;

Tubular padding bandage for padding under P.o.P. and synthetic casts, easy, practical, safe, economic.

tg[®] soft Tubular Padding Bandage

10 m roll in a cardboard dispenser	Size	REF	Shipping Units
for forearm, child's arm, child's leg	S	13 200	1 pc. / 8 pcs.
for lower leg, arm	M	13 201	1 pc. / 6 pcs.
for leg	L	13 202	1 pc. / 4 pcs.

Plain Stockinette

Very supple and skin-friendly undercast protection for P.o.P. and synthetic casts. The seamlessly knitted tube made from raw white cotton yarns is elastic, especially widthwise.

This stockinette is also used for crease-free skin protection on amputation stumps.



Plain Stockinette

100% cotton, raw white; sterilisable (A 121 °C).

Plain Stockinette			
	Width	REF	Shipping Units
25 m roll individually in a foil pouch	6 cm	14 820	1 roll / 10 rolls
	8 cm	14 821	1 roll / 10 rolls
	10 cm	14 822	1 roll / 5 rolls
	12 cm	14 823	1 roll / 5 rolls
	15 cm	14 824	1 roll / 6 rolls
	20 cm	14 825	1 roll / 3 rolls
	25 cm	14 826	1 roll / 3 rolls
	30 cm	14 827	1 roll / 3 rolls
	35 cm	14 828	1 roll / 3 rolls
4 m	4 cm	14 799	1 pc. / 110 pcs.
	6 cm	14 800	1 pc. / 80 pcs.
	8 cm	14 801	1 pc. / 48 pcs.
	10 cm	14 802	1 pc. / 40 pcs.
	12 cm	14 803	1 pc. / 32 pcs.
	15 cm	14 804	1 pc. / 24 pcs.
	20 cm	14 805	1 pc. / 20 pcs.
	25 cm	14 806	1 pc. / 18 pcs.

tg[®] Tubular Bandage

A fine, seamlessly knitted tube with widthwise extensibility for use as crease-free undercast protection – from the little finger through to body-size P.o.P. or synthetic casts; also suitable to hold padding material in place, and for dermatological indications.

tg[®] Tubular Bandage

67% bleached cotton, 33% viscose, sterilisable (A 121°C).



tg tubular bandage is also used for dermatological indications, over ointment dressings or zinc paste bandages or as undercast protection for P.o.P. or synthetic casts.

tg[®] Tubular Bandage				REF	Shipping Units
white, individually in a cardboard dispenser					
5 m roll	size 1	(1.4 cm)	for finger and toe	24 020	1 roll / 80 rolls
	size 2	(2.3 cm)	for fingers and toes	24 021	1 roll / 50 rolls
	size 3	(3 cm)	for several fingers, child's arm	24 022	1 roll / 40 rolls
	size 5	(5.5 cm)	for arm, lower leg, child's leg	24 023	1 roll / 40 rolls
	size 6	(6.5 cm)	for arm, leg (medium)	24 024	1 roll / 40 rolls
	size 7	(7 cm)	for arm (large), leg, child's head	24 025	1 roll / 40 rolls
	20 m roll	size 1	(1.4 cm)	for finger and toe	24 000
size 2		(2.3 cm)	for fingers and toes	24 001	1 roll / 20 rolls
size 3		(3 cm)	for several fingers, child's arm	24 002	1 roll / 10 rolls
size 5		(5.5 cm)	for arm, lower leg, child's leg	24 003	1 roll / 10 rolls
size 6		(6.5 cm)	for arm, leg (medium)	24 004	1 roll / 10 rolls
size 7		(7 cm)	for arm (large), leg, child's head	24 005	1 roll / 18 rolls
size 9		(8.5 cm)	for thigh (large), head, armpit	24 006	1 roll / 12 rolls
size 12		(12 cm)	for thigh (extra large), head (large), child's trunk	24 007	1 roll / 12 rolls
10 m roll	size K1	(16 cm)	for trunk (small) < M	24 030	1 roll / 10 rolls
	size K2	(21 cm)	for trunk (large) > L	24 031	1 roll / 10 rolls

tg® Vest

Short-sleeved vest, soft, supple and skin-friendly, as an easy-to-apply and time-saving under-cast protection for P.o.P. and synthetic casting tapes. The elastic knitted material has widthwise extensibility, allow-

ing for intermediate sizes. If two vests are applied – and a special casting technique is used – the vest that is in contact with the skin can be changed, an advantage for body hygiene under long-term casts.

tg vests are also used in dermatology to cover large areas of affected skin.

100% cotton, not degraded or optically brightened, sterilisable (A 121°C).



tg® Vest

	REF	Shipping Units
individually in a cardboard box		
small (ladies' size 36 - 44, men 42 - 50)	24 044	1 pc. / 14 pcs.
large (ladies' up to size 46, men up to size 56)	24 045	1 pc. / 9 pcs.

Haftan®

Crepe Paper Bandages

Haftan®

Extra-thin, supple bandage made from white polyurethane foam which does not slip and creates a smooth, even surface, holding padding material in place under casting tapes, es-

pecially in combination with synthetic material.

The fine, open-pored surface of Haftan allows the individual turns of the bandage to adhere well to each other, thus preventing slippage even on the skin.

The PU foam reduces any penetration of moisture.

Haftan is also used for joint supports (sports), combined with adhesive bandages or tape.



Crepe Paper Bandage

Fine, crimped paper bandage, lengthwise slightly stretchable. Allows crease-free application; for protecting the skin under unpadded plaster casts and for holding padding material in place. The bandage retains its strength even when wet.

100% cellulose

Haftan®		REF	Shipping Units
white,	7 cm x 27.5 m	22 315	6 pcs. / 60 pcs.
thickness: 0.8 mm	10 cm x 27.5 m	22 316	6 pcs. / 60 pcs.
Crepe Paper Bandages			
4 m	4 cm	14 250	20 pcs. / 200 pcs.
	6 cm	14 251	20 pcs. / 200 pcs.
	8 cm	14 252	20 pcs. / 200 pcs.
	10 cm	14 253	20 pcs. / 200 pcs.
	12 cm	14 254	20 pcs. / 200 pcs.



Cellona® Synthetic Undercast Padding

Undercast padding that retains its softness and elasticity to protect exposed bones and nerves under P.o.P. and synthetic cast-

ing tapes. The fine fibrous web is strong enough to permit application in circular turns, but is easy to tear and provides a

smooth surface. The material's self-adhesion facilitates application by simply smoothing down.

Cellona synthetic undercast padding is air-permeable and compensates for changes in temperature. The fibres do not absorb any moisture, the casting tape hardens faster, and the good padding effect is retained with long-term wear.



Cellona® Synthetic Undercast Padding

100% polyester fibres, no optical brighteners, in various sizes in roll form, with paper wrapper; sterilisable (A 121°C; or other methods).

Cellona synthetic undercast padding is also used for padding splints, zinc paste bandages (Unna's Boot) and adhesive bandages. For wound treatment – in combination with textile dressings or ointment tulle - Cellona padding acts as an air-permeable cushion.

Cellona® Synthetic Undercast Padding		
	REF	Shipping Units
6 cm x 3 m	10 686	6 rolls / 60 rolls
10 cm x 3 m	10 687	4 rolls / 40 rolls
40 cm x 10 m	10 697	1 roll / 4 rolls
4 cm x 3 m	10 690	108 rolls
6 cm x 3 m	10 692	72 rolls
10 cm x 3 m	10 694	48 rolls
15 cm x 3 m	10 695	36 rolls
20 cm x 3 m	10 696	24 rolls
5 cm x 2.7 m	92 840	96 rolls
7.5 cm x 2.7 m	92 841	72 rolls
10 cm x 2.7 m	92 842	48 rolls
15 cm x 2.7 m	92 843	36 rolls
20 cm x 2.7 m	92 844	24 rolls

Cotton Padding Padding Bandages

Padding in various combinations of fibres for partial, all-over or circular application under P.o.P. and synthetic casts.

Cotton Padding

Approx. 80% cotton, not degreased, approx. 20% viscose, on a roll, with paper layer.

Padding Bandages

Padding for medical purposes DAB (German Pharmacopoeia), 50% bleached cotton, 50% viscose, in rolls with paper wrapper.



Cotton Padding			
off-white	Size	REF	Shipping Units
individually in a pouch	40 cm x 2.40 m	10 603	1 pack / 20 packs
Padding Bandages			
Padding for medical purposes DAB	10 cm x 5 m 15 cm x 5 m	10 640 10 641	1 box (75 rolls) 1 box (50 rolls)



Cellona® Adhesive Padding Cast Edge Padding

Adhesive, very supple padding material for non-slip application under partially padded P.o.P. and synthetic casting tapes.

Cellona® Adhesive Padding

Preferred for individually cut shapes to protect exposed support zones when applying trunk

or pelvis/leg casts. The polyacrylate adhesive is applied in strips, thus ensuring good air permeability and high skin tolerance.



Cellona® Adhesive Cast Edge Padding

For extra padding at hard edges of casts, especially with synthetic casting tapes; preferred for forearm, upper arm and tutor casts and at all potential pressure points.

Cellona® Adhesive Padding Cellona® Adhesive Cast Edge Padding

Nonwoven padding or air-permeable needled nonwoven: 50% polyester, 30% polypropylene, 20% viscose; polyacrylate adhesive, high skin tolerance.

Cellona® Adhesive Padding			
adhesive, white		REF	Shipping Units
5 mm thick, in a foil pouch	19 cm x 38 cm	50 810	10 pcs. / 60 pcs.
2 mm thick, roll form	58 cm x 1 m 28 cm x 10 m	50 811 50 813	1 box / 10 boxes 1 roll / 11 rolls
Cellona® Adhesive Cast Edge Padding			
white			
2 mm thick, roll form, in a cardboard dispenser	8 cm x 5 m	50 812	1 box / 24 boxes

Cellona® Shoecast • Walking Heel

For physiologically correct movement with a walking cast and for protection against moisture and dirt, L&R offer two walking soles:



Cellona® Shoecast

Walking sole with a high degree of lengthwise and widthwise stability plus maximum ground contact. The anatomical shape permits natural heel-to-toe movement. The highly abrasion-resistant sole profile prevents slipping, while the foamed material cushions every step. A built-up heel cup and cross-over straps with infinitely adjustable hook & loop fastener ensure a secure fit.

Walking sole: 100% foamed polyurethane; straps: 100% polyamide.



Cellona® Walking Heel

Profiled walking heel cast from light grey synthetic rubber; lightweight, anti-slip and shock-absorbing; grooves and lugs for attachment either immediately or after application of P.o.P. or synthetic lower limb casts.



Cellona® Shoecast			
		REF	Shipping Units
individually in a foil pouch			
for children	Size 0, left	50 858	1 pc. / 8 pcs.
	Size 0, right	50 859	1 pc. / 8 pcs.
for adults	Size 1, small left	50 860	1 pc. / 10 pcs.
	Size 1, small right	50 861	1 pc. / 10 pcs.
	Size 2, medium left	50 862	1 pc. / 10 pcs.
	Size 2, medium right	50 863	1 pc. / 10 pcs.
	Size 3, large left	50 864	1 pc. / 8 pcs.
	Size 3, large right	50 865	1 pc. / 8 pcs.
	Size 4, extra large left	50 869	1 pc. / 6 pcs.
	Size 4, extra large right	50 870	1 pc. / 6 pcs.
	Replacement Strap with Hook & Loop Fastener	40 cm (for size 0)	50 868
	90 cm (for sizes 1 and 2)	50 873	10 pcs. / 100 pcs.
	100 cm (for sizes 3 and 4)	50 874	10 pcs. / 100 pcs.
Cellona® Walking Heel			
small		50 820	1 pc. / 20 pcs.
medium		50 821	1 pc. / 20 pcs.
large		50 822	1 pc. / 10 pcs.



Cellona® Shoe

Cellona Shoe is a dressing protection shoe which can be worn over dressings either as a closed shoe or a sandal. When closed, Cellona Shoe offers reliable protection all around: the back of the foot is covered by a softly padded flap which protects foot and toes against cold and wet conditions. If it is placed onto the inside sole before putting on the shoe, Cellona

Shoe is changed into a sandal and there is free air circulation to the foot.

Cellona Shoe provides the bandaged foot with a reliable hold while walking – no matter whether the shoe is worn open or closed. The integrated hook & loop fasteners allow an individual adjustment to the foot. A buckle in the ankle area sup-



ports perfect adjustment and fit. The non-slip special profile allows safe walking with Cellona Shoe even in wet weather conditions.

Cellona Shoe can be worn both on the right and left foot. Its adjustable hook & loop fasteners make it conformable to any dressing – like wound dressings, Plaster of Paris or synthetic casts or to phlebological and lymphological dressings.

Cellona Shoe is water- and dirt-repellent and can be hand-washed at 30°C (no tumble-drying).

Cellona® Shoe			
Size		REF	Shipping Units
for children individually packed in a foil pouch	P (shoe size 25-30)	16 471	1 pc./ 12 pcs.
	XS (shoe size 31-34)	16 472	1 pc./ 12 pcs.
for adults individually packed in a foil pouch	S (shoe size 35-38)	16 473	1 pc./ 12 pcs.
	M (shoe size 39-41)	16 474	1 pc./ 12 pcs.
	L (shoe size 42-43)	16 475	1 pc./ 12 pcs.
	XL (shoe size 44-47)	16 476	1 pc./ 12 pcs.

Cellona® Cream

Protective cream for skin that needs extra care. It is particularly suitable when the skin comes into contact with aqueous alkaline solutions or solutions containing surface-active agents, or with the use of substances which are insoluble in water, e.g. oils and dyes.

This oil-in-water emulsion contains ingredients to restore the oil balance of the skin; it can be easily and evenly spread and is rapidly absorbed. The cream contains a component which



creates a protective barrier against everyday environment influences. The skin is covered with an air-permeable protective film. Cellona Cream can also be used for full body care, even on sensitive skin.



Cellona® Cream

Ingredients:
(INCI)
Aqua (Water)
Glyceryl Stearate SE
Cetearyl Alcohol
Isopropyl Palmitate
Propylene Glycol
Caprylic/Capric Triglyceride
Phenoxyethanol
Methylparaben
Ethylparaben
Propylparaben
Butylparaben
Isobutylparaben
Glyceryl Ricinoleate
Sodium Cetearyl Sulfate
Perfume (Fragrance)

Tube: aluminium
Tube cap: polypropylene
Jar: polypropylene, aluminium foil
Lid: polystyrene



Cellona® Cream

		REF	Shipping Units
in aluminium tube	100 ml	50 814	1 pc. / 12 pcs.
in plastic jar	200 ml	50 805	1 pc. / 20 pcs.