



OR Set Systems

- **OR Drape System**
- **OR Clothing**
- **OR Dressing Materials System**
- **OR Custom Procedure Trays**



Protection & Safety

Prevention of intraoperative and postoperative infection is a top priority in any OR setting. With regard to the products used in the OR, this means that they have to meet the most exacting demands. The only way to provide safe protection for patients and staff is to interfere as fully as possible with all potential routes of transmission. OR products by Lohmann & Rauscher go a long way in this regard.

The Raucodrape sterile OR drape system by L&R frequently meets and even exceeds the requirements stipulated by the European EN 13795 series of standards, in which the basic requirements for OR drapes and

gowns are defined. Relevant specifications include minimum values for resistance against liquid and germ penetration. Raucodrape NEW provides a complete barrier and is an essential factor in effective infection prophylaxis in the OR setting. In combination with Sentinex OR clothing, the Kitpack OR Procedure Tray, and the Setpack OR Dressing Materials, L&R provides a comprehensive range of products for any type of surgical procedure.

Cost-effectiveness & Competence

Our medical product consultants, will help you analyse your specific work processes and procedures and develop individualised, economical solutions.

Monitoring & Quality

All products are specifically designed in function and material to meet applicable individual requirements, and are subject to continuous monitoring for quality in accordance with strict criteria. We attach utmost importance to in-house safety standards, some of which even exceed applicable legal requirements.

Raucodrape®

NEW

Safety and convenience for the OR team and the patient

Innovations in the

Raucodrape® OR Drape System

The new generation of Raucodrape OR drapes focuses on improved materials, innovative service, and a higher level of safety.

Special benefits of the new Raucodrape®

- substantially improved material performance
- wider choice of sets and individual products
- user-friendly folding styles
- self-explanatory pictograms
- double set of documentation labels for each product

The new materials

- are fully absorbent and impermeable to fluids
- meet the requirements of the high-performance standard in terms of EN 13795
- provide enhanced safety of use
- provide increased user convenience
- enable efficient infection prophylaxis in the OR setting
- show excellent draping behaviour

L&R uses the following names for improved materials and a systematically structured range of OR drapes:

Raucodrape®

a 2-layered material for surgical procedures involving a normal amount of fluids

Raucodrape® Plus

a 2-layered material with 2 reinforcing layers for use in procedures involving large amounts of fluids

Raucodrape® Drape Sheets

a 2-layered material for use in minor surgical procedures

Raucodrape® will continue to be available in a 3-layered version as well.





Raucodrape®

Standards defined by indications

From our new Raucodrape range, you choose the drapes that best suit your specific purposes. A large number of single-sheet drapes are available to provide convenient and handy all-in-one-solutions.

Raucodrape® OR standard sets for specific operations

Easier storage and precious time saved in preparing surgical procedures are two strong arguments in favour of specific procedure sets. The Raucodrape product range includes a wide variety of OR standard sets for specific indications.

Solutions for outpatient surgery

For the purposes of outpatient surgery, the L&R portfolio includes an expanded range of Raucodrape products with features tailored to the specific requirements involved.

The comprehensive range of 2-layered OR drapes meets the exacting requirements for convenience and safety. By using Raucodrape OR drapes, process flow can be optimised, and the cost-efficiency of outpatient surgical procedures can be increased.

Kitpack® OR Custom Procedure Tray

Kitpacks include OR drapes as well as all single-use OR components you need for your particular type of surgical procedure. Our indication-specific, process-oriented Kitpacks are thoroughly assembled in accordance with your individual demands and preferences.

Product tests

Product tests provide a sound and solid basis for decision making. L&R invites you to test Raucodrape OR drapes for yourself in a trial run under everyday conditions.

For details and specific enquiries, please feel free to contact our medical device consultants.

Training courses

Topical product training courses are the starting point for the optimised process implementation of Raucodrape OR Drape Systems in your own routine practice. This is why we make it a top priority to provide thorough, systematic training for your OR staff.

You are kindly invited to take advantage of our training course programme. We are looking forward to telling you more about the special benefits and advantages of the Raucodrape OR Drape System.




Raucodrape®

High Performance for your OR

The new Raucodrape OR drapes meet the highest standards in terms of EN 13795.

The safety features of Raucodrape®

- safe barrier function against liquid penetration
- safe barrier function against bacterial penetration
- superior tear resistance - even when wet
- high cleanliness both microbial and particulate
- extremely low linting




Raucodrape® with even higher tear-resistance and good absorption capacity

* PE: polyethylene, PP: polypropylene

Tear resistance


Raucodrape has excellent tear resistance. The improvement achieved in the wet material is outstanding.



all values given in Newton N; tested in accordance with ISO 9073-3
*as compared to the similar previous version of the material

Water absorption capacity

Raucodrape can absorb 250 ml of water per m².



tested in accordance with ISO 9073-6

Liquid penetration

Raucodrape is liquid-proof. Even under a water column of 125 cm, it still forms an impenetrable barrier.



tested in accordance with EN 20811

Raucodrape®


A “plus” in High Performance Standard

Raucodrape® Plus unfolds its strength right where it is needed most. The field of operation is framed by 2 additional layers of reinforcing material.

With a total of 4 layers in the area of maximum stress, the material meets the most exacting demands for performance in liquid-intensive surgical procedures. Raucodrape® Plus is far above the highest standard in terms of EN 13795.

Additional safety features of Raucodrape® Plus:

- reinforced field of operation with 4 layers of material
- increased tear resistance - both dry and wet
- excellent absorption capacity
- enhanced safety in surgical procedures involving large amounts of liquids



Raucodrape® Plus with 4 layers for increased safety under highest demands

* PE: polyethylene, PP: polypropylene

Tear resistance


Raucodrape Plus is strong - and right on target. Test results have confirmed its suitability or the highest demands in the OR.



all values given in Newton N; tested in accordance with ISO 9073-3
*as compared to the similar previous 3-layered version

Water absorption capacity


Raucodrape Plus absorbs 520 ml of water per m². Absorption takes place over the entire surface.



tested in accordance with ISO 9073-6

Liquid penetration

Raucodrape Plus is particularly liquid-proof. Even under a water column in excess of 150 cm, it still forms an impenetrable barrier.



tested in accordance with EN 20811

Raucodrape®

Strong solutions, right down to the last detail - that is what the new product range is all about

Standards are defined by actual needs in daily practice

The strong points of Raucodrape are apparent in the smallest details - a multitude of innovative design features have been added for easier unfolding of the drape sheets, a wider variety of applications, and enhanced safety and convenience in surgical procedures.

Strong solutions in detail

■ Raucodrape® makes unfolding faster and easier

Raucodrape OP drapes now feature optimised folding styles. Depending on the specific purposes for which it is used, each type of drape has its own unique folding style which ensures safety in use and makes the process of unfolding faster and easier.

Now all our instrument table covers come with an extra-wide folding style for better handling and application.

In all self-adhesive fenestrated drapes, the aperture is visible even before the sheet is unfolded. This makes it considerably easier to apply the product correctly to the patient

■ Self-explanatory pictograms

Safe and fast unfolding of the drapes is also facilitated by the easily comprehensible pictograms. At a glance, they show the first step to unfolding the sheet, the direction in which it unfolds the correct positioning of the drape on the patient, as well as any splits, integrated leggings, and lots more.

■ Expanded range of single-sheet drape systems

The Raucodrape product range provides a comprehensive choice of modern single-sheet drape systems. The special features integrated in the single-sheet drapes save considerable time and effort in preparing for a surgical procedure.

As a general rule, all single-sheet drapes have a three-dimensional cut in the area of intervention, which prevents any undesirable tension and tautness in the drape. Thus, there is always enough space, even when several surgeons are working on a patient.

All integrated leggings in the single-sheet drapes are extra-long (125 cm), i.e. they can be used even with padded leg rests.

■ New integrated special features

L&R has placed increased focus on adding special features to


the Raucodrape OR drapes. The portfolio now includes even more of the following product categories:

- fluid collection pouches
- incision films
- tube retainers
- instrument pouches

■ Three-piece adhesive strips

Adhesive strips are an essential factor in the overall efficiency of any process, in particular with regard to product and user safety:

- All split sheet have an ergonomically designed, three-piece adhesive strip, which makes them easier to apply.
- Many adhesive sheets feature an extra-long adhesive strip to accommodate adipose patients.
- All adhesive strips are equipped with a 2 cm finger-lift for a better and safer grip.



A small detail that makes work so much easier: the 2 cm finger-lifts on all adhesive strips








Raucodrape® OR Drape System

General Surgery

Raucodrape® Basic-Set Mini, 2-layered, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 300	10 sets / 20 sets
4 adhesive drapes, 2-layered	75 x <u>90</u> cm		
Raucodrape® Set for Minor Surgery, sterile			
1 table cover	150 x 190 cm	33 315	22 sets
2 towels	30 x 40 cm		
1 adhesive tape	9 x 50 cm		
1 fenestrated drape, 2-layered, selfadhesive, circular aperture (centered)	225 x 240 cm		
	12 cm		
Raucodrape® Basic-Set, 2-layered, sterile			
1 table cover	150 x 190 cm	33 301	13 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
2 adhesive drapes, 2-layered	75 x <u>90</u> cm		
1 adhesive drape, 2-layered adhesive strip	175 x <u>180</u> cm 80 cm		
1 adhesive drape, 2-layered adhesive strips	150 x <u>240</u> cm 15 + 70 + 15 cm		
Raucodrape® Universal-Set II, 3-layered, sterile			
1 table cover	150 x 190 cm	33 305	11 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
2 adhesive drapes, 3-layered	75 x <u>90</u> cm		
1 adhesive drape, 3-layered adhesive strip	175 x <u>180</u> cm 80 cm		
1 adhesive drape, 3-layered adhesive strips	150 x <u>240</u> cm 15 + 70 + 15 cm		
Raucodrape® Universal-Set II Special, 3-layered, sterile			
1 table cover	150 x 190 cm	33 306	9 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
2 adhesive drapes, 3-layered	75 x <u>90</u> cm		
1 adhesive drape, 3-layered adhesive strip	175 x <u>180</u> cm 80 cm		
1 adhesive drape, 3-layered adhesive strips	175 x 300 cm 20 + 80 + 20 cm		

— adhesive side




Raucodrape® Universal-Set I, 2-layered Plus, sterile

Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 307	13 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive strip	9 x 50 cm		
2 adhesive drapes, 2-layered Plus	75 x 90 cm		
1 adhesive drape, 2-layered Plus, adhesive strip	175 x 180 cm 80 cm		
1 adhesive drape, 2-layered Plus, adhesive strips	150 x 250 cm 15 + 70 + 15 cm		

— adhesive side

Raucodrape® OR Drape System General Surgery



Raucodrape® Universal-Set I Super, 2-layered Plus, sterile

Contents	Size	REF	Shipping Units		
1 table cover	150 x 190 cm	33 308	11 sets		
4 towels	30 x 40 cm				
1 instrument table cover	80 x 145 cm				
1 adhesive tape	9 x 50 cm				
2 adhesive drapes, 2-layered Plus	90 x <u>110</u> cm				
1 adhesive drape, 2-layered Plus, adhesive strip	190 x 225 cm 15 + 70 + 15 cm				
1 anaesthesia drape, 2-layered Plus, split, with armrest cover	175 x <u>270</u> cm 45 x <u>65</u> cm				
adhesive side					

Raucodrape® Universal Split Sheet Set I, sterile

Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 310	12 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 adhesive drape, 2-layered	75 x <u>90</u> cm		
1 split sheet, 2-layered Plus, split	<u>225</u> x <u>260</u> cm 10 x 100 cm		
1 adhesive drape, 2-layered adhesive strips	150 x <u>240</u> cm 15 + 70 + 15 cm		

Raucodrape® Universal Split Sheet Set II, sterile

1 table cover	150 x 190 cm	33 311	12 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 split sheet, 2-layered Plus split, adhesive strip	<u>225</u> x 300 cm 10 x 100 cm 70 cm		
1 anaesthesia drape, 2-layered Plus, split	175 x <u>270</u> cm		
with armrest cover	45 x <u>65</u> cm		

Raucodrape® Universal Split Sheet Set III, sterile

1 table cover	150 x 190 cm	33 312	12 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 split sheet, 2-layered Plus, split, with integrated tube retainers	<u>225</u> x 300 cm 30 x 80 cm		
1 anaesthesia drape, 2-layered Plus, split	175 x <u>270</u> cm		
with armrest cover	45 x <u>65</u> cm		


Raucodrape® Variable Split Sheet Set, sterile

1 table cover	150 x 190 cm	33 313	15 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 split sheet, 2-layered Plus, split	<u>225</u> x 240 cm 20 x 20 cm		
1 split sheet, 2-layered, split	150 x <u>260</u> cm 20 x <u>60</u> cm		

— adhesive side



Raucodrape® OR Drape System General Surgery




Raucodrape® LAP-Set I, sterile

Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm		
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 LAP drape I, for use in a level position, abdominal aperture, with armrest cover, with integrated tube retainers	260/200 x 335 cm 28 x 32 cm	33 320	14 sets


Raucodrape® LAP-Set II, sterile			
Contents	Size	REF	Shipping Units
1 table cover 2 towels 1 instrument table cover 1 LAP drape II, with integrated leggings, abdominal aperture, with armrest cover, with integrated tube retainers	150 x 190 cm 30 x 40 cm 80 x 145 cm 260 x 310 cm 125 cm 28 x 32 cm	33 321	12 sets
Raucodrape® LAP-Set III, sterile			
1 table cover 2 towels 1 instrument table cover 1 LAP drape III, with integrated leggings, abdominal aperture, with fenestrated incision film, with 2-chamber instrument pouches with armrest cover, with integrated tube retainers	150 x 190 cm 30 x 40 cm 80 x 145 cm 260 x 310 cm 125 cm 35 x 37 cm 28 x 32 cm	33 322	11 sets
Raucodrape® Rectum-Set I, sterile			
1 table cover 4 towels 1 instrument table cover 1 Rectum drape I, with integrated leggings, abdominal aperture, perineal aperture, with armrest cover, with integrated tube retainers	150 x 190 cm 30 x 40 cm 80 x 145 cm 260 x 310 cm 125 cm 28 x 32 cm 10 x 15 cm	33 325	11 sets
Raucodrape® Rectum-Set II, sterile			
1 table cover 4 towels 1 instrument table cover 1 bottom sheet 1 rectum drape II, with integrated leggings, abdominal aperture, with fenestrated incision film, perineal aperture, with 2-chamber instrument pouches, with armrest cover, with integrated tube retainers	150 x 190 cm 30 x 40 cm 80 x 145 cm 75 x 120 cm 260 x 310 cm 125 cm 35 x 37 cm 28 x 32 cm 10 x 15 cm	33 326	11 sets



Raucodrape® OR Drape System Paediatrics



REF.: 33 330



Raucodrape® Paediatric-Set II, 2-layered, sterile

Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 330	14 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape,	9 x 50 cm		
2 adhesive drapes, 2-layered	75 x 90 cm		
1 adhesive drape, 2-layered	175 x 180 cm		
adhesive strip	80 cm		
1 adhesive drape, 2-layered	150 x 240 cm		
adhesive strip	80 cm		

Raucodrape® Fenestrated Drape with Incision Film, sterile

Fenestrated drape, 2-layered, aperture (lengthwise, centered), with incision film	120 x 150 cm 7 x 10 cm	33 055	22 pcs. / 44 pcs.
---	---------------------------	--------	-------------------

— adhesive side

Raucodrape® OR Drape System

Gynaecology/Urology

Raucodrape® Pfannenstiel-Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 360	12 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 Pfannenstiel Drape	260 x 310 cm		
with integrated leggings,	125 cm		
abdominal aperture,	35 x 37 cm		
with fenestrated incision film,	24 x 29 cm		
with armrest cover,			
with integrated tube retainers			
Raucodrape® Rectum-Set I, sterile			
1 table cover	150 x 190 cm	33 325	11 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 rectum drape I,	260 x 310 cm		
with integrated leggings,	125 cm		
abdominal aperture,	28 x 32 cm		
perineal aperture,	10 x 15 cm		
with armrest cover,			
with integrated tube retainers			
Raucodrape® Rectum-Set II, sterile			
1 table cover	150 x 190 cm	33 326	11 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 bottom sheet	75 x 120 cm		
1 rectum drape II,	260 x 310 cm		
with integrated leggings,	125 cm		
abdominal aperture,	35 x 37 cm		
with fenestrated incision film,	28 x 32 cm		
perineal aperture,	10 x 15 cm		
with 2-chamber instrument			
pouches			
with armrest cover,			
with integrated tube retainers			



Raucodrape® OR Drape System

Gynaecology/Urology

Raucodrape® C-section Set I, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm		
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 OR drape sheet soft, fluid-repellent nonwoven material	90 x 110 cm		
1 C-section drape I for use in a level position, abdominal aperture, fenestrated incision film integrated fluid collection pouch with exit port, with armrest cover	260/200 x 335 cm 35 x 37 cm 24 x 29 cm	33 361	11 sets
Raucodrape® C-section Drape, sterile			
C section drape, for use in a level position, abdominal aperture, with fenestrated incision film, integrated fluid collection pouch with exit port, with armrest cover	260/200 x 335 cm 35 x 37 cm 24 x 29 cm	33 359	8 pcs. / 16 pcs.
Raucodrape® OR Drape Sheet Soft, sterile			
OR drape sheet soft fluid-repellent nonwoven material	90 x 110 cm	33 006	20 pcs. / 80 pcs.
Raucodrape® C-section Set II, sterile			
1 table cover	150 x 190 cm		
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 OR drape sheet soft, fluid-repellent nonwoven material	90 x 110 cm		
1 C-section drape II, with integrated leggings, abdominal aperture, with fenestrated incision film, integrated fluid collection pouch with exit port, with armrest cover,	260 x 310 cm 125 cm 35 x 37 cm 24 x 29 cm	33 362	11 sets

Raucodrape® Gynaecology/Cystoscopy Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover 2 leg protectors 1 gynaecology/cystoscopy drape, perineal aperture (centered)	120 x 140 cm 75 x 120 cm 90 x 175 cm 9 x 15 cm	33 363	8 sets / 16 sets
Raucodrape® Gynaecology/Cystoscopy Drape, sterile			
1 gynaecology/cystoscopy drape, perineal aperture (centered)	90 x 175 cm 9 x 15 cm	33 039	35 pcs. / 70 pcs.
Raucodrape® Lithotomy-Set, sterile			
1 table cover 2 towels 1 adhesive tape 1 bottom drape, 2-layered 1 lithotomy drape, with integrated leggings, perineal aperture	150 x 190 cm 30 x 40 cm 9 x 50 cm 75 x 120 cm 175/260 x 210 cm 125 cm 10 x 15 cm	33 366	21 sets
Raucodrape® Gynaecology Set I, sterile			
1 table cover 2 towels 1 instrument table cover 1 hook and loop fastener, 1 adhesive drape (PE film) 1 gynaecology drape I, with integrated leggings, perineal aperture,	150 x 190 cm 30 x 40 cm 80 x 145 cm 2.5 x 25 cm 50 x 50 cm 240/260 x 230 cm 125 cm 10 x 15 cm	33 364	14 sets
Raucodrape® Gynaecology Set II, sterile			
1 table cover 2 towels 1 instrument table cover 1 adhesive drape (PE film) 1 gynaecology drape II, with integrated leggings, perineal aperture, framed with Raucodrape Plus, integrated fluid collection pouch with sieve and exit port	150 x 190 cm 30 x 40 cm 80 x 145 cm 9 x 50 cm 50 x 50 cm 240/260 x 230 cm 125 cm 10 x 15 cm	33 365	14 sets
— adhesive side			



Raucodrape® OR Drape System

Gynaecology/Urology

Raucodrape® T.U.R. Set I, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 367	23 sets
2 towels	30 x 40 cm		
1 hook and loop fastener	2.5 x 25 cm		
1 T.U.R. drape I,	175/260 x 210 cm		
with integrated leggings,	125 cm		
perineal aperture,	5 cm		
suprapubic aperture,	8 cm		
integrated fluid collection			
pouch with sieve an exit			
port,			
with a latex-free finger stall			
Raucodrape® TUR Set II, sterile			
1 table cover	150 x 190 cm	33 368	18 sets
2 towels	30 x 40 cm		
1 hook and loop fastener	2.5 x 25 cm		
1 T.U.R. drape II,	175/260 x 210 cm		
with integrated leggings,	125 cm		
perineal aperture,	5 cm		
suprapubic aperture,	8 cm		
with integrated sieve,			
with a latex-free finger stall			

Raucodrape® OR Drape System

Orthopaedics/Traumatology

Raucodrape® Universal Split Sheet Set I, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 310	12 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 adhesive drape, 2-layered	75 x 90 cm		
1 split sheet, 2-layered Plus, split,	<u>225</u> x 260 cm 10 x 100 cm		
1 adhesive drape, 2-layered, adhesive strip	<u>150</u> x 240 cm 15 + 70 + 15 cm		
Raucodrape® Hip-Set I, sterile			
1 table cover	150 x 190 cm	33 340	8 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 adhesive drape, 2-layered	75 x 90 cm		
OR drape sheet, 2-layered	150 x 180 cm		
1 split sheet, 2-layered Plus, split	<u>225</u> x 260 cm 10 x 100 cm		
1 anaesthesia drape, 2-layered Plus, split with armrest cover	175 x <u>270</u> cm 45 x <u>65</u> cm		
Raucodrape® Hip Set II, sterile			
1 table cover	150 x 190 cm	33 341	8 sets
4 towels	30 x 40 cm		
1 instrument table cover Special	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 adhesive drape, 2-layered	75 x 90 cm		
1 OR drape sheet, 2-layered	150 x 180 cm		
1 orthopaedic stockinette,	33 x 110 cm		
2 adhesive film strips	10 x 50 cm		
1 split sheet, 2-layered Plus, split with integrated tube retainers	<u>225</u> x 280 cm 10 x 100 cm		
1 anaesthesia drape, 2-layered Plus, split with armrest cover	175 x <u>270</u> cm 45 x <u>65</u> cm		
— adhesive side			



Raucodrape® OR Drape System

Orthopaedics/Traumatology

Raucodrape® Extremity Set I, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 343	11 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive drape, 2-layered	75 x 75 cm		
1 OR drape sheet, 2-layered	150 x 180 cm		
1 extremity drape, 2-layered Plus, elastic aperture, with integrated tube retainers	225 x 320 cm 7 cm		
Raucodrape® Extremity-Set II, sterile			
1 table cover	150 x 190 cm	33 344	11 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 OR drape sheet, 2-layered	150 x 180 cm		
1 orthopaedic, stockinette	33 x 55 cm		
1 adhesive film strip	10 x 50 cm		
1 extremity drape, 2-layered Plus, elastic aperture, with integrated tube retainers	225 x 320 cm 7 cm		
Raucodrape® Extremity Set III, sterile			
1 table cover	150 x 190 cm	33 345	10 sets
2 towels	30 x 40 cm		
1 instrument table cover Special	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 OR drape sheet, 2-layered	150 x 180 cm		
1 extremity drape, 2-layered Plus, elastic aperture, with armrest cover, with integrated tube retainers	270/200 x 335 cm 7 cm		
Raucodrape® Extremity Drapes, sterile			
Extremity drape, 2-layered, elastic aperture (centered)	200 x 300 cm 7 cm	33 085	13 pcs. / 26 pcs.
Extremity drape, 2-layered, elastic aperture (centered)	200 x 300 cm 7 cm	33 083	12 pcs. / 24 pcs.
Extremity drape, 2-layered Plus, elastic aperture (centered), with integrated tube retainers	200 x 300 cm 7 cm	33 084	11 pcs. / 22 pcs.
— adhesive side			

Raucodrape® Extremity Set Hand/Foot, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm		
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 extremity drape Hand/Foot 2-layered Plus, elastic aperture, with integrated tube retainers	225 x 300 cm 4 cm	33 347	17 sets
Raucodrape® Hand Surgery Set, sterile			
1 table cover	150 x 190 cm		
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 extremity drape Hand/Foot 2-layered Plus, elastic aperture, with integrated tube retainers	225 x 300 cm 4 cm	33 349	14 sets
Raucodrape® Extremity Drapes Hand/Foot, sterile			
Extremity drape, 2-layered, elastic aperture (centered)	100 x 175 cm 4 cm	33 080	20 pcs. / 40 pcs.
Extremity drape, 2-layered, elastic aperture (centered)	150 x 180 cm 4 cm	33 081	18 pcs. / 36 pcs.
Extremity drape, 2-layered Plus, elastic aperture (off centered), with integrated tube retainers	225 x 300 cm 4 cm	33 082	10 pcs. / 20 pcs.
Raucodrape® Arthroscopy Set I, sterile			
1 wrapping paper	80 x 100 cm		
1 instrument table cover Lite	80 x 145 cm		
1 stockinette, short	24 x 80 cm		
2 adhesive tapes	9 x 50 cm		
1 extremity drape, 2-layered, elastic aperture	200 x 300 cm 7 cm	33 350	23 sets




Raucodrape® OR Drape System

Orthopaedics/Traumatology


Raucodrape® Arthroscopy Set II, sterile			
Contents	Size	REF	Shipping Units
1 wrapping paper	80 x 100 cm	33 351	21 sets
1 instrument table cover Lite	80 x 145 cm		
1 camera drape (reverse-fold)	13 x 250 cm		
1 orthopaedic stockinette,	33 x 55 cm		
1 adhesive film strip	10 x 50 cm		
1 adhesive tape	9 x 50 cm		
1 extremity drape, 2-layered, elastic aperture	200 x 300 cm 7 cm		
Raucodrape® Arthroscopy Set III, sterile			
1 wrapping paper	80 x 100 cm	33 352	11 sets
2 SMMS surgical gowns Standard	150 cm		
1 instrument table cover Lite	80 x 145 cm		
1 camera drape (reverse-fold)	13 x 250 cm		
1 orthopaedic stockinette,	33 x 55 cm		
1 adhesive film strip	10 x 50 cm		
1 adhesive tape	9 x 50 cm		
1 extremity drape, 2-layered, elastic aperture	200 x 300 cm 7 cm		
Raucodrape® Arthroscopy Set III, sterile			
1 wrapping paper	80 x 100 cm	33 353	15 sets
1 SMMS surgical gown Standard	150 cm		
1 instrument table cover Lite	80 x 145 cm		
1 camera drape (reverse-fold)	13 x 250 cm		
1 orthopaedic stockinette,	33 x 55 cm		
1 adhesive film strip	10 x 50 cm		
1 adhesive tape	9 x 50 cm		
1 extremity drape, 2-layered, elastic aperture	200 x 300 cm 7 cm		
Raucodrape® Shoulder Arthroscopy Set I, sterile			
1 wrapping paper	80 x 100 cm	33 354	20 sets
1 instrument table cover Lite	80 x 145 cm		
1 orthopaedic stockinette	33 x 55 cm		
1 adhesive film strip	10 x 50 cm		
1 split sheet, 2-layered, split	225 x 260 cm 10 x 100 cm		

— adhesive side

Raucodrape® Knee Arthroscopy Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 348	10 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 OR drape sheet, 2-layered	150 x 180 cm		
1 orthopaedic stockinette,	33 x 55 cm		
1 adhesive film strip	10 x 50 cm		
1 knee arthroscopy drape, elastic apertures, integrated fluid collection pouch with exit port, integrated tube retainers	225 x 320 cm 6 cm + 7 cm		
Raucodrape® Knee Arthroscopy Drape, sterile			
Knee arthroscopy drape, elastic apertures, integrated fluid collection pouch with exit port, integrated tube retainers	225 x 320 cm 6 cm + 7 cm	33 086	11 pcs. / 22 pcs.
Raucodrape® Shoulder Arthroscopy Set Beach Chair, sterile			
1 table cover	150 x 190 cm	33 355	14 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 stockinette, short	24 x 80 cm		
1 shoulder arthroscopy drape, elastic, pear-shaped aperture, integrated fluid collection pouch with exit port, with 3 hook and loop fasteners	225 x 380 cm 11.5 x 12.5 cm		
Raucodrape® Vertical Drape Set I, sterile			
1 table cover	150 x 190 cm	33 342	21 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 vertical drape (transparent PE film)	250 x <u>330</u> cm		
with integrated incision film, integrated fluid collection pouch with sieve and exit port, with 2 instrument pouches, with 2 hook and loop fasteners for retaining tubes and cords	80 x 30 cm		
 adhesive side			



Raucodrape® OR Drape System Orthopaedics/Traumatology



Raucodrape® Vertical Drape, sterile

Contents	Size	REF	Shipping Units
1 vertical drape (transparent PE film), with integrated incision film, integrated fluid collection pouch with sieve and exit port, with 2 instrument pouches, with 2 hook and loop fasteners for retaining tubes and cords	250 x 330 cm 80 x 30 cm	33 038	13 pcs. / 26 pcs.
— adhesive side			

Raucodrape® OR Drape System

Neurosurgery/ENT/CMF/ Ophthalmology

Raucodrape® Kranio-Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 370	14 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
3 adhesive drapes (PE film)	50 x 50 cm		
1 craniotomy drape, aperture,	225 x 280 cm		
with integrated incision film, integrated fluid collection pouch with lateral gutters, sieve, and exit port, with integrated tube retainers	19 x 25 cm		
Raucodrape® LAP/Spinal Drape, sterile			
LAP/Spinal drape, 3-layered aperture (oval, lengthwise, off-centre, self-adhesive)	225 x 320 cm 11 x 22 cm	33 051	8 pcs. / 16 pcs.
Raucodrape® Head Set I, sterile			
1 table cover	150 x 190 cm	33 371	16 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 OR drape sheet, 2-layered	75 x 90 cm		
1 split sheet, 2-layered, split	150 x 250 cm 10 x 80 cm		
1 adhesive drape, 2-layered	75 x 90 cm		
Raucodrape® Head-Set II, sterile			
1 table cover	150 x 190 cm	33 372	13 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 turban drape (2-piece) fluid-repellent nonwoven material,	130 x 100 cm 225 x 300 cm		
2-layered,	14 x 80 cm		
1 split sheet, 2-layered Plus, split (rectangular), with integrated tube retainers			
— adhesive side			



Raucodrape® OR Drape System

Neurosurgery/ENT/CMF/ Ophthalmology

Raucodrape® Variable Split Sheet Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 313	15 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 split sheet, 2-layered Plus, split	225 x 240 cm 20 x 20 cm		
1 split sheet, 2-layered Plus, split	150 x 260 cm 20 x 60 cm		
Raucodrape® Split Sheet, sterile			
Split Sheet, 2-layered split	150 x 250 cm 10 x 80 cm	33 073	14 pcs. / 28 pcs.
Raucodrape® Fenestrated Face Drape, sterile			
Fenestrated face drape, 2-layered, aperture, with integrated incision film, (oval, widthwise, off-centre)	150 x 240 cm 8 x 19 cm	33 054	15 pcs. / 30 pcs.
Raucodrape® Fenestrated Ear Drape, sterile			
Fenestrated ear drape, 2-layered, aperture (oval, lengthwise, off-centre, self-adhesive)	175 x 240 cm 5 x 7.5 cm	33 050	11 pcs. / 22 pcs.
Raucodrape® Eye Drape, sterile			
Eye drape (fluid-repellent nonwoven material), central aperture, with integrated incision film, integrated fluid collection pouch	150 x 150 cm 7 x 10 cm	33 056	17 pcs. / 34 pcs.
— adhesive side			

Raucodrape® OR Drape System

Cardiology/Angiography

Raucodrape® Cardiovascular Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 380	7 sets
4 towels	30 x 40 cm		
2 instrument table cover	80 x 145 cm		
2 adhesive tapes	9 x 50 cm		
2 foot protectors	37.5 x 24 cm		
1 OR drape sheet	150 x 180 cm		
1 split sheet bilateral	<u>175</u> x 260 cm		
2-layered Plus, split, with integrated tube retainer	20 x 100 cm		
1 cardiovascular drape	260/200 x 335 cm		
2-/3-layered, aperture, with integrated incision film, with 3-chamber instrument pouches on both sides, with integrated tube retainer and a cord organiser	33 x 38 cm		
Raucodrape® Angiography Drape I, sterile			
1 wrapping paper	80 x 100 cm	33 381	11 sets / 22 sets
4 towel	30 x 40 cm		
1 adhesive tape	9 x 25 cm		
1 angiography drape, 3-layered	150 x 240 cm		
2 apertures (circular, off-centre)	12 cm		
Raucodrape® Angiography Drape II, sterile			
1 wrapping paper	80 x 100 cm	33 382	8 sets / 16 sets
1 towel	30 x 40 cm		
1 adhesive tape	9 x 25 cm		
1 angiography drape, 3-layered, PE film on both sides of drape,	227 x 330 cm		
2 apertures (off-centre) with fenestrated incision film	12 x 12 cm 8 cm		
— adhesive side			



Raucodrape® OR Drape System

Drape Sheets and Accessories

Raucodrape® Drape Sheets, sterile			
Contents	Size	REF	Shipping Units
Drape sheet, 2-layered	38 x 45 cm	33 001	220 pcs. / 880 pcs.
Drape sheet, 2-layered	45 x 75 cm	33 002	120 pcs. / 480 pcs.
Drape sheet, 2-layered	75 x 90 cm	33 003	65 pcs. / 260 pcs.
Raucodrape® OR Drape Sheets, sterile			
OR Drape sheet, 2-layered	75 x 90 cm	33 005	40 pcs. / 160 pcs.
OR Drape sheet, 2-layered	100 x 150 cm	33 008	27 pcs. / 54 pcs.
OR Drape sheet, 2-layered	150 x 180 cm	33 009	13 pcs. / 26 pcs.
OR Drape sheet, 2-layered	150 x 240 cm	33 004	13 pcs. / 26 pcs.
OR Drape sheet soft (fluid-repellent nonwoven material)	90 x 110 cm	33 006	20 pcs. / 80 pcs.
Raucodrape® OR Table Covers, sterile			
OR Table cover	120 x 140 cm	33 090	26 pcs. / 52 pcs.
OR Table cover	150 x 190 cm	33 091	15 pcs. / 30 pcs.
OR Table cover Special	150 x 190 cm	33 092	11 pcs. / 22 pcs.
Instrument Table cover Standard	80 x 145 cm	33 093	35 pcs. / 70 pcs.
Instrument Table cover Lite	80 x 145 cm	33 094	40 pcs. / 80 pcs.
Instrument Table cover Special	80 x 145 cm	33 095	25 pcs. / 50 pcs.
Raucodrape® OR Adhesive Drapes, sterile			
Adhesive drape, 2-layered	50 x <u>50</u> cm	33 010	80 pcs. / 320 pcs.
Adhesive drape, 2-layered	50 x <u>75</u> cm	33 011	60 pcs. / 240 pcs.
Adhesive drape, 2-layered	75 x <u>75</u> cm	33 012	40 pcs. / 160 pcs.
Adhesive drape, 2-layered	75 x <u>90</u> cm	33 013	38 pcs. / 152 pcs.
Adhesive drape, 2-layered	90 x <u>110</u> cm	33 020	25 pcs. / 100 pcs.
Adhesive drape, 2-layered	150 x <u>180</u> cm	33 025	20 pcs. / 40 pcs.
Adhesive drape, 2-layered Plus	150 x <u>180</u> cm	33 015	18 pcs. / 36 pcs.
Adhesive drape, 2-layered	175 x <u>180</u> cm	33 016	17 pcs. / 34 pcs.
Adhesive drape, 2-layered Plus	175 x <u>180</u> cm	33 017	15 pcs. / 30 pcs.
Adhesive drape, 2-layered	150 x <u>240</u> cm	33 018	15 pcs. / 30 pcs.
Adhesive drape, 2-layered Plus	150 x <u>250</u> cm	33 019	14 pcs. / 28 pcs.
Adhesive drape, 2-layered Plus	<u>190</u> x 225 cm	33 021	11 pcs. / 22 pcs.
Adhesive drape, 3-layered	175 x <u>300</u> cm	33 023	6 pcs. / 12 pcs.
Instrument Collection Drape, 2-layered	<u>90</u> x 150 cm	33 022	19 pcs. / 76 pcs.
— adhesive side			

Raucodrape® OR Fenestrated Drapes, sterile				
	Size/position of aperture	Size	REF	Shipping Units
Fenestrated drape, 2-layered, self-adhesive	variable	50 x 75 cm	33 058	50 pcs. / 200 pcs.
Fenestrated drape, 2-layered, self-adhesive	variable	75 x 90 cm	33 059	30 pcs. / 120 pcs.
Fenestrated drape, 2-layered	6 cm, circular, centered	45 x 75 cm	33 040	130 pcs. / 520 pcs.
Fenestrated drape, 2-layered, self-adhesive	6 cm, circular, centered	50 x 60 cm	33 041	70 pcs. / 280 pcs.
Fenestrated drape, 2-layered	8 cm, circular, centered	75 x 90 cm	33 042	70 pcs. / 280 pcs.
Fenestrated drape, 2-layered, self-adhesive	8 cm, circular, centered	75 x 90 cm	33 043	40 pcs. / 160 pcs.
Fenestrated drape, 2-layered, self-adhesive	7 x 10 cm, oval, widthwise, off-centre	120 x 150 cm	33 044	25 pcs. / 50 pcs.
Fenestrated drape, 2-layered, self-adhesive	12 cm, circular, off-centre	150 x 200 cm	33 046	18 pcs. / 36 pcs.
Raucodrape® OR Fenestrated Drapes, sterile				
Epidural drape, 2-layered, self-adhesive	10 cm, circular, off-centre	<u>60</u> x 75 cm	33 047	68 pcs. / 272 pcs.
Epidural drape, 2-layered, self-adhesive	10 x 15 cm, rectangular, lengthwise, off-centre	<u>55</u> x 90 cm	33 049	45 pcs. / 180 pcs.
Fenestrated Ear drape, 2-layered, self-adhesive	5 x 7.5 cm, oval, lengthwise, off-centre	175 x 240 cm	33 050	11 pcs. / 22 pcs.
Fenestrated drape with incision film, 2-layered	7 x 10 cm, oval, lengthwise, centered	120 x 150 cm	33 055	22 pcs. / 44 pcs.
Fenestrated Face drape, 2-layered, with incision film	8 x 19 cm, oval, widthwise, off-centre	150 x 240 cm	33 054	15 pcs. / 30 pcs.
LAP/Spinal drape, 3-layered, self-adhesive	11 x 22 cm, oval, lengthwise, off-centre	225 x 320 cm	33 051	8 pcs. / 16 pcs.
— adhesive side				



Raucodrape® OR Drape System

Drape Sheets and Accessories

Raucodrape® OR Split Sheets, sterile				
	Split Size	Size	REF	Shipping Units
Split Sheet, 2-layered	7 x 30 cm	75 x <u>90</u> cm	33 062	35 pcs. / 140 pcs.
Split Sheet, 2-layered	10 x 45 cm	125 x <u>150</u> cm	33 064	28 pcs. / 56 pcs.
Split Sheet, (PE film)	14 x 50 cm	120 x <u>150</u> cm	33 065	60 pcs. / 120 pcs.
Split Sheet, (PE film)	14 x 75 cm	<u>150</u> x 210 cm	33 066	40 pcs. / 80 pcs.
Split Sheet, 2-layered	10 x 65 cm	150 x <u>200</u> cm	33 067	14 pcs. / 28 pcs.
Split Sheet, 3-layered	10 x 65 cm	150 x <u>200</u> cm	33 068	10 pcs. / 20 pcs.
Split Sheet, 2-layered	10 x 80 cm	<u>150</u> x 250 cm	33 073	14 pcs. / 28 pcs.
Split Sheet, 2-layered	10 x 100 cm	<u>225</u> x 260 cm	33 069	13 pcs. / 26 pcs.
Split Sheet, 2-layered Plus	10 x 100 cm	<u>225</u> x 260 cm	33 071	11 pcs. / 22 pcs.
Split Sheet, 2-layered Plus, with integrated tube retainers	10 x 100 cm	<u>225</u> x 280 cm	33 060	10 pcs. / 20 pcs.
Split sheet bilateral, 2-layered Plus, with integrated tube retainers	20 x 100 cm	<u>175</u> x 260 cm	33 072	13 pcs. / 26 pcs.

— adhesive side

Raucodrape® Accessories, sterile			
	Size	REF	Shipping Units
Arm Protector	50 cm	33 103	45 pcs. / 180 pcs.
Leg Protector (2 pcs.)	50 x 110 cm	33 104	26 pcs. / 104 pcs.
Leg Protector (2 pcs.)	75 x 120 cm	33 105	30 pcs. / 60 pcs.
Orthopaedic Stockinette, with 1 adhesive film strip	33 x 55 cm	33 106	65 pcs. / 130 pcs.
Orthopaedic Stockinette, with 2 adhesive film strips	33 x 110 cm	33 107	38 pcs. / 76 pcs.
Stockinette, short	24 x 80 cm	18 302	14 pcs. / 56 pcs.
Stockinette, small	18 x 50 cm	18 303	25 pcs. / 100 pcs.
Stockinette, medium	24 x 110 cm	18 304	12 pcs. / 48 pcs.
Stockinette, large	32 x 120 cm	18 305	9 pcs. / 36 pcs.
Towels (in pairs)	30 x 40 cm	33 108	120 pcs. / 480 pcs.
Bottom Drape I	75 x 120 cm	33 101	28 pcs. / 112 pcs.
Bottom Drape II with collection pouch	75 x 100 cm	33 102	35 pcs. / 70 pcs.

Raucodrape® Accessories, sterile

	Size	REF	Shipping Units
Absorption Pad (2 pcs.)	35 x <u>50</u> cm	33 125	80 pcs. / 320 pcs.
Absorption Pad	50 x <u>75</u> cm	33 126	33 pcs. / 132 pcs.
Adhesive Tape	9 x <u>50</u> cm	33 116	250 pcs. / 750 pcs.
Hook and Loop Fastener	2.5 x 25 cm	17 746	130 pcs. / 520 pcs.
Instrument Pouch, 1 chamber	38 x 40 cm	33 109	110 pcs. / 440 pcs.
Instrument Pouch, 2 chambers	38 x 40 cm	33 110	110 pcs. / 440 pcs.
Perineal Drape	20 x 40 cm	33 100	140 pcs. / 560 pcs.
Rectal Shield	60 x 90 cm	33 111	70 pcs. / 280 pcs.

— adhesive side

Raucodrape® Instrument Covers, sterile

	Size	REF	Shipping Units
Chair Cover	∅ = 50 cm	33 124	30 pcs. / 120 pcs.
Control Panel Cover, round	80 x 90 cm	33 113	40 pcs. / 160 pcs.
Control Panel Cover, round	∅ = 90 cm	33 120	20 pcs. / 80 pcs.
Control Panel Cover, round	∅ = 120 cm	33 121	15 pcs. / 60 pcs.
Control Panel Cover, rectangular	80 x 145 cm	33 112	25 pcs. / 100 pcs.
Image Converter Set, 3-piece	80 x 45 cm	33 122	30 pcs. / 120 pcs.
Image Converter Set, one-piece	-	33 114	36 pcs. / 72 pcs.
Image Converter Set, one-piece	-	33 115	19 pcs. / 38 pcs.
Camera Drape, telescopic fold	13 x 250 cm	15 965	40 pcs. / 160 pcs.
Camera Drape, reverse-fold	13 x 250 cm	15 966	40 pcs. / 160 pcs.
Camera Drape with ring	15 x 250 cm	33 119	20 pcs. / 80 pcs.
Camera Drape	16 x 250 cm	33 118	25 pcs. / 50 pcs.
Armrest Cover (2 pcs.)	35 x 73 cm	33 123	46 pcs. / 184 pcs.

Raucodrape® Surgical Incision Drape, sterile

	Size	REF	Shipping Units
Individually sealed in a special pouch	12 x 17 cm	25 440	10 pcs. / 80 pcs.
	15 x 20 cm	25 441	10 pcs. / 60 pcs.
	30 x 20 cm	25 442	10 pcs. / 80 pcs.
	45 x 20 cm	25 443	10 pcs. / 80 pcs.
	40 x 35 cm	25 444	11 pcs. / 80 pcs.
	45 x 50 cm	25 445	12 pcs. / 80 pcs.
	45 x 65 cm	25 446	10 pcs. / 40 pcs.
	55 x 80 cm	25 447	10 pcs. / 40 pcs.
	72 x 88 cm	25 448	10 pcs. / 40 pcs.

— adhesive side



Sentinex[®] Surgical Face Masks

Sentinex[®] Surgical Caps



Sentinex lite



Sentinex CLASSIC



Sentinex SuperSoft



Sentinex Extra Touch



Sentinex Safety Shield with anti-fog visor



Sentinex Ultra



Sentinex AntiFog



Sentinex Soft



Sentinex Fancy



Sentinex Astro



Sentinex Astro Extra



Sentinex Color



Sentinelx Nova



Sentinelx Trend



Sentinelx Flex



Sentinelx Easy

Sentinelx Surgical Face Masks in different versions, made as if “to measure”, e.g. with optimum wearer comfort and proven BFE value* for efficient prophylaxis against infections for patients and staff, in the OR and doctor’s surgery, in intensive care units and in areas with high infection tendency. The extremely long nose bow is made of plastic covered steel

and is integrated into soft, skin-friendly seams and offers good fit around the nose. Filtration material, inner and outer layers and the ties are ultrasound sealed – a hygienically perfect and high quality finish. All Sentinelx Surgical Face Masks and Sentinelx Surgical Caps are free of latex and fiberglass.

Sentinelx® Surgical Face Masks

	REF	Shipping Units
Sentinelx Lite, blue	11 983	50 pcs. / 300 pcs.
Sentinelx Lite, green	11 991	50 pcs. / 300 pcs.
Sentinelx Lite, green, with ear loops, latex-free	11 992	50 pcs. / 500 pcs.
Sentinelx CLASSIC, green	11 996	50 pcs. / 300 pcs.
Sentinelx SuperSoft	11 997	50 pcs. / 300 pcs.
Sentinelx Extra Touch	11 998	50 pcs. / 300 pcs.
Sentinelx Safety Shield	11 989	25 pcs. / 100 pcs.
Sentinelx AntiFog	11 986	50 pcs. / 300 pcs.
Sentinelx Ultra	11 999	50 pcs. / 300 pcs.

Sentinelx® Surgical Caps

Sentinelx Soft	Beret cap, green, breathable	17 420	100 pcs. / 600 pcs.
Sentinelx Astro	Astro cap, green	17 421	100 pcs. / 600 pcs.
Sentinelx Astro Extra	Astro cap with sweatband, green	17 422	100 pcs. / 600 pcs.
Sentinelx Fancy	Nurse’s cap with center seam, green	17 423	100 pcs. / 600 pcs.
Sentinelx Color	Nurse’s cap with center seam, colored	17 430	100 pcs. / 600 pcs.
Sentinelx Nova	Nurse’s cap with strings, green	17 424	100 pcs. / 600 pcs.
Sentinelx Trend	Nurse’s cap, long, green	17 425	100 pcs. / 600 pcs.
Sentinelx Flex	Surgeon’ cap, green	17 426	100 pcs. / 600 pcs.
Sentinelx Easy, green	Beret cap, green	17 427	100 pcs. / 600 pcs.
Sentinelx Easy, blue	Beret cap, blue	17 428	100 pcs. / 600 pcs.
Sentinelx Easy, white	Beret cap, white	17 429	100 pcs. / 600 pcs.



Sentinex® Surgical Gowns • OR Stockings

Sentinex Surgical gowns by L&R represent a clearly structured concept:

- Spunlace and SMMS
- Basic, Standard, Special versions
- With and without sterile wrapping (“direct”)

Sentinex OR gowns offer high wearer comfort and protection for users and patients.




Two materials: Spunlace and SMMS

With the introduction of EN 13795 uniform requirements have been defined for OR gowns. The central requirement of this standard is the compliance with a barrier function against fluids and germs in order to ensure protection against infections, for the patient as well as for the user.

Both materials fulfill the requirements of the standard and are suitable for all types of surgical procedures – if reinforced through additional reinforcements in the critical areas.

Spunlace

- excellent wearing comfort
- skin-friendly and breathable
- ideally suited for physical demanding operations
- best possible temperature exchange




Surgical gown with sterile wrapping



- ideally suited for sterile preparation
- with sterile wrapping
- 2 towels

SMMS

- light
- protects reliably against bacteria and fluids
- guarantees the circulation of air, heat, and water vapour
- reduces thermo-physiological stress



Surgical gown “direct”



- for immediate sterile handover
- for fast removal
- without sterile wrapping and towels

Sentinex[®] Surgical Gown Basic

Gown with tie-closure for small operations and interventions, with 3 ties for a safe closure in the back.



Sentinex[®] Surgical Gown Standard

Wraparound surgical gown with a sterile back and a specially wide overlapping area. Ideally suited for operations with small amounts of fluids.



Sentinex[®] Surgical Gown Special

Wraparound surgical gown on the basis of the Standard version, with additional impermeable reinforcement in the front and sleeve area. Ideally suited for operations with a large amount of fluids.



Sentinex[®] Surgical Gown Spunlace, sterile

Contents	Size	REF	Shipping Units
Standard with sterile wrapping	120 cm	19 300	44 pcs.
	130 cm	19 301	40 pcs.
	150 cm	19 302	36 pcs.
	170 cm	19 303	32 pcs.
	150 cm large	19 304	32 pcs.
Special with sterile wrapping	120 cm	19 310	40 pcs.
	130 cm	19 311	36 pcs.
	150 cm	19 312	32 pcs.
	170 cm	19 313	28 pcs.
	150 cm large	19 314	28 pcs.
Standard "direct"	130 cm	19 321	40 pcs.
	150 cm	19 322	36 pcs.
Special "direct"	130 cm	19 331	36 pcs.
	150 cm	19 332	32 pcs.

Sentinex[®] Surgical Gowns SMMS, sterile

Standard with sterile wrapping	120 cm	19 340	52 pcs.
	130 cm	19 341	48 pcs.
	150 cm	19 342	44 pcs.
	170 cm	19 343	40 pcs.
	150 cm large	19 344	40 pcs.
Special with sterile wrapping	120 cm	19 350	48 pcs.
	130 cm	19 351	44 pcs.
	150 cm	19 352	40 pcs.
	170 cm	19 353	36 pcs.
	150 cm large	19 354	36 pcs.
Standard "direct"	120 cm	19 360	52 pcs.
	130 cm	19 361	48 pcs.
	150 cm	19 362	44 pcs.
Special "direct"	120 cm	19 380	48 pcs.
	130 cm	19 381	44 pcs.
	150 cm	19 382	40 pcs.

Sentinex[®] Surgical Stockings

Sentinex Surgical Stocking	20 225	200 pcs. / 200 pcs.
Sentinex Surgical Sock	20 224	200 pcs. / 400 pcs.

Sentina Examination Gloves

Sentina CP Examination Gloves

Made of tear-resistant, transparent sealed copolymer foil; hard-wearing.

For the hygienic and economical everyday requirements in hospitals and doctor's surgeries.



Sentina CP Examination Gloves			
	Size	REF	Shipping Units
sterile individually sealed	S	51 049	100 pcs. / 1000 pcs.
	M	51 046	100 pcs. / 1000 pcs.
	L	51 043	100 pcs. / 1000 pcs.
sterile sealed in pairs	S	51 048	50 pairs / 500 pairs
	M	51 045	50 pairs / 500 pairs
	L	51 042	50 pairs / 500 pairs
Examination Fingerstalls, made of Latex			
non-sterile, 100 pieces in a foil pouch	2	52 080	100 pcs. / 25000 pcs.
	3	52 081	100 pcs. / 25000 pcs.
	4	52 082	100 pcs. / 25000 pcs.
	5	52 083	100 pcs. / 25000 pcs.
	2-5	52 043	100 pcs. / 25000 pcs.
Cotton surgical gloves, non-sterile			
PE bags with 1 pair	7,5	17 574	20 pcs. / 1200 pcs.
	8	17 575	20 pcs. / 1200 pcs.
	8.5	17 576	20 pcs. / 1200 pcs.
	9	17 577	20 pcs. / 1000 pcs.
	9.5	17 578	20 pcs. / 1000 pcs.
	10	17 579	20 pcs. / 1000 pcs.
Sentina Ambidextrous Latex Examination Gloves			
non-sterile powdered with rolled rim, light	6-7 S	51 120	100 pcs. / 1000 pcs.
	7-8 M	51 121	100 pcs. / 1000 pcs.
	8-9 L	51 122	100 pcs. / 1000 pcs.
powder-free with rolled rim, light	6-7 S	51 125	100 pcs. / 1000 pcs.
	7-8 M	51 126	100 pcs. / 1000 pcs.
	8-9 L	51 127	100 pcs. / 1000 pcs.

Setpack® OR Dressing Materials

Flexible combination of standardised components

Thorough preparation takes its time – also in emergency situations. What is required is immediate accessibility and safety.

The Setpack Standard System offers ready-for-use sets which contain the **surgical dressing materials** most commonly used in operations – sterile and ready at any time.

Ready-for-use Standard Sets contain **wound dressings, dissecting swabs, gauze balls or laparotomy swabs** in quantities as they are needed. All Setpack products are provided with an X-ray detectable marker (min. 60% barium sulfate). Should a case of doubt arise, an X-ray control will give immediate



clarification. Depending on the requirements for your specific operations, you can make your individual choice. Determined quantities offer maximum economy.



Patented Double Counting Card System

- dual check system by having each counting card both inside the pack and in the peel bag
- risk-free product handover during the surgery
- clear documentation by the use of self-adhesive counting cards



Special covering of the inside pack

- a handy flap makes it easy to open the package
- fast product removal through special packaging design



Patented packaging box for gauze balls

- increased safety by optimized counting control before and after the surgery
- special developed packaging design for handsome removal of dissecting swabs



Setpack® OR Dressing Materials Standard Sets

Setpack® Gauze Dressings, X-ray detectable, sterile					
		Size	REF	Shipping Units	
Ready-for-use dressings made of absorbent gauze type 17 (17 threads/cm ²) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	Sets of 10 pcs.	12-ply	7.5 x 7.5 cm	15 114	28 / 84 sets
			10 x 10 cm	15 100	20 / 60 sets
			10 x 12.5 cm	15 116	20 / 60 sets
			10 x 20 cm	15 102	16 / 48 sets
	16-ply	5 x 5 cm	15 110	34 / 102 sets	
		7.5 x 7.5 cm	15 106	24 / 72 sets	
		10 x 10 cm	15 101	20 / 60 sets	
		10 x 12.5 cm	15 111	18 / 54 sets	
Absorbent gauze type 20 (20 threads/cm ²) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	Sets of 10 pcs.	32-ply	10 x 10 cm	15 104	12 / 36 sets
Type 17 (17 threads/cm ²) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	Sets of 20 pcs.	12-ply	7.5 x 7.5 cm	15 119	20 / 60 sets
			10 x 10 cm	15 103	18 / 54 sets
			10 x 12.5 cm	15 109	10 / 30 sets
			10 x 20 cm	15 108	8 / 24 sets
	16-ply	5 x 5 cm	15 117	26 / 78 sets	
		7.5 x 7.5 cm	15 118	16 / 48 sets	
		10 x 10 cm	15 105	14 / 42 sets	
		10 x 12.5 cm	15 112	8 / 24 sets	
Setpack® Gauze Balls, X-ray detectable, sterile					
Ready-for-use gauze balls made of absorbent gauze type 20 (20 threads/cm ²) with elastic, for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card	Sets of 5 pcs.	fist size	7	22 800	16 / 48 sets
	Sets of 10 pcs.	walnut size	2	22 801	72 / 144 sets
			3	22 802	38 / 114 sets
			4	22 803	26 / 78 sets
			6	22 804	18 / 54 sets
	Sets of 15 pcs.	plum size	3	22 805	30 / 90 sets
			4	22 806	20 / 60 sets
			6	22 807	12 / 36 sets
	Sets of 20 pcs.	plum size	3	22 808	26 / 78 sets
			4	22 809	16 / 48 sets
			6	22 810	6 / 18 sets
Setpack® Dissecting Swabs, sterile					
Ready-for-use dissecting swabs made of absorbent gauze type 24 (24 threads/cm ²) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	Sets of 10 pcs.	small	2	12 780	40 / 120 sets
		medium	3	12 781	40 / 120 sets
		large	4	12 782	40 / 120 sets
	Sets of 20 pcs.	small	2	12 783	40 / 120 sets
		medium	3	12 784	40 / 120 sets
		large	4	12 785	40 / 120 sets

Setpack® OR Dressing Materials

Setpack® Toptex® lite Laparotomy Swabs, X-ray detectable

Ready-for-use laparotomy swabs made of absorbent gauze type 20 (20 threads/cm²) EN 14079, prewashed, for use in the OR, with X-ray detectable chip and self-adhesive, double, separate counting card.



Toptex® lite Laparotomy Swabs, X-ray detectable, sterile

Ready-for-use laparotomy swabs made of absorbent gauze type 20 (20 threads/cm²), prewashed, for use in the OR, with X-ray detectable chip and double, separated, adhesive counter card.

			Size	REF	Shipping Units	
White	set of 2 pcs.	4 layers	40 x 40 cm	15 004	14 / 42 sets	
	set of 5 pcs.	4 layers		15 001	6 / 18 sets	
	set of 2 pcs.	6 layers		15 012	12 / 36 sets	
	set of 5 pcs.	6 layers	20 x 30 cm	15 003	6 / 18 sets	
	set of 2 pcs.	4 layers		15 011	28 / 84 sets	
	set of 5 pcs.	4 layers		15 007	14 / 42 sets	
	set of 2 pcs.	4 layers		50 x 60 cm	15 014	8 / 24 sets
	set of 2 pcs.	4 layers			15 010	26 / 78 sets
	set of 2 pcs.	4 layers			15 010	26 / 78 sets
	Green	set of 2 pcs.	4 layers	40 x 40 cm	15 006	14 / 42 sets
set of 5 pcs.		4 layers	15 002		6 / 18 sets	
set of 2 pcs.		6 layers	15 008		12 / 36 sets	
set of 5 pcs.		6 layers	20 x 30 cm	15 000	6 / 18 sets	
set of 2 pcs.		4 layers		15 017	28 / 84 sets	
set of 5 pcs.		4 layers		15 005	14 / 42 sets	
set of 2 pcs.		4 layers		50 x 60 cm	15 013	8 / 24 sets
set of 2 pcs.		6 layers			15 018	6 / 18 sets
set of 2 pcs.		4 layers			15 009	26 / 78 sets
Blue		set of 2 pcs.	4 layers	40 x 40 cm	15 016	14 / 42 sets
	set of 5 pcs.	4 layers	15 015		6 / 18 sets	

Toptex® lite Disposable Laparotomy Swabs, X-ray detectable, non-sterile

made of absorbent gauze Ph. Eur. type 20 (20 threads/cm²) with 30 cm loop, prewashed, with X-ray detectable chip.

4 layers	White	8 x 90 cm	10 921	100 / 600 pcs.	
		20 x 30 cm	10 922		
		40 x 40 cm	10 923		
		50 x 60 cm	10 925		
	Green	8 x 90 cm	10 916	100 / 600 pcs.	
		20 x 30 cm	10 917		
		40 x 40 cm	10 918		
		50 x 60 cm	10 920		
	Blue	40 x 40 cm	10 915	90 / 180 pcs.	
	6 layers	White	40 x 40 cm	10 924	75 / 150 pcs.
		Green	40 x 40 cm	10 919	75 / 150 pcs.



Setpack[®] OR Dressing Materials Individual Sets

The Setpack Individual Range offers the possibility to create OR Dressing Material Sets according to your requirements.

Together with our medical consultants you decide on quantities and sizes of wound and absorbent dressings, gauze balls, dissecting swabs, cotton tipped

applicators, ribbon gauze, laparotomy swabs and bandages. We take care of the professional production of your "made to measure" set.

The Individual Set should form the basis for the specific operation. Additional requirements will be supplied from the Standard Set.



Kitpack[®] OR Custom Procedure Tray

Kitpack is an efficient OR system made up from various ready-for-use OR components. Apart from dressing materials of various types and quantities, OR drape systems and clothing, it contains additional disposable medical products. For example drains, disposable syringes and cannulae, suture materials and other disposable components. In compact packaging, sterile and ready-for-use.

Sets are compiled standardized and yet individually: determined for the requirements in your hospital, selected for the surgical intervention for which it is intended.

By using prefabricated Kitpacks, you will reduce preparation time and procedure costs – with full safety and the necessary flexibility.



Summary of components:

Lohmann & Rauscher products:

- Gauze dressings, X-ray detectable
- Gauze balls, X-ray detectable
- Dissecting swabs, X-ray detectable
- Toptex lite disposable laparotomy swabs, X-ray detectable
- Absorbent dressings and wound Dressings
- Raucodrape OR drape system
- Sentinex OR gowns
- Gloves
- Bandages
- Hygiene products
- Adhesive dressings
- Raucodrape surgical incision drapes

Disposable OR components:

- Drainage systems
- Disposable syringes and cannulae
- Disposable scalpels and blades
- Infusion systems
- Catheters
- Suture materials
- Other disposable components that can be integrated into the sterilization cycle



Kitpack® – designed for specific indications, following your specifications:

- Individually compiled OR set
- Contains all required dressing materials and medical products
- Sterile and ready-for-use
- Compact packaging to suit OR requirements