



Products for Dermatology



Curatest® Curatest® F

Curatest® Nonwoven Patch Test Strips

The patch test strip with 10 test discs is as simple to apply as a first-aid dressing. The detachable side edges can be written on to assist diagnosis.

Nonwoven carrier material: 100% polypropylene
Frame of test discs: polyester film
Absorbent test discs: 60% cotton, 40% viscose
Polyacrylate adhesive: free of colophonium and colophonium derivatives, hypoallergenic.

Curatest® F Film Patch Test Strips

This patch test strip with 5 or 10 test discs has a semi-permeable, transparent carrier film that is permeable to water vapour but waterproof. It adheres reliably, like a second skin, even when taking a shower. Easy application of the test substances to the discs and of the test strip to the skin.

Carrier film: 100% polyurethane.
Frame of test discs: polyester film.
Absorbent test discs: 60% cotton, 40% viscose.
Polyacrylate adhesive: free of colophonium and colophonium derivatives, hypoallergenic.



Curatest® Nonwoven Patch Test Strips

Size	REF	Shipping Units
each strip has 10 test discs 7.5 x 12.5 cm	30 054	50 pcs. / 1200 pcs.

Curatest® F Film Patch Test Strips

each strip has 10 test discs 7.5 x 12.5 cm	30 062	50 pcs. / 1200 pcs.
each strip has 5 test discs 3.2 x 12.5 cm	30 061	100 pcs. / 2000 pcs.

tg® Ready-to-wear Dressings

These ready-to-wear dressings are ideal for covering body parts affected by dermatological disorders, e. g. contact dermatitis, photodermatitis and medically treated skin.

For such indications, ready-to-wear dressings are quick and easy to apply. The material used was selected for its excellent hypoallergenic characteristics.



tg® Glove

Soft, supple hand/finger dressing shaped for a comfortable fit. Also for preventing superinfections.

100% cotton, washable, sterilisable (A 121°C).

tg® Mitten

Thumbless hand dressing with integrated tie strings for securing at the wrist. As protective covering for dermatological disorders, and to prevent children with itchy skin conditions from scratching themselves.

100% cotton, washable, sterilisable (A 121°C).

tg® Glove		REF	Shipping Units
each pair individually in a plastic bag			
for children		24 749	1 pair / 100 pairs
for adults	Gr. 6-7 (small)	24 750	1 pair / 50 pairs
	Gr. 7½-8½ (medium)	24 751	1 pair / 50 pairs
	Gr. 9-10 (large)	24 752	1 pair / 50 pairs
tg® Mitten			
each pair individually small (infants)		24 746	1 pair / 100 pairs
in a plastic bag	large (small children)	24 747	1 pair / 100 pairs

tg® Vest

Short-sleeved vest, soft, supple and hypoallergenic, as a practical ready-to-wear dressing for covering extensive skin disorders; also as a time-saving under-cast protection for P.o.P. and synthetic casts.

The elastic knitted material has widthwise extensibility allowing for intermediate sizes.

Ribbed material with widthwise extensibility, raw white, 100% cotton, not degreased or optically brightened, particularly hypoallergenic, washable, sterilisable (A 121°C).



tg® Vest

	REF	Shipping Units
individually in a cardboard box		
small (ladies' size 36 - 44, men 42 - 50)	24 044	1 pc. / 14 pcs.
large (ladies' up to size 46, men up to size 56)	24 045	1 pc. / 9 pcs.

Curaderm® All-in-one Neurodermatitis Suit

Single-piece suit made from unbleached and undyed cotton for protecting infants and children with neurodermatitis, psoriasis and other itching diseases. Toes and fingernails are securely enclosed in the integrated socks and mittens whilst the elasticated band adjustable

by parents and carers prevents hands and feet from slipping out.

New skin defects are thus avoided, while itching is steadily reduced. Children are able to move freely and feel comfortable in this breathable all-in-one suit.



Curaderm® All-in-one Neurodermatitis Suit

Made from 100% cotton, unbleached and undyed, zip fasteners and externally closable elastic draw strings have no skin contact; seams are on the outside; washable up to 95°C.

Curaderm® All-in-one Neurodermatitis Suit

	Size	REF	Shipping Units
individually in a plastic bag	74 – 80 cm	24 300	1 pc. / 28 pcs.
	86 – 92 cm	24 301	1 pc. / 24 pcs.
	98 – 104 cm	24 302	1 pc. / 20 pcs.
	110 – 116 cm	24 303	1 pc. / 20 pcs.
	122 – 128 cm	24 304	1 pc. / 18 pcs.
	134 – 140 cm	24 305	1 pc. / 15 pcs.