



People.Health.Care.

VENOSAN®

Medical support and compression stockings



Our prevention and treatment solutions for lower limbs.



Product overview

Medical support stockings.

Travel and leisure socks

VENOSAN® COTTON SUPPORT TRAVEL & LEISURE

- For non heavy legs in daily life and while travelling
- Promotes blood circulation
- Efficient thrombosis prophylaxis
- 100% cotton on the skin
- With soft band that does not constrict
- Can be worn with jeans and suits

Use

- For people who sit or stand a lot at work
- Long journeys by plane, train, and car



15–20
mmHg

Composition
65 % polyamide
20 % cotton
15 % elastane

Colours*

White	Olive	Silver
Beige	Jeans	Anthracite
Wood	Navy	Black

Sizes

	S	M	L	XL
EUR	36–39	39–42	42–45	45–48
GB	3–6	6–8	8–10.5	10.5–13
USA / CAN	4–6.5	6.5–8.5	8.5–11	11–13

Style



VENOSAN® SILK SUPPORT TRAVEL & LEISURE

- For non heavy legs in daily life and while travelling
- Promotes blood circulation
- Efficient thrombosis prophylaxis
- With soft band that does not constrict
- Soft, smooth silk fibres on the skin
- Elegant appearance

Use

- For people who sit or stand a lot at work
- Long journeys by plane, train, and car



15–20
mmHg

Composition
65 % polyamide
20 % elastane
15 % silk

Colours*

White	Sisko (Violet)	Marine
Beige	Silver	Black
Apricot	Jeans	

Sizes

	S	M	L
EUR	36–39	39–42	42–45
GB	3–6	6–8	8–10.5
USA / CAN**	5.5–8	8–10	10–12.5

Style



* Actual product colours may vary slightly from the colours in this catalogue due to the printing process.

** Women's sizes

VENOSAN® JET LEGS®

TRAVEL WITH COMFORT SOCKS

- Promotes blood circulation and prevents blood congestion
- Efficient thrombosis prophylaxis

Use

- Wear for plane, train, and car travel
- Long journeys



15–20 mmHg

Composition
86 % nylon polyamide
14 % LYCRA® elastane

Colours*



Black

Sizes

S-M
L-XL

Style



Anti-embolism stockings

VENOSAN® AES

- Perfect fit thanks to the anatomical shape
- Colour coding for simple and reliable identification of the size
- Top inspection window for safety

Indications

- To prevent thrombosis and embolism before, during and after operations
- General immobility, confinement to bed
- Before, during and after childbirth



16–18 mmHg

Composition
AD + AGH:
91 % nylon polyamide
9 % LYCRA® elastane

Colours*



White

Sizes

		S	M	L	XL	XXL
AD	cB (cm)	18–20	20–22	23–25	25–27	28–30
	cC (cm)	to 32	32–37	37–42	42–47	47–52
	ID (cm)	S regular to 41	M regular to 43		Regular to 45	
		S long over 41	M long over 43		Long over 45	
AGH	cB (cm)	18–20	20–22	23–25	25–27	–
	cC (cm)	to 32	32–37	37–42	42–47	–
	cG (cm)	45–64	53–71	61–80	70–84	–
	IG (cm)	Short to 74 / Regular 74–79 / Long over 79				

Style



Product overview

Medical compression stockings and pantyhose.

VENOSAN® 3000 COTTON

NATURAL LOOK

- 100% cotton on the skin
- Absorbs moisture
- Durable
- Easy to don

Indications

- Chronic venous insufficiency CVI 2–3
- Mild lymphoedema
- Lipoedema stage 2–3
- Following vein surgery
- Post-thrombotic syndrome PTS



Composition

- CCL 2 49% cotton
27% nylon polyamide
24% elastane
- CCL 3 40% cotton
27% nylon polyamide
33% elastane

VENOSAN® 4000

CLIMATE EFFECT

- Highly elastic and comfortable thanks to fine microfibre
- Climate effect
- Semi-transparent

Indications

- Chronic venous insufficiency CVI 1
- Pregnancy-induced varicose veins
- Following vein surgery
- Swollen, heavy legs
- Spider veins



Composition

- 71% TACTEL® polyamide
29% elastane

VENOSAN® 5000

THE MAGIC OF SILVER

- Odour-inhibiting
- Temperature-equalising
- Anti-static

Indications

- Chronic venous insufficiency CVI 1–2
- Lipoedema stage 1
- Following vein surgery
- Varicose related mild skin problems



Composition

- 40% TACTEL® polyamide
30% elastane
20% cotton
10% polyamide/silver (X-Static®)

VENOSAN® 6000 CLASSIC

- Slight transparency
- Made of high elastic yarns for a perfect fit
- Cost-efficient stocking that satisfies all requirements
- Breathable
- Ease of care

Indications

- Chronic venous insufficiency CVI 1-2
- Varicose veins
- Pregnancy-induced varicose veins
- Following vein surgery



Composition
70% nylon polyamide
30% elastane

VENOSAN® 7000 NATURAL TOUCH

- High stiffness (resistance)
- Robust
- Temperature-equalising
- Skin-friendly

Indications


- Chronic venous insufficiency CVI 2-3
- Mild lymphoedema
- Lipoedema stage 2-3
- Following vein surgery
- Following acute thrombosis (after bandaging)
- Post-thrombotic syndrome PTS
- Marked orthostatic stress



Composition
62 % TACTEL® polyamide
26 % elastane
8 % cotton
4 % Lyocell/SeaCell®

Colours and styles

Medical compression stockings and pantyhose.




VENOSAN® 3000

	AD	ADH	AG	AGH	AGH	AGG	AT	ATU	ATH
CCL 2, 23–32 mmHg									
Standard size	1	1	–	1	–	–	–	–	–
Custom size	1	1	1	1	1	1	1	–	1
CCL 3, 34–46 mmHg									
Standard size	1	–	–	–	–	–	–	–	–
Custom size	1	1	1	1	1	1	1	–	1


VENOSAN® 4000

	AD	ADH	AG	AGH	AGH	AGG	AT	ATU	ATH
CCL 1, 18–21 mmHg									
Standard size	2 3 4 7	2 3 7	–	2 3 4 7	2 3 4 7	–	2 3 4 7	2 7	–
Custom size	2 7	2 7	2 7	2 7	2 7	2 7	2 7	2 7	–
CCL 2, 23–32 mmHg									
Standard size	2 3 4 5 6 7	2 3 7	–	2 3 4 5 6 7	2 3 4 5 6 7	2	2 3 4 7	2 7	–
Custom size	2 3 4 7	2 3 4 7	2 3 4 7	2 3 4 7	2 3 4 7	2 7	2 3 4 7	2 3 4 7	–

Colours*



* Actual product colours may vary slightly from the colours in this catalogue due to the printing process.



VENOSAN® 5000

	AD	ADH	AG	AGH	AGH ☺	AGG	AT	ATU	ATH
CCL 1, 18–21 mmHg									
Standard size	2 3 4 7	2 3 7	–	2 3 4 7	2 3 4 7	2	2 3 4 7	–	–
Custom size	2 7	2 7	2 7	2 7	2 7	2 7	2 7	2 7	–
CCL 2, 23–32 mmHg									
Standard size	2 3 4 7	2 3 7	–	2 3 4 7	2 3 4 7	2	2 3 4 7	–	–
Custom size	2 3 4 7	2 3 4 7	2 3 4 7	2 3 4 7	2 3 4 7	2 7	2 3 4 7	2 3 4 7	–
CCL 3, 34–46 mmHg									
Standard size	2 3 7	–	–	–	–	–	–	–	–
Custom size	2	2	2	2	2	2	2	2	–

VENOSAN® 6000

	AD	ADH	AG	AGH	AGH ☺	AGG	AT	ATU	ATH
CCL 1, 18–21 mmHg									
Standard size	1 7	–	–	1 7	–	–	1 7	–	–
CCL 2, 23–32 mmHg									
Standard size	1 7	1 7	–	1 7	–	–	1 7	–	–
CCL 3, 34–46 mmHg									
Standard size	1	–	–	1	–	–	1	–	–

VENOSAN® 7000

	AD	ADH	AG	AGH	AGH ☺	AGG	AT	ATU	ATH
CCL 1, 18–21 mmHg									
Standard size	1 7	1 7	–	1 7	–	1	–	–	–
Custom size	1 7	1 7	1 7	1 7	1 7	1 7	1 7	–	1 7
CCL 2, 23–32 mmHg									
Standard size	1 7	1 7	–	1 7	–	1	–	–	–
Custom size	1 7	1 7	1 7	1 7	1 7	1 7	1 7	–	1 7
CCL 3, 34–46 mmHg									
Standard size	1 7	1 7	–	1 7	–	–	–	–	–
Custom size	1 7	1 7	1 7	1 7	1 7	1 7	1 7	–	1 7



People.Health.Care.



Publisher

Lohmann & Rauscher GmbH & Co. KG · Irlicher Straße 55 · D-56567 Neuwied
Lohmann & Rauscher GmbH · Kirchengasse 17 · A-2525 Schönan an der Triesting

Contact addresses

Lohmann & Rauscher GmbH & Co. KG · Westerwaldstraße 4 · D-56579 Rengsdorf
Lohmann & Rauscher GmbH · Johann-Schorsch-Gasse 4 · A-1141 Vienna

VENOSAN® is a registered trademark of
SWISSLASTIC AG St. Gallen, Switzerland – an  Company