



People.Health.Care.

Rosidal® 1C **NEW**

The new, faster
compression solution.

win

effective

time-saving

easy

Rosidal® 1C NEW

The 3-in-1 short stretch bandage with skin layer, padding layer and compression layer.



Proven effectiveness:

- reliable, strong compression: high working and low resting pressure (up to 40–60 mmHg)
- the therapeutically relevant pressure is maintained for several days*
- can be worn day and night for up to seven days*
- good conformability and easy to apply due to cross elasticity

Outstanding efficiency:

- **easy** and safe application with the innovative Safe-Loc system
- 3-in-1 system enables quick application and **saves time**

Feel the comfort:

- textile terry structure
- slim design allows it to be worn in most types of shoes in comparison to multi component systems

Sources:

* Based on preliminary results of observational study NCT06728137 (www.clinicaltrials.gov)

Rosidal® 1C NEW

Maximum benefit with just one short stretch bandage.

Intended use:

For compression, retention of dressings, support and relief and for immobilisation of limbs.

- Indications:**
- varicosis
 - venous insufficiency
 - thrombophlebitis
 - phlebothrombosis
 - post-thrombotic syndrome
 - stage I to III b chronic venous insufficiency according to Widmer
 - venous and mixed ulcers
 - oedema (phleboedema, post-traumatic, post-operative, during pregnancy)
 - after sclerotherapy and venous surgery
 - contusions, sprains, tendon injuries and joint dislocations



Available in two widths:
8 cm & 10 cm

Trust your strength with the Safe-Loc system:
pull to full stretch, apply and you're done.

The Safe-Loc system makes Rosidal 1C safe, quick and easy to use.



- Rosidal 1C is applied at full stretch – with **full tension**
- 
- it reduces overstretching of the bandage and excessive compression pressure, so is **very user-friendly and offers treatment safety**
 - the bandage can improve blood flow in the venous blood vessels, promoting **successful treatment outcomes**



Rosidal® 1C **NEW**

The 3-in-1 short stretch bandage.

Skin layer

Padding layer

Compression layer



Rosidal® 1C

Dimensions (cm x m stretched length)*	REF	Quantity	Shipping Units (pcs.)
8 x 6	162 154	1 roll	1/32
8 x 6	162 155	6 rolls	6/24
10 x 6	162 156	1 roll	1/60
10 x 6	162 157	6 rolls	6/24

* The width of the bandage should match the diameter of the extremity being bandaged.

Material

Polyamide 48%, cotton 42%, polyester 6%, elastane 4%



Does not contain latex

Natural rubber latex is not an intended ingredient in the production or packaging of the above products. We cannot completely exclude any potential contamination with traces of latex from production or the environment.

