



People.Health.Care.

Suprasorb® P non-adhesive **NEW**

Soft, conformable and highly absorbent.

With improved
breathability,
allowing excess
moisture to evaporate
from the dressing



The perforated wound
contact layer prevents
adhesion to the wound



Maintains a moist
wound environment*



Conforms
to the wound bed

Suprasorb® P non-adhesive NEW

Soft, conformable and highly absorbent.



Three-layer material:

Your benefit:

Waterproof film layer

breathable and permeable to vapour

Absorbent PU foam


efficient exudate management

Perforated non-adherent wound contact layer

prevents adhesion to the wound

Indications:

- acute and chronic wound such as: venous leg ulcers, pressure ulcers, diabetic foot ulcers, donor sites, abrasions and surgical wounds
- for moderate to highly exuding wounds; superficial wounds
- can be used as a secondary dressing for infected wounds and wounds at risk of infection as recommended by a physician



Suprasorb® P non-adhesive

Size (cm)	REF	Shipping unit (pcs.)
5x5	401 180	10/400
7.5x7.5	401 181	10/300
10x10	401 182	10/300
10x12	401 184	10/300
12.5x12.5	401 185	10/300
15x15	401 183	5/200
15x20	401 179	5/200
20x20	401 186	5/100

