



People.Health.Care.

Suprasorb® P Heel

The anatomical shape makes it easy to apply.
Soft, flexible and conformable.



Enables atraumatic
dressing changes



Minimizes the risk
of maceration



Great conformability
in awkward areas

Suprasorb® P Heel

Offers comfort and a cushioning protection, allows patient comfort.

Three-layer material:

Your benefit:

Waterproof film layer

breathable and permeable to vapour

Absorbent PU foam

efficient exudate management

Perforated non-adherent wound contact layer

enables atraumatic dressing changes

Indications:

- acute and chronic wound such as: venous leg ulcers, pressure ulcers, diabetic foot ulcers, donor sites, abrasions and surgical wounds
- for superficial wounds and in combination with a suitable wound filler for deep wounds
- for moderate to highly exuding wounds



Suprasorb® P Heel non-adhesive

Size (cm)	REF	Shipping unit (pcs.)
10 x 13	180040	5/120

