




People.Health.Care.

Advieskaart

Confectietabel VENOSAN®

Prestige 2 en 3 en Prestige 2 Katoen rondbrei beenkousen



Type	Confectie									Maatwerk
	Lengte	Kort (Alleen AD)	Normaal	Lang	Omvang	Small	Medium	Large	Extra large	Omvang
AT	/GT		29-35	29-35	cT cH	60-100 80-120	60-100 80-120	60-100 80-120	60-100 80-120	≤160 ≤160
AGT	/GT		30 ⁵ -31 ⁵	31 ⁵ -32 ⁵	cT	60-120	60-120	60-120	60-1200	≤160
AG met sil. Band			65-76	76-88	cG breed cG norm.	57-66 48-56	61-72 50-60	66-80 53-66	74-90 57-73	41-90
AG	/G		65-76	76-88	cG	48-66	50-72	53-80	57-90	41-90
	/F		55-64	64-74	cF	42-58	44-62	47-68	51-77	37-85
	/E		44-49	49-55	cE	32-41	34-44	37-48	41-54	32-75
AD	/D	35-38	38-42	42-47	cD	30-37	32-39	34-42	37-47	28-73
	/C	27-30	30-33	33-37	cC	31-39	33-41	35-44	39-48	28-65
	/B1	19-21	21-23	23-25	cB1	24-30	26-32	29-35	33-39	21-50
	/B	10-11	11-12	12-14	cB	19-21	21-23 ⁵	23-27	27-32	18-43
					cY	28-34	29-35	31-37	34-40	27-50
	/A	16-17	17-19	19-21	cA	20-24	21-25	22-27	24-30	18-38
	/Z	22-24	23-26	25-28						



VENOSAN®
Compression Stockings



People.Health.Care.



Confectie assortiment VENOSAN® Prestige 2 en 3 en Prestige 2 Katoen rondbrei beenkousen

	kleur	teenstuk	AD	AD silb. 3 cm / 5 cm	AG	AG silb. 5 cm	AT
VENOSAN® Prestige 2	● beige / ● zwart	open	● ●	● ●	●	●	-
		gesloten	● ●	● ●	-	●	●
VENOSAN® Prestige 2 Katoen	● beige	open	●	●	-	●	-
		gesloten	●	●	-	-	●
VENOSAN® Prestige 3	● beige	open	●	●	-	●	-
		gesloten	-	-	-	-	-

AD: kniekous
AG: dijbeen kous
AT: panty
silb.: siliconenband

	lengte AD (cm)	lengte AG (cm)	lengte GT voor AT (cm)	voetlengte IA / IZ (cm)
kort	35 - 38			16 - 17 / 22 - 24
normaal	38 - 42	65 - 76	29 - 35	17 - 19 / 23 - 26
lang	42 - 47	76 - 88		19 - 21 / 25 - 28



VENOSAN®
Compression Stockings