



People.Health.Care.

Compression Bandaging Product Catalog



Contents

Multi-layer Bandaging
Products for Lymphedema 4

Lymphedema Bandaging

Skin Protection

tg® tubular bandage	9
tg® cotton stockinette	10
tg® soft	11
tg® grip	12
tg® shape	14

Fingers/Toes

Transelast® Classic	15
Mollelast®	16
Raicolast®	17

Padding

Cellona®	18
Rosidal® soft	19

Focal Compression

Komprex®	20
Komprex® II	21

Bandaging

Rosidal® K	22
Idealbinde™	24
Lenkelast®	25
Durelast®	26
Dauerbinde® K	27

Tape

Silkafix®	28
-----------	----

Bandaging Systems

Rosidal® Lymphset	29
Rosidal® TCS	30
Rosidal® Sys	32

Specialty Bandages

Porelast®	33
Panelast®	34
tg® fix	35

Accessories

Bandage Roller	36
----------------	----

Multi-layer Bandaging Products for Lymphedema

Compression bandaging is the gold standard of lymphedema treatment, and L&R offers a comprehensive range of multi-layer bandaging products.

STEP 1 Skin Protection  tg[®] tubular stockinette  tg[®] cotton stockinette  tg[®] soft stockinette	STEP 2 Fingers/Toes  Transelast[®] Classic conforming bandage  Mollelast[®] conforming bandage  Raucolast[®] conforming bandage	STEP 3 Padding  Cellona[®] synthetic padding  Rosidal[®] soft foam padding
 Caresia[™] foam bandage liners		

STEP 4

Focal Compression



Komprex®
medical-grade foam
pads, rolls, and sheet



Komprex® II
channel foam sheet



SwellSpots™
edema and fibrosis pads



STEP 5

Bandaging



Rosidal® K
short-stretch bandage



Idealbinde®
wide short-stretch bandage



Lenkelast®
medium-stretch bandage



Durelast®
very-short-stretch bandage

STEP 6

Tape



Silkafix®
adhesive tape

Recommended Bandaging Products for Lymphedema

	STEP 1 Skin Protection	STEP 2 Fingers/Toes	STEP 3 Padding
	tg [®] tubular, tg [®] cotton, or tg [®] soft	Transelast [®] Classic, Mollelast [®] or Raucolast [®]	Cellona [®] or Rosidal [®] soft
Arm Fingertips/ Wrist to Axilla	1 piece cut slightly longer than the length of the arm. Select size based on patient need. OR  Caresia™ 1 garment. Select arm and/or glove/gauntlet style and size from sizing chart	6 cm: 2 rolls ■ Typically folded in half lengthwise before use	10 cm: 2–3 rolls
Below Knee Toes/Foot to just below Knee	1 piece cut slightly longer than the length of the lower leg. Select size based on patient need. OR  Caresia™ 1 garment. Select below knee style and size from sizing chart based on patient need	4 cm: 1 roll ■ Typically folded in half lengthwise before use	10 cm: 2–3 rolls
Full Leg Toes/Foot to Groin	2 pieces: one cut slightly longer than the length of the lower leg, one cut slightly longer than the length of the knee to the thigh. ■ These two pieces will need to overlap for full coverage. Select size based on patient need. OR  Caresia™ 2 garments. Select below knee and thigh style and size from sizing chart	4 cm: 1 roll ■ Typically folded in half lengthwise before use	10 cm: 2–3 rolls 15 cm: 2 rolls

This list is provided for your guidance, but please note the sizes and quantities may vary based on the size and shape of the limb. Patients are typically recommended to have two sets of multi-layer bandaging products — one to wash and one to wear.

STEP 4 Focal Compression	STEP 5 Bandaging	STEP 6 Tape
Komprex® or SwellSpots™	Rosidal® K	Silkafix®
Optional. Select type, size, and amount based on patient need.	6 cm: 1 roll 8 cm: 1 roll 10 cm: 2–3 rolls	1.25 cm Use enough to secure bandaging.
g chart based on patient need. ■ Replaces Steps 1–4.		
2 pads: Komprex 0/1 kidney-shaped pad, or SwellSpots Curved pad. ■ Used for the medial and lateral Malleoli ■ Additional focal compression can be used based on patient need	6 cm: 1 roll 8 cm: 1 roll 10 cm: 2–3 rolls	1.25 cm Use enough to secure bandaging.
on patient need. ■ Replaces Steps 1–4.		
2 pads: Komprex 0/1 kidney-shaped pad, or SwellSpots Curved pad. ■ Used for the medial and lateral Malleoli ■ Additional focal compression can be used based on patient need	6 cm: 1 roll 8 cm: 1 roll 10 cm: 2–3 rolls 12 cm: 3–4 rolls Idealbinde® 15 cm or 20 cm can also be used for larger thighs.	1.25 cm Use enough to secure bandaging.
hart based on patient need. ■ Replaces Steps 1–4.		

This page intentionally left blank.

tg[®] tubular

Tubular bandage



Features/Benefits

- Seamlessly knit tubular protection bandage
- Absorbs excess perspiration and oils
- Lightweight and breathable with close-mesh construction
- Conforms comfortably to body contours with crease-free application
- Stretches approximately four times its original width
- Retains shape; can be narrowed down by pulling lengthwise

- Comes in easy to use dispenser box
- Can be used to enclose chip foam for fibrosis pads

Indications

- Lymphedema
- Edema
- Venous leg ulcer (VLU)

Product composition

Cotton, viscose

Not made with natural rubber latex.



tg[®] Tubular bandage

individually-packed in a cardboard dispenser

Size	Width	REF	HCPCS Code	Shipping Units (per box/case)
20 m (22 yd) roll				
1 for finger and toe	1.4 cm (0.55 in.)	24000	A6607	1/20
2 for fingers, toes, child's hand	2.3 cm (0.91 in.)	24001	A6607	1/20
3 for several fingers, child's arm	3 cm (1.2 in.)	24002	A6607	1/10
5 for arm, lower leg, child's leg	5.5 cm (2.2 in.)	24003	A6607	1/10
6 for arm, leg (medium)	6.5 cm (2.6 in.)	24004	A6607	1/10
7 for arm (large), leg, child's head	7 cm (2.8 in.)	24005	A6607	1/18
9 for thigh (large), head, armpit	8.5 cm (3.3 in.)	24006	A6607	1/12
12 for thigh (extra large), head (large), child's trunk	12 cm (4.7 in.)	24007	A6607	1/12
10 m (11 yd) roll				
K1 for trunk (small)	16 cm (6.3 in.)	24030	A6607	1/10
K2 for trunk (large)	21 cm (8.27 in.)	24031	A6607	1/10

tg[®] cotton

Plain stockinette / Tubular bandage

Features/Benefits

- Seamlessly knit tubular protection bandage
- Absorbs excess perspiration and oils
- Lightweight and breathable with close-mesh construction
- Excellent for patients with sensitive skin
 - ◆ Skin-friendly due to its 100% cotton material composition
 - ◆ Soft and supple touch
- Conforms comfortably to body contours with crease-free application

- Stretches approximately two times its original width
- Retains shape; can be narrowed down by pulling lengthwise

Indications

- Lymphedema
- Edema
- Venous leg ulcer (VLU)

Product composition

100% raw cotton (unbleached)

Not made with natural rubber latex.



Plain stockinette/Tubular bandage Cotton stockinette

individually-wrapped in a plastic bag

Size	REF	HCPCS Code	Shipping Units (per box/case)
25 m (27.5 yd) roll			
6 cm (2.4 in.)	14820	A6607	1/10
8 cm (3.2 in.)	14821	A6607	1/10
10 cm (4 in.)	14822	A6607	1/5
12 cm (4.7 in.)	14823	A6607	1/5
15 cm (6 in.)	14824	A6607	1/6
20 cm (7.9 in.)	14825	A6607	1/3
25 cm (9.8 in.)	14826	A6607	1/3
30 cm (11.8 in.)	14827	A6607	1/3

tg® soft

Tubular padding bandage



Features/Benefits

- Seamlessly knit tubular protection padding bandage
- Absorbs high amounts of excess perspiration and oils
- Soft and comfortable terry cloth padding
- Skin-friendly with high unbleached cotton content
- Excellent for patients with small limbs, fragile skin, over boney prominences
- Conforms comfortably to body contours with crease-free application
- Stretches approximately five times its original width

- Retains shape; can be narrowed down by pulling lengthwise
- Comes in easy-to-use dispenser box

Indications

- Lymphedema
- Edema
- Venous leg ulcer (VLU)

Product composition

85% raw cotton (unbleached), 15% spandex

Not made with natural rubber latex.



tg® soft Tubular padding bandage

individually-packed in a cardboard dispenser

Size	Width		REF	HCPCS Code	Shipping Units (per box/case)
	stretched	unstretched			
10 m (11 yd) unstretched length					
S for forearm, child's arm, child's leg	30–32 cm (11.8–12.6 in.)	5.5–6.5 cm (2.2–2.6 in.)	13200	A6608	1/8
M for lower leg, arm	41–48 cm (16.1–19 in.)	8–10 cm (3.2–4 in.)	13201	A6608	1/6
L for leg	56–66 cm (22.1–26 in.)	13–15 cm (5.1–5.9 in.)	13202	A6608	1/4

tg® grip

Elasticated tubular support bandage

**Field of Application**

- Light compression for effective management of strains, sprains, edema, and soft tissue injuries
- Support and relief of ligaments, tendons and joints
- Dressing retention

Features/Benefits

- Delivers uniform light compression and support (<15mmHg)
- Seamless tubular design
- Lightly cushioned cotton construction for comfort
- Skin-friendly high cotton content
- Conforms well to the body with virtually no creases
- Flexible and extensible for freedom of movement
- Air-permeable to provide cool comfort and moisture management

- Widthwise elasticity to stay in place
- Practical quick and easy application and removal (without need of clips or tape)
- Available in the widest range of sizes available in the market

Tip

- Compression can be enhanced with a second layer of a smaller size

Care instructions

- Washable (60°C/140°F) and Reusable
- Steam sterilization (134°C/273°F)

Product composition

Cotton, elastodiene (natural rubber latex), polyamide yarns



tg® grip Elasticated tubular support bandage

individually-packed in a cardboard dispenser

Size	Width	REF	HCPSC Code	Shipping Units (per box/case)
10 m (11 yd) roll				
A for children's hand/forearm	4.5 cm (1.75 in.)	24120	A6457	1/24
B for children's arms/lower leg, small hand/wrist	6 cm (2.33 in.)	24121	A6457	1/16
C for slim arms, children's leg, medium hand/wrist, small elbow/ankle	6.5 cm (2.75 in.)	24122	A6457	1/16
D for medium arm, lower leg, large hand/wrist, medium elbow/knee	7.5 cm (3 in.)	24123	A6457	1/16
E for large arm/elbow/ankle, medium leg/knee, small thigh, XL hand/wrist	8.75 cm (3.5 in.)	24124	A6457	1/12
F for large leg/knee, medium thigh, XL elbow/ankle	10 cm (4 in.)	24125	A6457	1/12
G for XL leg/knee, large thigh	12 cm (4.75 in.)	24126	A6457	1/8
J for XL thigh, small torso/trunk	17.5 cm (7 in.)	24127	A6457	1/6
K for medium torso/trunk	21.5 cm (8.5 in.)	24128	A6457	1/4
L for large torso/trunk	32.5 cm (13 in.)	24129	A6457	1/4
M for XL torso/trunk	37.5 cm (15 in.)	24130	A6457	1/2

tg® shape

Elasticated Tubular Support Bandage



Field of Application

- Light compression for effective management of strains, sprains, edema, and soft tissue injuries
- Support and relief of ligaments, tendons and joints
- Dressing retention

Features/Benefits

- Delivers uniform light compression and support (<15 mmHg)
- Seamless tubular design
- Lightly cushioned cotton construction for comfort
- Skin-friendly high cotton content
- Conforms well to the body with virtually no creases
- Provides graduated compression for lower limb
- Flexible and extensible for freedom of movement
- Air-permeable to provide cool comfort and moisture management

- Widthwise elasticity to stay in place
- Practical quick and easy application and removal (without need of clips or tape)

Tip

- Compression can be enhanced with a second layer or smaller size of tg shape

Care instructions

- Washable (60°C/140°F) and Reusable
- Steam sterilization (134°C/273°F)

Product composition

Cotton, elastodiene (natural rubber latex), polyamide yarns



tg® shape Shaped elasticated tubular support bandage

individually-packed plastic bag

Size	Calf Circumference	REF	HPCPS Code	Shipping Units (per box)
Full Leg - S	32–36 cm (12.5–14.25 in.)	88900	A6457	10
Full Leg - M	35–39 cm (13.75–15.25 in.)	88901	A6457	10
Full Leg - L	38–42 cm (15–16.5 in.)	88902	A6457	10
Full Leg - XL	41–45 cm (16.25–17.75 in.)	88903	A6457	10
Below Knee - S	32–36 cm (12.5–14.25 in.)	88904	A6457	10
Below Knee - M	35–39 cm (13.75–15.25 in.)	88905	A6457	10
Below Knee - L	38–42 cm (15–16.5 in.)	88906	A6457	10

Transelast® Classic

Conforming bandage



Field of Application

- Compression bandage for lymphedema and edema
- Intended for bandaging of fingers, toes, and genitals
- Dressing retention

Features/Benefits

- Soft, comfortable, close-weave gauze bandage
- Ideal narrow width for bandaging fingers and toes
- Conforms well to body contours
- Air-permeable
- Extensibility is approximately 130%

Product composition

59% polyamide, 41% viscose
Not made with natural rubber latex.

Transelast® Classic Conforming bandage

rolled, bulk packed in a cardboard box



Width	REF	Shipping Units (per box/case)
4 m (4.4 yd) stretched, 1.7 m (1.9 yd) unstretched length		
4 cm (1.6 in.)	153374	20/660
6 cm (2.4 in.)	153375	20/420

Mollelast®

Conforming Bandage



Field of Application

- Compression bandaging for
- lymphedema
- Fingers, Toes and Genital lymphedema
- Fixing foam padding
- Dressing retention

Features/Benefits

- Comfortable, open-weave gauze
- bandage
- Ideal for bandaging fingers and toes
- Conforms well to body contours
- Air permeable
- Wider sizes are excellent to secure
- foam padding from moving under bandages
- Extensibility is approximately 100%

Product composition

56% viscose, 44% polyamide

Not made with natural rubber latex.



Mollelast® Conforming bandage

rolled, bulk packed in a cardboard box

Width	REF	Shipping Units (per box/case)
4 m (4.4 yd) stretched length		
4 cm (1.6 in.)	19410	20/520
6 cm (2.4 in.)	19411	20/360
8 cm (3.2 in.)	19412	20/500
10 cm (4 in.)	19413	20/360
12 cm (4.7 in.)	19414	20/320

Raucolast®

Conforming bandage



Field of Application

- Compression bandaging for lymphedema
- Fingers, Toes and Genital lymphedema
- Fixing foam padding
- Dressing retention

Product composition

100% polyester
Not made with natural rubber latex.

Features/Benefits

- Soft, comfortable, closed-weave gauze bandage
- Ideal for bandaging fingers and toes
- Conforms well to body contours
- Air permeable
- Wider size is excellent to secure foam padding from moving under bandages
- Extensibility up to approximately 80% lengthwise



Raucolast® Conforming bandage

rolled, bulk packed in a cardboard box

Width	REF	Shipping Units (per box/case)
4 m (4.4 yd) stretched length		
4 cm (1.6 in.)	11493	20/320
6 cm (2.4 in.)	11496	20/200

Cellona® Synthetic Padding



Field of Application

- Padding under compression bandaging (lymphology and phlebology)

Features/Benefits

- Cushions the limb
- Prevents sharp indentations and/or irritations of the skin
- Ensures even distribution of pressure from compression bandages
- Self-adhesive to facilitate application and prevent slippage

- Soft and elastic; provides comfort and excellent conformability around body contours
- Helps form a smooth surface
- Air permeable
- Retains fullness with long-term wear

Product composition

100% polyester fibers
Not made with natural rubber latex.



Cellona® Synthetic padding individually-wrapped in paper

Width	REF	Shipping Units (per box/case)
3 m (3.3 yd) length		
4 cm (1.6 in.)	136338	9/90
	10690	108
6 cm (2.4 in.)	136339	6/60
	10692	72
10 cm (4 in.)	136336	4/40
	10694	48
15 cm (6 in.)	136337	3/30
	10695	36

Rosidal® soft Foam Padding Bandage



Field of Application

- Padding under compression bandaging (lymphology and phlebology)

Features/Benefits

- Reusable; washable up to 50 times
- Ensures even distribution of pressure from compression bandages
- Interlocking open pores make compression bandages slip resistant
- Durable and tear resistant
- Provides some compression and shaping of the limb when pulled slightly
- Highly breathable and permeable to moisture

- Washing temperature up to 95°C/203°F (temperatures from 40°–60°C / 104°–140°F are recommended)
- Maintains elasticity and cushioning effect even after washing

Product composition

100% polyurethane foam

Product and clips are not made with natural rubber latex.



Rosidal® soft Foam padding bandage

individually-packed in a cardboard box

Width	Thickness	REF	Shipping Units (per box/case)
2.5 m (2.7 yd) length			
10 cm (4 in.)	0.3 cm (0.12 in.)	23110	1/30
10 cm (4 in.)	0.4 cm (0.16 in.)	23111	1/24
12 cm (4.7 in.)	0.4 cm (0.16 in.)	23112	1/16
15 cm (6 in.)	0.4 cm (0.16 in.)	23113	1/14

Komprex®

Foam Rubber Pads, Bandage, Sheet



Field of Application

- Compression in treatment of edema (lymphology and phlebology)
- Increase local pressure (in combination with other compression products)
- Soften fibrosis
- Fill anatomical depressions

- Available in ready-to-use shapes, rolls and sheet
- Permeable to air and water vapor

Product composition

100% synthetic latex

Not made with natural rubber latex.

Features/Benefits

- Reusable, closed-cell, dense foam
- Easily cuts to required shape and size

Komprex® Foam rubber pads

individually-packed in a cardboard box

Style	Size	Thickness	REF	Shipping Units (per box/case)
00 Oval	9 x 6.5 cm (3.5 x 2.6 in.)	10 mm (0.4 in.)	22300	1/60
0 Kidney-shaped	9 x 5 cm (3.5 x 2 in.)	10 mm (0.4 in.)	22301	1/75
1 Kidney-shaped	12 x 6.5 cm (4.8 x 2.6 in.)	12 mm (0.5 in.)	22302	1/30

Komprex® Foam rubber bandage

rolled, individually-packed in a cardboard box

Width	Thickness	REF	Shipping Units (per box/case)
2 m (2.2 yd) length			
8 cm	5 mm (0.2 in.)	22313	1/30
8 cm	10 mm (0.4 in.)	22314	1/15

Komprex® Foam rubber sheet

rolled, individually-packed in a plastic bag

Size	Thickness	REF	Shipping Units (per box/case)
100 x 50 cm (39.4 x 20 in.)	10 mm (0.4 in.)	22320	1/10

Komprex® II

Channel Foam Sheet



Field of Application

- Compression in treatment of edema (lymphology and phlebology)
- Increase local pressure (in combination with other compression products)
- Soften fibrosis

Features/Benefits

- Special surface creates a firm, yet gentle micro-massaging effect
- Narrow strips of open cell foam material are surrounded by a skin-friendly, non-woven, breathable, non-slip padding
- Provides channeling effect
- Loosens hardened and fibrotic tissue
- Better air flow for greater patient comfort and compliance
- Washable (up to 30°C/86°F) and reusable.

Product composition

Narrow strips of polyurethane foam rubber embedded in a polyester, non-woven fabric

Not made with natural rubber latex.

To note

To avoid obstructions of lymphatic flow, the foam line must always be applied in a direction lengthwise to the extremity with the soft, comfortable side of the material facing the body.

Komprex® II Foam sheet

individually-wrapped in paper inside a plastic bag

Size	Thickness	REF	Shipping Units (per box/case)
65 x 65 cm (26 x 26 in.)	10 mm (0.4 in.)	16985	1/6

Rosidal® K

Short stretch bandage



Field of Application

- Management of venous hypertension, edema, venous leg ulceration
- Strong compression bandaging for management of lymphedema and varicose veins
- Support Bandage for immobilization and prevention or relief of swelling due to sprains and strains (sports injury)
- Immobilization of limbs

Features/Benefits

- 100% cotton; cool, skin-friendly comfort and natural moisture wicking
- Reusable; washable up to 50 times
- Provides high working and low resting pressures, safe and comfortable for day and night wear
- Textured textile structure prevent bandage slipping
- Air permeable
- Available in double width to simplify application on larger circumference thigh or lower leg
- Lengthwise extensibility (approximately 90%)
- Color: Beige
- Steam sterilizable (max. 134 C / 273.2 F)

Care instructions

- Machine washable on delicate cycle and dryable up to 50 times (max. 95 C / 203 F)
- Use mild detergents
- Do not use bleach or fabric softeners

Product composition

100% cotton

Product and bandage clips are not made with natural rubber latex.



Rosidal® K Short stretch bandage

individually-rolled and secured with a clip, packed in a cardboard box

Width	HCPCS Code	REF	Shipping Units (per box/case)
5 m (5.5 yd) length			
4 cm (1.6 in.)	A6448	90684	1/20
		22199	125
6 cm (2.4 in.)	A6448	90685	1/20
		22200	108
8 cm (3.2 in.)	A6449	90686	1/20
		22201	108
10 cm (4 in.)	A6449	90687	1/20
		22202	90
12 cm (4.7 in.)	A6449	90688	1/20
		22203	54
10 m (11 yd) length			
10 cm (4 in.)	A6449	90689	1/20
		22206	60
12 cm (4.7 in.)	A6449	90690	1/20
		22205	40

Idealbinde®

Light Compression Short Stretch Bandage



Field of Application

- Venous leg ulcer
- Venous hypertension
- Edema
- Lymphedema
- Swelling
- Used as second layer in multi-layer bandaging systems

Product composition

100% Cotton

Product and clips are not made with natural rubber latex.

Features/Benefits

- 100% cotton
- Short stretch bandage
- Textured fabric structure ensures bandage layers do not slip
- Skin-friendly
- Air permeable
- Extensibility approximately 90%
- Washable up to 95°C/203°F



Idealbinde® Short Stretch Bandage

individually-wrapped

Size	REF	Suggested HCPCS Code	Shipping Units (per box/case)
5 m (5.5 yd) stretched, 2.7 m (3 yd) unstretched length			
10 cm (4 in.)	14103	A6449	10/60
15 cm (6 in.)	14105	A6450	10/60
20 cm (8 in.)	14106	A6450	1*/35

Lenkelast®

Medium Stretch Bandage



Field of Application

- Compression bandaging (lymphology and phlebology)
- Genital lymphedema
- Pediatric lymphedema
- Breast bandaging
- Fixing layer over underlying foam padding
- Dressing retention
- Immobilization of limbs

Features/Benefits

- Soft, comfortable, easy to use all purpose bandage
- Conforms well to all body contours
- Very skin-friendly and breathable
- Durable elasticity for a wide variety of applications
- Extensibility is approximately 130%

Product composition

64% cotton, 35% polyamide, 1% spandex
Not made with natural rubber latex.



Lenkelast® Medium stretch bandage
individually-wrapped in plastic packed in a cardboard box

Width	REF	Shipping Units (per box/case)
5 m (5.5 yd) stretched length		
6 cm (2.4 in.)	39430	10/250
8 cm (3.2 in.)	39431	10/200
10 cm (4 in.)	39432	10/240
12 cm (4.7 in.)	39433	10/200
15 cm (6 in.)	39434	10/160
20 cm (7.9 in.)	39435	10/100

Durelast®

Very Short Stretch Bandage



Field of Application

- Very strong compression bandaging (lymphology and phlebology)
- Support and stress relief in acute care and sports medicine

Features/Benefits

- Firm grip fabric texture ensures bandage layers hold securely
- Very skin-friendly
- Air permeable
- Extensibility is approximately 45%
- Washable up to 95°C/203°F

Product composition

66% cotton, 34% polyamide

Product and clips are not made with natural rubber latex.



Durelast® Very short stretch bandage
individually-packed in a cardboard box

Width	REF	Shipping Units (per box/case)
5 m (5.5 yd) stretched length		
6 cm (2.4 in.)	22100	1/108
8 cm (3.2 in.)	22101	1/108
10 cm (4 in.)	22102	1/90
12 cm (4.7 in.)	22103	1/54

Dauerbinde® K

Long Stretch Bandage



Field of Application

- Compression bandaging (phlebology)
- Management of venous leg ulcers
- Support and relief of musculoskeletal system
- Management of contusions, sprains, tendon injuries and joint dislocations
- Dressing retention
- Limb immobilization
- Reduction of hematoma

- Ribbed surface in combination with high stretch-regain prevents slippage
- Compressive power can be easily adjusted due to its long stretch properties
- Durable textured fabric structure
- Air permeable
- Extensible up to 170%
- Washable up to 15 times at temperatures up to 95°C / 203°F

Features/Benefits

- Easy to apply and conforms well without creases or folds to body contours – even for those who have little bandaging experience
- Soft and skin-friendly

Product composition

58% Cotton, 24% Polyamide, 8% Elastane, 10% Viscose

Not made with natural rubber latex.

Dauerbinde® K Long stretch bandage

individually-packed in a cardboard box



Width	REF	Shipping Units (per box/case)
7 m (7.7 yd) stretched length		
6 cm (2.4 in.)	105904	1/80
8 cm (3.2 in.)	105905	1/100
10 cm (4 in.)	105900	1/80
12 cm (4.7 in.)	105906	1/80
20 cm (7.9 in.)	105901	1/30

Silkafix®

Adhesive Tape



Field of Application

- Fastens bandages, dressings, catheters, cannulas, drains or probes

Features/Benefits

- Strong, durable, satiny tape adheres well, yet removes easily
- High performance adhesive leaves no residue on skin or bandages
- Easily tears lengthwise and width-wise for customized application; no scissors needed
- Particularly suitable for sensitive skin

- Skin-friendly, hypoallergenic
- White artificial silk, non-elastic tape
- Water and dirt repellent
- Quick and easy application; dispenses from a quality, protected, plastic spool

Product composition

Rayon taffeta, white, 100% acetate; skin friendly polyacrylate adhesive

Not made with natural rubber latex.



Silkafix® Adhesive tape
rolled on a plastic dispensing reel

Width	REF	Shipping Units (per box)
5 m (5.5 yd) roll		
1.25 cm (0.5 in.)	143853	18/450

Rosidal® Lymphset

Complete kit for Lymphedema treatment



Field of Application

- Professional lymphological compression bandaging of legs and arms

Features/Benefits

- Convenient, prepackaged kit saves time and simplifies the order process and inventory management
- Contains all the necessary component required for

a professional multi-layer compression bandage

- Kits contain two complete sets of bandaging material - one to wear and one to wash for uninterrupted usage
- Requires no preparation time for specific patient care. It's ready for use at any time.

Rosidal® Lymphset Complete kit for lymphedema treatment

packed together in a cardboard box

Qty	Item	Width	Length	REF	Shipping Units
Deluxe Double Arm					
1	tg®	6.5 cm (2.6 in.)	20 m (22 yd)	136672	1 kit
20	Mollelast®	6 cm (2.4 in.)	4 m (4.4 yd)		
4	Cellona®	10 cm (4 in.)	3 m (3.3 yd)		
2	Rosidal® K	6 cm (2.4 in.)	5 m (5.5 yd)		
2	Rosidal® K	8 cm (3.2 in.)	5 m (5.5 yd)		
4	Rosidal® K	10 cm (4 in.)	5 m (5.5 yd)		
1	Silkafix®	1.25 cm (0.5 in.)	5 m (5.5 yd)		
Deluxe Double Leg					
1	tg®	8.5 cm (3.3 in.)	20 m (22 yd)	136673	1 kit
20	Mollelast®	4 cm (1.6 in.)	4 m (4.4 yd)		
4	Cellona®	10 cm (4 in.)	3 m (3.3 yd)		
6	Cellona®	15 cm (6 in.)	3 m (3.3 yd)		
2	Rosidal® K	6 cm (2.4 in.)	5 m (5.5 yd)		
2	Rosidal® K	8 cm (3.2 in.)	5 m (5.5 yd)		
4	Rosidal® K	10 cm (4 in.)	5 m (5.5 yd)		
4	Rosidal® K	12 cm (4.7 in.)	5 m (5.5 yd)		
2	Idealbinde™	15 cm (6 in.)	5 m (5.5 yd)		
2	Komprex®	9 cm (3.5 in.)	5 cm (2 in.)		
1	Silkafix®	1.25 cm (0.5 in.)	5 m (5.5 yd)		

Rosidal® TCS

Two Component System



Field of Application

- Venous leg ulcer
- Edema
- Lymphedema
- ABPI above 0.8 or with an ABPI of 0.5–0.8 under medical supervision

Features/Benefits

- Two components work together to deliver effective levels of working compression and low levels of resting compression

- Safe-Loc system provides built in safety when applied at full stretch
- Thin enough to allow the use of normal footwear to encourage compliance
- Rosidal SC (soft compression)
 - ◆ Thin, absorbent layer combines padding and skin protection in a single bandage
 - ◆ Acts as a cushion to evenly distribute pressure
 - ◆ Integrated fabric lining protects sensitive skin from irritation
 - ◆ Absorbs exudate and perspiration
- Rosidal CC (cohesive compression)
 - ◆ Durable short stretch cohesive bandage layer
 - ◆ Short stretch for effective levels of working compression
 - ◆ Cohesive to secure the compression bandage system for up to 7 days

Product composition

Rosidal SC: Polyurethane foam with integrated woven textile material made from 60% polyamide and 40% cotton

Not made with natural rubber latex.

Rosidal CC: 51% cotton, 39% polyamide, 10% polyurethane, with a hypo-allergenic synthetic cohesive micro-coating

Not made with natural rubber latex.



Rosidal® TCS Two-Component Compression System

individually-packed in a cardboard box

Component	Width	Length		REF	HCPCS Code	Shipping Units (per box/case)
		stretched	unstretched			
Rosidal® SC layer 1	10 cm (4 in.)	n/a	3.5 m (3.8 yd)	26484	A6441*	1/14
Rosidal® CC layer 2	10 cm (4 in.)	6 m (6.6 yd)	3.9 m (4.2 yd)		A6452*	

Rosidal® Sys

Short Stretch Compression System



Field of Application

- Compression therapy of venous leg ulcers

Features/Benefits

- Ready-to-use products combined in one handy set
- Cost effective with reusable padding (Rosidal® soft) and compression bandages (Rosidal® K)
- Rosidal® K short stretch bandage provides high working pressure and low resting pressure making the system tolerable at night
- System can be worn continuously for up to 7 days

- Set saves handling time and eliminates administration of individual components
- Rosidal® soft and Rosidal® K are washable up to 50 times at up to 95°C/203°F
- Mollelast® haft contains natural rubber latex

To Note

For wound treatment, a suitable dressing can be easily added so the wound is optimally cared for.



Rosidal® Sys Short-stretch compression system
packed together in a cardboard box

Qty	Item	Width	Length	REF	Shipping Units (per set/case)
1	tg®	8 cm (4 in., Size 7)	20 m (22 yd)	16984	1/4
4	Rosidal® soft (0.2 cm / 0.08 in. thickness)	10 cm (4 in.)	2 m (2.2 yd)		
2	Rosidal® K	8 cm (3.2 in.)	5 m (5.5 yd)		
2	Rosidal® K	10 cm (4 in.)	5 m (5.5 yd)		
1	Porofix®	2.5 cm (3.2 in.)	5 m (5.5 yd)		
1	Mollelast® haft	8 cm (3.2 in.)	20 m (22 yd)		

Porelast®

Short Stretch Adhesive Bandage



Field of Application

- Compression bandaging in phlebology
- Functional bandaging in sports medicine
- Support and relief with distortions, contusions and dislocations

Features/Benefits

- Lengthwise stretchable woven, 100% cotton bandage
- Fabric with porously applied adhesive is permeable to air and water vapor
- X-ray transparent
- Extensibility approximately 80%

Product composition

100% cotton with synthetic rubber adhesive (contains zinc oxide)

Not made with natural rubber latex.



Porelast® Short-stretch adhesive bandage

individually-packed in a cardboard box

Size	REF	Shipping Units (per box/case)
2.5 m (2.7 yd) unstretched length		
6 cm (2.4 in.)	30210	1/50
8 cm (3.2 in.)	30211	1/50
10 cm (4.0 in.)	30212	1/50

Panelast®

Short Stretch Adhesive Bandage



Field of Application

- Compression bandaging in phlebology
- Base layer in tape applications in sports medicine
- Support and relief with distortions, contusions and dislocations
- Follow-up treatment of fractures

Features/Benefits

- Lengthwise and widthwise stretchable woven, 100% cotton bandage
- Excellent conformability around body contours
- Fabric with porously applied adhesive is permeable to air and water vapor
- X-ray transparent
- Extensibility widthwise approximately 40%
- Extensibility lengthwise approximately 80%

Product composition

100% cotton with synthetic rubber adhesive (contains zinc oxide)

Not made with natural rubber latex.



Panelast® Short-stretch adhesive bandage

individually-packed in a cardboard box

Size	REF	Shipping Units (per box/case)
2.5 m (2.7 yd) unstretched length		
8 cm (3.2 in.)	30271	1/50

tg[®] fix

Tubular Net Bandage



Field of Application

- Retention of all types of wound dressings

Features/Benefits

- Wide-mesh, highly elastic in both directions tubular net
- Fast and easy dressing inspection
- Material is ladderproof; it can be cut open at any point without the cut tearing
- Conformability can be further increased by placing incisions in the right places without incurring the risk of fraying
- Excellent conformability around body contours
- Soft, comfortable and does not restrict movement
- Washable (up to 60°C/140°F) and reusable

Product composition

71% polyamide,
29% elastodiene (natural rubber latex)



tg[®] fix Tubular net bandage

individually-packed in a cardboard box

Size	Width		REF	Shipping Units (per box/case)
	Untretched	Stretched		
25 m (27.3 yd) stretched length				
A for finger	1.4 cm (0.5 in.)	17.5 cm (6.9 in.)	24250	1/66
B for small extremities	2.7 cm (1 in.)	40 cm (15.7 in.)	24251	1/20
C for large extremities	4.2 cm (1.7 in.)	65 cm (25.6 in.)	24252	1/16
D for large head, small trunk	7.0 cm (3 in.)	105 cm (41.3 in.)	24253	1/16
E for large trunk, hip, armpit	9.1 cm (3.6 in.)	140 cm (55.1 in.)	24254	1/12

Bandage Roller



- Quick bandage re-rolling with manually operated device
- Bandage is placed into the slot and draped over the tension bar
- Easy to use; free hand guides bandage while controlling tension
- Roller can be mounted to any flat surface with enclosed bracket
- After usage, bandage roller can be taken out of the bracket holder and stored

Bandage Roller

individually-packed in a plastic bag

REF	Shipping Units
30272	1



People.Health.Care.

Support for Life 

L&R USA INC. · 3880 W Wheelhouse Road, Milwaukee, WI 53208 USA
855-892-4140 · Fax: 414-892-4150 · inquiries@us.LRmed.com · www.Lohmann-Rauscher.com
L&R CAN INC. · 2820 Argentinia Road, Unit 2, Mississauga, ON L5N 8G4 Canada
888-836-6726 · Fax: 888-434-5347 · inquiries@ca.LRmed.com · www.Lohmann-Rauscher.ca