



People.Health.Care.

# ExoSoft™

## Comfortable circular-knit compression.



Gauntlet

Glove  
Quarter Finger

Glove  
Half Finger

Glove  
Full Finger

ExoSoft™ is a full line of off-the-shelf circular knit compression garments designed for the self-management of edema or mild lymphedema during the day or high activity. Best known for its soft, high-containment seamless gauntlet and gloves, ExoSoft offers soft, breathable comfort to support patient compliance.



### Features

- Wide off-the-shelf selection of seamless, high-containment gloves and gauntlets
- Gloves available in 3 finger lengths: Quarter, Half, and Full
- Made in USA
- Ships quickly<sup>1</sup>
- Supports mobility with its comfortable, seamless design
- Quick and easy fitting with only 2–5 measurements needed to select size for each style
- Easy garment care in washing machine and dryer
- 6-month material and manufacturing warranty

<sup>1</sup>. Orders placed before 5pm Central ship same day.



People.Health.Care.



No Seams  
Gauntlet & Glove

### ExoSoft™ Gauntlet

sold individually

<b>Compression</b>	±↑ 20–30 mmHg	
<b>Color</b>	<input type="checkbox"/> BEIGE	<input type="checkbox"/> BLACK
<b>Size</b>	<b>Circumference</b>	
	B <sup>c</sup>	C <sup>c</sup>
Small	17–19	15–17
Medium	19–21	17–19
Large	22–24	20–22
X-Large	25–27	23–25

### ExoSoft™ Glove

sold individually

<b>Compression</b>	±↑ 20–30 mmHg	
<b>Color</b>	<input type="checkbox"/> BEIGE	<input type="checkbox"/> BLACK
<b>Finger length</b>	<input type="checkbox"/> QUARTER	<input type="checkbox"/> HALF FULL
<b>Size</b>	<b>Circumference</b>	
	B <sup>c</sup>	C <sup>c</sup>
Small	17–19	15–17
Medium	19–21	17–19
Large	22–24	20–22
X-Large	25–27	23–25

### ExoSoft™ Arm

sold individually

<b>Compression</b>	±↑ 20–30 mmHg			
<b>Color</b>	<input type="checkbox"/> BEIGE	<input type="checkbox"/> BLACK		
<b>Top style</b>	<input type="checkbox"/> KNIT TOP	<input type="checkbox"/> SILICONE TOP		
<b>Size</b>	<b>Circumference</b>		<b>Length</b>	
	C <sup>c</sup>	E <sup>c</sup> G <sup>c</sup>		CG <sup>t</sup>
<b>Short</b>				
Small	15–17	24–27	27–30	< 43
Medium	17–19	27–31	31–35	
Large	20–22	31–35	35–40	
X-Large	23–25	36–41	40–45	
<b>Long</b>				
Small	15–17	24–27	27–30	> 43
Medium	17–19	27–31	31–35	
Large	20–22	31–35	35–40	
X-Large	23–25	36–41	40–45	

### ExoSoft™ Below Knee


sold in pairs

<b>Compression</b>	±↑ 15–20 mmHg 20–30 mmHg		
<b>Color</b>	<input type="checkbox"/> BEIGE	<input type="checkbox"/> BLACK	
<b>Top style</b>	<input type="checkbox"/> KNIT TOP		
<b>Toe style</b>	<input type="checkbox"/> OPEN TOE	<input checked="" type="checkbox"/> CLOSED TOE	
<b>Size</b>	<b>Circumference</b>		<b>Length</b>
	B <sup>c</sup>	C <sup>c</sup>	
Small	18–20	27–37	up to 40
Medium	20–24	29–39	up to 43
Large	24–28	32–43	up to 45
X-Large	28–31	35–45	up to 48

### ExoSoft™ Thigh High

sold in pairs

<b>Compression</b>	±↑ 15–20 mmHg 20–30 mmHg		
<b>Color</b>	<input type="checkbox"/> BEIGE	<input type="checkbox"/> BLACK	
<b>Top style</b>	<input checked="" type="checkbox"/> SILICONE TOP		
<b>Toe style</b>	<input checked="" type="checkbox"/> CLOSED TOE		
<b>Size</b>	<b>Circumference</b>		<b>Length</b>
	B <sup>c</sup>	C <sup>c</sup> G <sup>c</sup>	
Small	18–20	27–37	41–56 up to 73
Medium	20–24	29–39	45–61 up to 76
Large	24–28	32–43	50–66 up to 78
X-Large	28–31	35–45	56–71 up to 81



All measurements in centimeters.



L&R USA INC. · 3880 W Wheelhouse Road, Milwaukee, WI 53208 USA · Phone: 855-892-4140 · Fax: 414-892-4150 · Email: inquiries@us.LRmed.com