



People.Health.Care.

ReadyWrap® Full Leg **NEW**

Complete coverage, right out of the bag.



Same garments, new package

The same garments you know and trust, now available in a single package.



Easier ordering

Reduce ordering time and potential errors with five garments together in one single package.



Cost-conscious

Reimbursable under HCPCS code A6586 to improve patient access.



Each package includes:

- 1 ReadyWrap® Calf
- 1 ReadyWrap® Knee
- 1 ReadyWrap® Thigh
- 2 ReadyWrap® Liners, Thigh High



Solaris Collection™



People.Health.Care.

ReadyWrap® Full Leg

Adjustable gradient
compression wrap system


ReadyWrap® Full Leg is an adjustable gradient compression wrap system made up of the same highly-adjustable ReadyWrap garments available separately — ReadyWrap Full Leg provides an easy way to ensure your patients get the therapeutic garments they need from a single easy-to-order product.

Indications

- Lymphedema
- Lipedema
- Edema
- Venous ulcers
- Venous insufficiency
- Leg fatigue and/or tiredness

Compression level

- 30-50mmHg



Off-the-Shelf



Day



High or Low Activity

Available in

Beige

Black

- 6-month material and manufacturing warranty
- Hand wash, lay flat to dry


ReadyWrap® Full Leg

sold as a complete system including a single ReadyWrap Calf, ReadyWrap Knee, and ReadyWrap Thigh garments with two Thigh High ReadyWrap Liners

Size	S	M	L	XL	XXL	
G ^c	48–63 cm	58–73 cm	68–83 cm	78–93 cm	88–103 cm	
E ^c	38–48 cm	45–58 cm	52–68 cm	59–78 cm	66–88 cm	
D ^c	28–38 cm	33–43 cm	38–48 cm	43–58 cm	53–68 cm	
C ^c	28–38 cm	33–43 cm	38–48 cm	43–58 cm	53–68 cm	
B ^c	18–23 cm	20–28 cm	25–33 cm	30–38 cm	35–43 cm	

Length	regular	long
YG ^L	up to 75 cm	up to 85 cm


Each package includes:



1 × ReadyWrap® Calf



1 × ReadyWrap® Knee



1 × ReadyWrap® Thigh



2 × ReadyWrap® Liners,
Thigh High