



People.Health.Care.

ReadyWrap®

Clinically-effective and easy-to-use adjustable compression wrap.




ReadyWrap®



Clinically-effective and easy-to-use.

ReadyWrap® is an adjustable compression wrap for transitioning out of therapy into self-management during the daytime, nighttime, and during high or low activity.

 Off-the-Shelf or Custom*

 Day or Night

 High or Low Activity





Low-stretch material

Mimics short-stretch bandaging — the gold standard of treatment.



Color-coded straps

Provides easy step-by-step donning guidance.



50% Overlapping straps[§]

Helps prevent gapping.¹



Blocked support spine

Helps maintain garment form to support all-day wear.¹



Foam Padding

Cushioning at elbow and wrist provides comfort and evenly distributes pressure.



Modular design

Garments work effectively on their own or in combination.

Why choose ReadyWrap®?

ReadyWrap® features low-stretch straps that mimic short stretch bandaging — delivering effective compression therapy with its high working and low resting pressure.

To support patient compliance, ReadyWrap offers a unique color-coded visual guide for application and alternating strapping system to support one handed donning as needed.

With ReadyWrap, you can be certain your patients will have an effective and comfortable garment that easily adjusts to their limb and donning abilities.

Indications

- Lymphedema
- Edema
- Leg ulcers*
- Venous insufficiency*
- Leg fatigue and/or tiredness*

Available in



- 6-month warranty
- Hand wash, lay flat to dry



“To manage the swelling reduction, I placed him in bi-lateral ReadyWrap garments. He and his spouse loved that they could easily don/doff it on their own.”

Vicki Ralph, MPPA, OTR/L, CLT-LANA

Lower Extremity Garments and Accessories



ReadyWrap Toe



ReadyWrap Foot



ReadyWrap Foot CT **NEW**



ReadyWrap Foot SL

ReadyWrap® Toe

left or right orientation

Size	Circumference	
	A ^c	
Small	18.0–20.4	
Medium	20.5–23.5	
Large	23.6–26.0	
X-Large	26.1–28.5	
XX-Large	28.6–31.0	

ReadyWrap Foot

left or right orientation,
includes single Below Knee ReadyWrap Liner

Size	Circumference		Length AY ^L
	Y ^c	B ^c	
Small	28–33	18–23	up to 21
Medium	30–38	20–28	up to 23
Large	35–43	25–33	up to 26
X-Large	40–48	30–38	up to 29

ReadyWrap Foot CT **NEW**

left or right orientation,
includes single Below Knee ReadyWrap Liner

Size	Circumference			Length AY ^L
	A ^c	A ^{1c}	B ^c	
Regular				
Small	18–20.4	20–25	18–23	16–19
Medium	20.5–23.5	24–29	20–28	
Large	23.6–26	28–33	25–33	
X-Large	26.1–28.5	32–37	30–38	
XX-Large	28.6–31	36–41	35–43	
Long				
Small	18–20.4	20–25	18–23	19–22
Medium	20.5–23.5	24–29	20–28	
Large	23.6–26	28–33	25–33	
X-Large	26.1–28.5	32–37	30–38	
XX-Large	28.6–31	36–41	35–43	

ReadyWrap Foot SL

left or right orientation,
includes single Below Knee ReadyWrap Liner

Size	Circumference		Length AY ^L
	A ^{1c}	B ^c	
Regular			
Small	21–24	18–23	16–19
Medium	25–28	20–28	
Large	29–32	25–33	
X-Large	33–36	30–38	
Long			
Small	21–24	18–23	19–22
Medium	25–28	20–28	
Large	29–32	25–33	
X-Large	33–36	30–38	

ReadyWrap Extender Strap

sold individually

Size	Length
One size	10 cm



Don't see the size or configuration needed for your patient?

Our Custom Design Center can work with you to create a custom lower extremity ReadyWrap to meet your patient's unique needs.

The Custom Design Center is available
Monday–Friday, 8:00AM–6:00PM Central Time

☎ 1-414-892-5158

✉ CustomDesignCenter@us.LRmed.com

📺 Video conference

ReadyWrap
Extender Strap





ReadyWrap Calf



ReadyWrap Fusion Kit

Full Kit
Includes
liner with
compressive
foot & ankle



ReadyWrap Knee



ReadyWrap Thigh

ReadyWrap Calf

includes single Below Knee ReadyWrap Liner

Size	Circumference		Length BD ^L
	B ^C	C ^C	
Regular (Average)			
Small	18–23	28–38	up to 30
Medium	20–28	33–43	
Large	25–33	38–48	
X-Large	30–38	43–58	
XX-Large	35–43	53–68	

Long (Tall)

Small	18–23	28–38	up to 35
Medium	20–28	33–43	
Large	25–33	38–48	
X-Large	30–38	43–58	
XX-Large	35–43	53–68	

ReadyWrap Liners

sold in pairs

Size	Circumference		
	B ^C	C ^C	G ^C
Below Knee			
S–L	18–33	28–48	
XL–XXL	30–43	43–68	
Thigh High			
S–L	18–33	28–48	48–83
XL–XXL	30–43	43–68	78–103

ReadyWrap Fusion Kit **NEW**

kit includes ReadyWrap Calf garment and single ReadyWrap Fusion Liner

Size	Circumference			Length	
	A ^C	B ^C	C ^C	AY ^L	BD ^L
Regular					
Small	21–24	18–23	28–38	< 25	up to 30
Medium	25–28	20–28	33–43		
Large	29–32	25–33	38–48		
X-Large	33–36	30–38	43–58		
XX-Large	37–40	35–43	53–68		

Long

Small	21–24	18–23	28–38	≥ 25	up to 35
Medium	25–28	20–28	33–43		
Large	29–32	25–33	38–48		
X-Large	33–36	30–38	43–58		
XX-Large	37–40	35–43	53–68		

ReadyWrap Fusion Liners **NEW**

Below Knee liner with focused 8–15 mmHg in foot and ankle, sold in pairs

Size	Circumference			Length	
	A ^C	B ^C	C ^C	AY ^L	YD ^L
Regular					
Small	21–24	18–23	28–38	< 25	< 40
Medium	25–28	20–28	33–43		
Large	29–32	25–33	38–48		
X-Large	33–36	30–38	43–58		
XX-Large	37–40	35–43	53–68		

Long

Small	21–24	18–23	28–38	≥ 25	≥ 40
Medium	25–28	20–28	33–43		
Large	29–32	25–33	38–48		
X-Large	33–36	30–38	43–58		
XX-Large	37–40	35–43	53–68		

ReadyWrap Knee

includes single Thigh High ReadyWrap Liner

Size	Circumference	
	D ^C	E ^C
Small	28–38	38–48
Medium	33–43	45–58
Large	38–48	52–68
X-Large	43–58	59–78
XX-Large	53–68	66–88

ReadyWrap Thigh

includes single Thigh High ReadyWrap Liner

Size	Circumference		Length EG ^L
	E ^C	G ^C	
Regular (Average)			
Small	38–48	48–63	up to 25
Medium	45–58	58–73	
Large	52–68	68–83	
X-Large	59–78	78–93	
XX-Large	66–88	88–103	

Long (Tall)

Small	38–48	48–63	up to 30
Medium	45–58	58–73	
Large	52–68	68–83	
X-Large	59–78	78–93	
XX-Large	66–88	88–103	

Medicare Reimbursable

ReadyWrap Calf and ReadyWrap Fusion Kit are covered under code **A6545**¹.

Upper Extremity Garments



ReadyWrap
Gauntlet



**Built-in
Liner**
Supports
one-handed
donning!

ReadyWrap Arm

ReadyWrap Gauntlet

left or right orientation

Size	Length BC ^L	Circumference	
		B [°]	C [°]
Small	up to 12.5	18–20	13.5–17
Medium	up to 13.5	20–22	15.5–20
Large	up to 14.5	22–24	17.5–23
X-Large	up to 15.5	24–26	19.5–26

ReadyWrap Arm

Size	Circumference			Length CG ^L
	C [°]	E [°]	G [°]	
Short				
Small	13.5–17	19–26	22–29	up to 43
Medium	15.5–20	24–33	26–36	
Large	17.5–23	29–40	31–43	
X-Large	19.5–26	30–43	36–50	
Regular (Average)				
Small	13.5–17	19–26	22–29	up to 46
Medium	15.5–20	24–33	26–36	
Large	17.5–23	29–40	31–43	
X-Large	19.5–26	30–43	36–50	
Long (Tall)				
Small	13.5–17	19–26	22–29	up to 49.5
Medium	15.5–20	24–33	26–36	
Large	17.5–23	29–40	31–43	
X-Large	19.5–26	30–43	36–50	





People.Health.Care.

How to Order

1 Select a style

Variety of garments available for both upper and lower extremity.

2 Measure

Only 1–5 measurements needed per style.

3 Order


Available through your preferred garment supplier:

- Local DME
- National provider
- Online distributor

www.lohmann-rauscher.us/where-to-buy

Dealers can purchase direct from L&R USA via phone, fax, or online store.

Measurement reference



*The "A" measurement represents the distal end of the garment and varies by style.

* Lower extremity styles only.

† For presence of open wound. Refer to LCD requirement for more information.

§ ReadyWrap Calf, Knee and Thigh styles only.

Sources

1. Ehmann, S., and Bock, K. (Jan. 2018), "Case report to demonstrate the need for selection criteria for optimal adjustable Velcro wrap prescription", The Journal of Advanced Wound Care, North American Supplement, vol. 27, no. 1

